



Tapes & Labels

FOR THE BLOOD BANK





Shamrock's Mission

Our mission is to supply quality, innovative labeling systems to the healthcare industry that help protect patient safety, fulfill modern technological needs and remain cost effective.



Environmental Policy

It is Shamrock's policy to provide a safe and healthy workplace for our employees and to be an environmentally responsible neighbor in our community. From product manufacture through shipping, in both office and production facilities, Shamrock supports comprehensive programs designed to positively impact the safety of our population and natural resources.



Joint Commission National Patient Safety Goals for Nursing Sustained with Shamrock Labels

- Goal 1.** Improve the accuracy of patient identification.
- Goal 2.** Improve the effectiveness of communication among caregivers.
- Goal 3.** Improve the safety of using medications.
- Goal 7.** Reduce the risk of health care-associated infections.
- Goal 15.** Identify safety risks inherent in the patient population.



Health Insurance Portability and Accountability

Prevent unauthorized access to individually identifiable health information.

Blood Bank

ISBT 128 LABELS FOR THE BLOOD BAG	Blank & Pre-Printed Overlay Labels	BB-2
	Prep Labels for Frozen Blood Bags	BB-2
ISBT 128 BLANK THERMAL LABELS	Blank Thermal Labels	BB-3
	Thermal Label Materials	BB-4
	Ribbons for Thermal Transfer Printers	BB-5
	Thermal Transfer Label/Ribbon Kits	BB-5
ISBT 128 PRE-PRINTED LABELS	Primary & Satellite Container Labels	BB-6
	Blood Product Labels	
	Whole Blood	BB-6
	Red Blood Cells	BB-7
	Plasma	BB-11
	Platelets	BB-13
	Cryoprecipitated AHF	BB-14
	Blood Group Labels	BB-15
	Auxiliary Labels	
	Facility Labels	BB-16
	Intended Recipient Labels	BB-17
	Donation ID Labels	BB-17
	Date & Time Labels	BB-17
	Special Testing Labels	BB-17
ISBT 128 DIN NUMBER LABEL SETS ISBT	Pre-Printed & Blank Piggyback Labels	BB-18
128 CELLULAR THERAPY LABELS	Blank Thermal Labels	BB-20
	Pre-Printed Special Message Label	BB-21
ISBT 128 TISSUE BANK LABELS	Blank Thermal Labels	BB-21
CODABAR LABELS	Blood Product Labels	BB-22
	Base Labels	BB-24
	Auxiliary Labels	BB-24
CONSECUTIVELY NUMBERED LABELS	Bar Coded Number Label Sets	BB-25
	Piggyback Number Labels	BB-26
	Custom Numbering Systems	BB-26
	Tape & Pre-Cut Number Labels	BB-26
COMMUNICATION LABELS	Pre-Printed Informational Labels	BB-27
	Biohazard Warning Labels	BB-33
	Shipping Labels	BB-33
AWARD STICKERS	Award Stickers	BB-34
BLOOD TAGS	Blood Tags	BB-34
PRODUCTS FOR COLOR CODING	Dot Labels	BB-35
	Labeling Tape & Tape Dispensers	BB-36

Specialty Labels for ISBT 128



Blood Bag Labels Blank and Pre-Printed

Stock blank and pre-printed labels are designed to meet exacting ICCBBA guidelines. Label materials are glove friendly and compositionally compliant with FDA 21 CFR 175.105.

5/8" x 2"
SB128-SCR Pre-printed
SB128-582 Blank

Accurate Blood Center Anywhere, USA
FDA Registration Number _____
US License Number _____

Properly Identify Intended Recipient. See Circular of Information for indications, contraindications, cautions and methods of infusion. This product may transmit infectious agents.
Rx Only

VOLUNTEER DONOR

1 1/8" x 2"
SB128-6 Pre-printed
SB128-182 Blank

RED BLOOD CELLS
ADENINE-SALINE (AS-5) ADDED

From 500 mL CPD Whole Blood Store at 1 to 6 C

2" x 2"
SB128-20-59 Pre-printed
SB128-22 Blank

PLACE DONATION IDENTIFICATION NUMBER HERE

DO NOT TRANSFUSE UNLESS ABO LABEL APPLIED HERE

ANTICOAGULANT CITRATE PHOSPHATE DEXTROSE SOLUTION, USP

63 mL Anticoagulant Citrate Phosphate Dextrose Solution USP for collection of 450 mL Whole Blood. Each 63 mL of anticoagulant contains 1.66 g Sodium Citrate (dihydrate) USP, 1.61 g Dextrose (monohydrate) USP, 1.88 mg Citric Acid (anhydrous) USP and 140 mg Monobasic Sodium Phosphate USP. pH may have been adjusted with Sodium Hydroxide.
CAUTION: Refer to instructions for use

Container Makers, Inc
Somewhere, USA

2" x 2"
SB128-AP62 Pre-printed
SB128-22 Blank

Expiration Date
0180312359
31 JAN 2016

5/8" x 2"
SB128-DATE Pre-printed
SB128-582 Blank

4 1/4" x 4" Manufacturer's Base Label

5/8" x 2"
SB128-TEST-8 Pre-printed
SB128-582 Blank

5/8" x 2"
SB128-TEST-8 Pre-printed
SB128-582 Blank

N0008
Negative for antibodies to CMV

Prep Labels to Facilitate Over-Labeling of Frozen Bags

Translucent Prep Labels are applied to the base label before freezing. Upon removal a special adhesive remains to facilitate application of overlay label.

Prep Label is placed over base label, leaving DIN number uncovered.

Pull Tab ensures easy removal of 4" x 4" prep label when ready to over label.

Prep Label won't flag or fall off in extreme temperatures.

Special Adhesive layer remains when Prep Label is removed, allowing over labeling as low as -65°F (-54°C).

Base Prep Label
SB128-44-TPIRG
4" x 4" (no liner slit)
500 labels/roll

Base Prep Label
SB128-44-TPIRGS
4" x 4" (with liner slit for pull tab)
500 labels/pack

Double Quadrant Prep Label
SB128-24-TTP
2" x 4"
500 labels/roll

Single Quadrant Prep Label
SB128-22-TTP
2" x 2"
500 labels/roll

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




Specialty Thermal Labels for ISBT 128



Shamrock offers top quality blank labels precisely formatted for ISBT 128 labeling software. Standard labels are available for on-demand Thermal Direct Printing using a heated printhead. Labels may also be ordered in Thermal Transfer material for printers using ribbons and in a premium material for more demanding procedures. See the chart below for material specifications and ordering information.

LINER SLIT OPTION



Base Label 4" x 4"
SB128-44-P (perfs thru label)
SB128-44-PL (perfs thru label & liner)
1,400 per roll


Base Label 4" x 4"
SB128-44 (no liner slit)
SB128-44-S (liner slit)
800 per roll

Base Label 4" x 4"
SB128-44-DIN-P (perfs thru label)
SB128-44-DINPL (perfs thru label & liner)
1,400 per roll


Base Label 4" x 4"
SB128-44-DIN (no liner slit)
SB128-44-DINS (liner slit)
1,400 per roll




SB128-3252 3 1/4" x 2"
1,400 labels per roll
For Facility and Product overlay in 2 left quadrants under DIN.




Base Label 4" x 4"
SB128-44-IRGP (perfs between labels)
SB128-44-IRGPS (perfs between labels, liner slit)
800 per roll



Base Label 4" x 4"
SB128-44-IRG (no liner slit)
SB128-44-IRGS (liner slit)
800 per roll



SB128-24 (no liner slit) 2" x 4"
SB128-24-S (liner slit)
1,600 labels per roll
For ABO, Date and Testing overlay in 2 right quadrants or for Product, Date and Testing overlay in 2 bottom quadrants.




SB128-12 1" x 2"
3,000 labels per roll
For overlay on 1/2 of any quadrant.




SB128-582 5/8" x 2"
4,500 labels per roll
For overlay of Donation ID No., Date or Testing sections.



SB128-22 2" x 2"
1,600 labels per roll
For ABO or for Product quadrant overlay.



THERM-466R 1" x 2 1/2"
3,000 per roll
Intended Recipient Label for labeling bag.



SB128-182 1 1/8" x 2"
2,750 labels per roll
For overlay on 1/2 of any quadrant.

SPECIALTY BLANK THERMAL LABEL MATERIALS FOR THE BLOOD BANK					
Material	Description	Application Temperature	Temperature Range	Cat. # Suffix	Cat. # Example
FOR DIRECT PRINTING WITH HEATED PRINTHEAD					
Thermal Direct	Paper, Permanent Adhesive	0°C	-54°C to 93°C	no suffix	SB128-44
Kimdura Direct	Film, Permanent Adhesive Solvent Resistant-Won't Smear	-10°C	-54°C to 148°C	-KD	SB128-44-KD
FOR TRANSFER PRINTING USING A RIBBON					
Thermal Transfer	Paper, Permanent Adhesive	0°C	-54°C to 200°C	-T	SB128-44-T
Kimdura Transfer	Film, Permanent Adhesive Solvent Resistant-Won't Smear	-12°C	-54°C to 148°C	-KT	SB128-44-KT

For additional premium materials see page BB-4 or contact Shamrock for complete information.

Specialty Thermal Label Materials



Standard Thermal Label Materials

MATERIAL	DESCRIPTION
Thermal Direct	Permanent adhesive won't stick to gloves. May be applied at 32°F (0°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C).
Thermal Transfer	Bright white in color and smudge resistant. Permanent adhesive won't stick to gloves. May be applied at 32°F (0°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C). Sticks to most surfaces including plastic and curved.
<p>Stock Thermal Labels are supplied on 3" cores for Thermal Direct printing unless otherwise specified. To order stock labels in Thermal Transfer material, please add the suffix "-T" to the catalog number.</p>	

Thermal Label Materials for Special Applications

MATERIAL	DESCRIPTION
Thermal Direct IR Readable	Smudge resistant, perfect for curved and plastic surfaces. Glove free permanent adhesive. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).
Thermal Direct Tight Release	Smudge resistant, ideal for serpentine feed printers. Glove free permanent adhesive. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).
Thermal Direct Removable	Smudge resistant material. Glove free removable adhesive leaves no residue. May be applied at 32°F (0°C), holds in temperatures ranging from -40°F to 200°F (-40°C to 93°C).
Thermal Direct Kimdura®	Solvent resistant, imprint won't smear or flake. Glove free permanent adhesive. May be applied at 10°F (-12°C), holds in temperatures ranging from -65°F to 300°F (-54°C to 148°C).
Thermal Direct Piggyback	Holds tightly to a wide variety of surfaces, won't stick to gloves. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).
Transfer Smudgeproof Kimdura®	Glove free, solvent resistant, aggressive permanent adhesive. May be applied at 10°F (-12°C), holds in temperatures ranging from -65°F to 300°F (-54°C to 148°C).
Thermal Transfer Matte Cold-Stick	Glove free, permanent adhesive. May be applied at -10°F (-23°C), holds in temperatures ranging from -65°F to 200°F (-53°C to 93°C). Good initial adhesion that becomes permanent after 24 hours.
Thermal Transfer Cryo-Stick	Glove free, permanent adhesive for cryogenic testing and storage applications. May be applied at 50°F (10°C), holds in temperatures ranging from -112°F to 176°F (-80°C to 80°C).
Thermal Transfer Piggyback	Bright white, smudge resistant material won't stick to gloves, holds tightly to most surfaces. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).
Transfer PolyFlex Piggyback	High performance, clear carrier sheet, excellent flexibility. Glove free permanent adhesive. May be applied at 40°F (5°C), holds in temperatures ranging from -40°F to 302°F (-40°C to 150°C).
Transfer Smudgeproof Kimdura® Piggyback	Resists moisture and chemicals. Glove free permanent adhesive. May be applied at -20°F (-29°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C).

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Specialty Ribbons for Thermal Printers

WAX/RESIN RIBBONS

- Bar Coded Labels
- Extended Use Labels

Wax/resin ribbons provide superior durability and smudge resistance for softer synthetic materials.

Catalog No.	Printer	W x L	Rolls/Ctn
RIBBON-209	Brady IP300	3.15" x 984'	24
RIBBON-210	Brady IP300	3.54" x 1476'	6

RIBBON-12	Datamax	3" x 1345'	12
RIBBON-144	Datamax	4.25" x 1181'	6
RIBBON-192	Datamax	4.5" x 1476'	24
RIBBON-120	Datamax I Class	4" x 1181'	6
RIBBON-177	Datamax I Class	4.33" x 1476'	24

RIBBON-178	Intermec Easy Coder 3100	4.09" x 509'	6
RIBBON-181	Intermec PM4i	4.17" x 1499'	6

RIBBON-208	Wasp	4.33" x 984'	24
------------	------	--------------	----

RIBBON-214	Sato-M84	4.33" x 1345'	6
------------	----------	---------------	---

RIBBON-149	Zebra	4.33" x 1476'	6
RIBBON-197	Zebra	4.33" x 1476'	6
RIBBON-102	Zebra 105Se, 90Xi	3.27" x 1476'	6

RIBBON-185	Zebra GX420t	4.33" x 243'	24
RIBBON-200	Zebra GX430t	4.33" x 243'	1
RIBBON-218	Zebra GX430t	3.29" x 243'	24
RIBBON-217	Zebra S4M	2.5" x 984'	36
RIBBON-54	Zebra S500, S600	3.15" x 984'	6

RIBBON-97	Zebra S600	4.33" x 1476'	6
RIBBON-187	Zebra TLP 2824	2.2" x 243'	36
RIBBON-161	Zebra TLP 2844	2.5" x 243'	12
RIBBON-164	Zebra TLP 2844	2.5" x 243'	12
RIBBON-212	Zebra TLP 3844-Z	2.5" x 243'	12
RIBBON-107	Zebra TLP 2742, 2242	3.29" x 299'	12

WAX/RESIN RIBBONS Cont.

Catalog No.	Printer	W x L	Rolls/Ctn
RIBBON-189	Zebra TLP 2844	4.33" x 243'	12
RIBBON-143	Zebra TLP 2844	4.33" x 243'	12
RIBBON-110	Zebra TLP 2844, TLP 3844-Z	4.33" x 243'	36

RIBBON-203	Zebra TLP 2844	4.33" x 243'	12
RIBBON-109	Zebra Z4M	3.5" x 1476'	6
RIBBON-112	Zebra Z4M	4" x 1476'	6
RIBBON-169	Zebra ZM400	1.57" x 1476'	6
RIBBON-173	Zebra ZM400	4" x 1476'	6

RESIN RIBBONS


- Extreme Environments
- Asset Labels

Resin ribbons provide superior smudge and smear resistance for printing on polyester, vinyl and other synthetic materials.

Catalog No.	Printer	W x L	Rolls/Ctn
RIBBON-70	Intermec 3400, 8646, 3440	4.09" x 509'	24
RIBBON-215	Intermec PM4i	4.33" x 984'	12

RIBBON-196	Zebra	2" x 1476'	36
RIBBON-156	Zebra	4.33" x 1476'	24
RIBBON-206	Zebra GX420t	4.33" x 243'	24
RIBBON-221	Zebra GX420t	4.33" x 243'	24
RIBBON-195	Zebra GX420t, GX-430t	2.5" x 243'	36

RIBBON-204	Zebra GX430t	2.2" x 243'	48
RIBBON-198	Zebra S4M	3.15" x 984'	24
RIBBON-103	Zebra S600	4.33" x 1476'	6
RIBBON-67	Zebra S500	4.33" x 984'	6
RIBBON-184	Zebra TLP 2844	2.5" x 243'	48
RIBBON-50	Zebra TLP 2844	3.29" x 243'	12
RIBBON-160	Zebra TLP 2844	4.33" x 243'	36
RIBBON-170	Zebra TLP 2844	2.5" x 243'	48
RIBBON-186	Zebra Z4M	4.33" x 1476'	6
RIBBON-219	Zebra ZM400	2.52" x 1476'	6



WAX-BASED RIBBONS

- Bar Coded Labels
- Patient Charge Labels

Premium wax-based film ensures high performance on a wide variety of facestocks.

Catalog No.	Printer	W x L	Rolls/Ctn
RIBBON-145	Brady	2.36" x 984'	2
RIBBON-213	Brady BBP-11	4.33" x 984'	1

RIBBON-211	Cognitive Cxi	4.25" x 459'	36
------------	---------------	--------------	----

RIBBON-49	Datamax	2.52" x 1181'	36
RIBBON-201	Datamax	3.5" x 1181'	6
RIBBON-36	Datamax	4" x 1346'	24
RIBBON-142	Datamax M Class	4.25" x 1181'	6

RIBBON-72	Fargo Prodigy	4" x 1108'	6
-----------	---------------	------------	---

RIBBON-44	Intermec 3400, 8646	4.09" x 509'	6
-----------	---------------------	--------------	---


RIBBON-205	Intermec PC41	2.5" x 299'	12
------------	---------------	-------------	----

RIBBON-42	Zebra	3.54" x 1476'	6
RIBBON-105	Zebra	4.33" x 1476'	6
RIBBON-122	Zebra 5105 SL	4.17" x 1476'	6
RIBBON-180	Zebra GX420t	4.33" x 243'	24
RIBBON-207	Zebra GX420t	4.33" x 243'	12
RIBBON-41	Zebra S600	3.5" x 1476'	6
RIBBON-194	Zebra Z4M	4" x 1476'	24
RIBBON-168	Zebra TLP 2844	3.29" x 243'	36

Specialty Label Kits for Thermal Transfer Printing

- ISBT 128 series Labels • Black Ink Printer Ribbons • Clean Card

Shamrock can supply any Thermal Transfer ISBT 128 Label matched with the Printer Ribbon that will best meet the demands of that label's application. Let us know the label needed, the number of rolls to be ordered and the label's intended use. Shamrock will provide the correct ribbon to match your application, in the quantity needed for the number of labels ordered, plus a handy Clean Card to help maintain your printer's guide paths, rollers and printheads. Contact us for complete information and ordering assistance.



Specialty ISBT 128 Pre-Printed Labels



Shamrock Blood Bank Labels for ISBT 128 are printed in compliance with the US Industry consensus Standard for the Uniform Labeling of Blood Components Using ISBT 128. They feature adhesive that is compositionally compliant with FDA 21 CFR 175.105. Labels permanent adhesive holds securely during use and stands up to storage in refrigerator or freezer.

Primary and Satellite Container Labels

Blood Container Base Labels measure 4" x 4" and are supplied 500 labels per roll. Labels shown at 70% of actual size.

PLACE DONATION IDENTIFICATION NUMBER HERE **DO NOT TRANSFUSE UNLESS ABO LABEL APPLIED HERE**

Properly Identify Intended Recipient.
See Circular of Information for indications, contraindications, cautions and methods of infusion. This product may transmit infectious agents.
Rx Only
VOLUNTEER DONOR

CAUTION: Refer to instructions for use

SB128-BAOS 4" x 4"

PLACE DONATION IDENTIFICATION NUMBER HERE **DO NOT TRANSFUSE UNLESS ABO LABEL APPLIED HERE**

Accurate Blood Center
Anywhere, USA
FDA Registration Number _____
US License Number _____

Properly Identify Intended Recipient.
See Circular of Information for indications, contraindications, cautions and methods of infusion. This product may transmit infectious agents.
Rx Only
VOLUNTEER DONOR

CAUTION: Refer to instructions for use

Second Blood Center
Anywhere, USA
FDA Registration Number _____
US License Number _____

SB128-BAOS (Volunteer) 4" x 4"
SB128-BAO-P (Paid Donor) 4" x 4"

Product Labels for ISBT 128

Product Labels are placed in the lower left quadrant of the Base Label. Labels, measuring 2" x 2" and supplied 500 labels per roll, are shown below at 70% of actual size. These are only a sampling of the Product Labels available.

WHOLE BLOOD

E0053100 AUTOLOGOUS

WHOLE BLOOD

Approx 450 mL plus 63 mL CPDA-1
Store at 1 to 6 C

SB128-10-1-1 Autologous

E0009V00

WHOLE BLOOD

Approx 450 mL plus 63 mL CPD
Store at 1 to 6 C

SB128-10-6

E0023V00

WHOLE BLOOD

Approx 500 mL plus 70 mL CPD
Store at 1 to 6 C

SB128-10-2

E0063X00 AUTOLOGOUS

**WHOLE BLOOD
LEUKOCYTES REDUCED**

From 450 mL plus 63 mL CPDA-1
Store at 1 to 6 C

SB128-10-7-X Autologous

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




Shamrock Blood Bank Labels for ISBT 128 are printed in compliance with the US Industry Consensus Standard for the Uniform Labeling of Blood and Blood Components Using ISBT 128. They feature adhesive that is compositionally compliant with FDA 21 CFR 175.105.




RED BLOOD CELLS




SB128-20-3




SB128-20-3-1 Autologous




SB128-20-20




SB128-20-49




SB128-20-17




SB128-20-4




SB128-20-21




SB128-20-22




SB128-20-59




SB128-20D-60A Part A0




SB128-20D-57A Part A0




SB128-20D-57B Part B0




SB128-20D-61A Open System, Part A0



SB128-20D-26A Part A0



SB128-20D-26B Part B0



SB128-20D-52A Part A0

Specialty Product Labels for ISBT 128



Blood Product Labels measure 2" x 2" and are supplied 500 labels per roll.
Shown at 70% of actual size. If a scannable proof is needed, please contact Shamrock.

RED BLOOD CELLS


 RED BLOOD CELLS ADENINE-SALINE (AS-5) ADDED DIVIDED LEUKOCYTES REDUCED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20D-52B Part B0	 RED BLOOD CELLS ADENINE-SALINE (AS-5) ADDED IRRADIATED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20-64	 RED BLOOD CELLS ADENINE-SALINE (AS-1) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20-60	 RED BLOOD CELLS ADENINE-SALINE (AS-1) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20-65 Open System
 RED BLOOD CELLS ADENINE-SALINE (AS-3) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL From 450 mL CP2D Whole Blood Store at 1 to 6 C SB128-20-35	 APHERESIS RED BLOOD CELLS ADENINE-SALINE (AS-3) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL containing approx _____ mL CP2D Store at 1 to 6 C SB128-20-55	 RED BLOOD CELLS ADENINE-SALINE (AS-5) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL From 450 mL CPD Whole Blood Store at 1 to 6 C SB128-20-36	 RED BLOOD CELLS ADENINE-SALINE (AS-5) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20-44
 RED BLOOD CELLS ADENINE-SALINE (AS-3) ADDED IRRADIATED LEUKOCYTES REDUCED _____ mL From 500 mL CP2D Whole Blood Store at 1 to 6 C SB128-20-57	 RED BLOOD CELLS ADENINE-SALINE (AS-1) ADDED LEUKOCYTES REDUCED _____ mL From 450 mL CPD Whole Blood Store at 1 to 6 C SB128-20-2	 RED BLOOD CELLS ADENINE-SALINE (AS-1) ADDED LEUKOCYTES REDUCED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20-26	 RED BLOOD CELLS ADENINE-SALINE (AS-3) ADDED LEUKOCYTES REDUCED _____ mL From 450 mL CP2D Whole Blood Store at 1 to 6 C SB128-20-27
 RED BLOOD CELLS ADENINE-SALINE (AS-3) ADDED LEUKOCYTES REDUCED _____ mL From 500 mL CP2D Whole Blood Store at 1 to 6 C SB128-20-48	 RED BLOOD CELLS ADENINE-SALINE (AS-5) ADDED LEUKOCYTES REDUCED _____ mL From 450 mL CPD Whole Blood Store at 1 to 6 C SB128-20-28	 RED BLOOD CELLS ADENINE-SALINE (AS-5) ADDED LEUKOCYTES REDUCED _____ mL From 500 mL CPD Whole Blood Store at 1 to 6 C SB128-20-52	 RED BLOOD CELLS ADENINE-SALINE (AS-1) ADDED LOW VOLUME Approx _____ mL from _____ mL CPD Whole Blood Store at 1 to 6 C SB128-20-47

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




RED BLOOD CELLS




SB128-20D-67A Part A0




SB128-20D-67B Part B0




SB128-20D-6A Part A0




SB128-20D-6AA Part Aa




SB128-20D-6AB Part Ab




SB128-20D-6AC Part Ac




SB128-20D-6B Part B0




SB128-20D-6C Part C0




SB128-20D-6D Part D0




SB128-20D-24A Part A0




SB128-20D-56A Part A0




SB128-20D-56B Part B0




SB128-20D-56C Part C0



SB128-20-6



SB128-20-68




SB128-20-70

Specialty Product Labels for ISBT 128



Blood Product Labels measure 2" x 2" and are supplied 500 labels per roll.
Shown at 70% of actual size. If a scannable proof is needed, please contact Shamrock.


RED BLOOD CELLS


E0274V00

RED BLOOD CELLS
IRRADIATED
LEUKOCYTES REDUCED

From 500 mL CP2D Whole Blood
Store at 1 to 6 C


SB128-20-71


E0209V00

RED BLOOD CELLS
LEUKOCYTES REDUCED

From 450 mL CPDA-1 Whole Blood
Store at 1 to 6 C


SB128-20-24


E0209D00 DIRECTED

RED BLOOD CELLS
LEUKOCYTES REDUCED

From 450 mL CPDA-1 Whole Blood
Store at 1 to 6 C


SB128-24-D Directed


E0164V00

RED BLOOD CELLS
LEUKOCYTES REDUCED

From 450 mL CPD Whole Blood
Store at 1 to 6 C


SB128-20-23


E0161V00

RED BLOOD CELLS
LEUKOCYTES REDUCED

From 500 mL CPD Whole Blood
Store at 1 to 6 C


SB128-20-72


E0259V00

RED BLOOD CELLS
LEUKOCYTES REDUCED

From 450 mL CP2D Whole Blood
Store at 1 to 6 C


SB128-20-25


E5017V00

RED BLOOD CELLS
LOW VOLUME

Approx. _____ mL from
_____ mL CPD Whole Blood
Store at 1 to 6 C


SB128-20-18


E4536V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED

_____ mL containing approx. _____ mL
ACD-A
Store at 1 to 6 C


SB128-20-75


E4546V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED

_____ mL containing approx. _____ mL
ACD-A
Store at 1 to 6 C
1st Container


SB128-20-76 1st Container


E4547V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED

_____ mL containing approx. _____ mL
ACD-A
Store at 1 to 6 C
2nd Container


SB128-20-77 2nd Container


E0668V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
DIVIDED
IRRADIATED
LEUKOCYTES REDUCED

_____ mL containing approx. _____ mL
_____ mL CP2D
Store at 1 to 6 C
1st Container


SB128-20D-51A Part A0
1st Container


E4526V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-1) ADDED
IRRADIATED
LEUKOCYTES REDUCED

_____ mL containing approx. _____ mL
ACD-A
Store at 1 to 6 C


SB128-20-80


E4538V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
IRRADIATED
LEUKOCYTES REDUCED

_____ mL containing approx. _____ mL
ACD-A
Store at 1 to 6 C


SB128-20-82


E0668V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
IRRADIATED
LEUKOCYTES REDUCED

_____ mL containing approx. _____ mL
CP2D
Store at 1 to 6 C
1st Container


SB128-20-51 1st Container


E0669V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
IRRADIATED
LEUKOCYTES REDUCED

_____ mL containing approx. _____ mL
CP2D
Store at 1 to 6 C
2nd Container

SB128-20-53 2nd Container


E4532V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-1) ADDED
LEUKOCYTES REDUCED

_____ mL containing approx. _____ mL
ACD-A
Store at 1 to 6 C
1st Container

SB128-20-45 1st Container

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Shamrock Blood Bank Labels for ISBT 128 are printed in compliance with the US Industry Consensus Standard for the Uniform Labeling of Blood and Blood Components Using ISBT 128. They feature adhesive that is compositionally compliant with FDA 21 CFR 175.105.



RED BLOOD CELLS



E4533V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-1) ADDED
LEUKOCYTES REDUCED

_____ mL containing approx _____ mL

ACD-A
Store at 1 to 6 C
2nd Container

SB128-20-46 2nd Container



E4531V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-1) ADDED
LEUKOCYTES REDUCED

_____ mL containing approx _____ mL

ACD-A
Store at 1 to 6 C

SB128-20-81



E0678V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
LEUKOCYTES REDUCED

_____ mL containing approx _____ mL

CP2D
Store at 1 to 6 C

SB128-20-79



E0685V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
LEUKOCYTES REDUCED

_____ mL containing approx _____ mL

CP2D
Store at 1 to 6 C
1st Container

SB128-20-54 1st Container



E0686V00

**APHERESIS
RED BLOOD CELLS**
ADENINE-SALINE (AS-3) ADDED
LEUKOCYTES REDUCED

_____ mL containing approx _____ mL

CP2D
Store at 1 to 6 C
2nd Container

SB128-20-83 2nd Container



E5160V00

**WASHED
RED BLOOD CELLS**

_____ mL

Store at 1 to 6 C

SB128-20-29 Open System



E0472V00

**WASHED
RED BLOOD CELLS**
IRRADIATED

From 450 mL Whole Blood

Store at 1 to 6 C

SB128-20-84 Open System


E5169V00


**WASHED
RED BLOOD CELLS**
LEUKOCYTES REDUCED

_____ mL

Store at 1 to 6 C

SB128-20-37 Open System

PLASMA



E2553V00

**PLASMA
CRYOPRECIPITATE
REDUCED**

_____ mL from CPD Whole Blood

Store at -18 C or colder

SB128-40-11



E2555V00

**PLASMA
FROZEN WITHIN 24 HOURS
AFTER PHLEBOTOMY**

_____ mL from CPD Whole Blood

Store at -18 C or colder

SB128-40-5



E2619V00

**PLASMA
FROZEN WITHIN 24 HOURS
AFTER PHLEBOTOMY**

_____ mL from CP2D Whole Blood

Store at -18 C or colder

SB128-40-10



E0707V00

FRESH FROZEN PLASMA

_____ mL from CPDA-1 Whole Blood

Store at -18 C or colder

SB128-41-2



E0701V00

FRESH FROZEN PLASMA

_____ mL from CPD Whole Blood

Store at -18 C or colder

SB128-41-6



E0707VA0

**FRESH FROZEN PLASMA
DIVIDED**

_____ mL from CPDA-1 Whole Blood

Store at -18 C or colder

SB128-41D-2A Part A0



E0701VA0

**FRESH FROZEN PLASMA
DIVIDED**

_____ mL from CPD Whole Blood

Store at -18 C or colder

SB128-41D-6A Part A0


E0701VB0

**FRESH FROZEN PLASMA
DIVIDED**

_____ mL from CPD Whole Blood

Store at -18 C or colder


SB128-41D-6B Part B0

Specialty Product Labels for ISBT 128




Blood Product Labels measure 2" x 2" and are supplied 500 labels per roll.
Shown at 70% of actual size. If a scannable proof is needed, please contact Shamrock.


PLASMA




SB128-41D-6C Part C0




SB128-41D-6D Part D0




SB128-41-9




SB128-41-11 1st Container




SB128-41-12 2nd Container




SB128-41D-9A Part A0




SB128-41D-9B Part B0




SB128-41D-9C Part C0




SB128-41-5




SB128-41-3




SB128-40-8




SB128-40-24




SB128-40-17



SB128-40-2



SB128-40-14




SB128-40-21

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




PLASMA


 E2700V00
THAWED PLASMA
CRYOPRECIPITATE REDUCED
 _____ mL from CPD Whole Blood
 Store at 1 to 6 C


SB128-40-20


 E2736V00
THAWED PLASMA
CRYOPRECIPITATE REDUCED
 _____ mL from CP2D Whole Blood
 Store at 1 to 6 C


SB128-40-22


 E2696V00
THAWED PLASMA
FROZEN WITHIN 24 HOURS AFTER PHLEBOTOMY
 _____ mL from CPD Whole Blood
 Store at 1 to 6 C


SB128-40-9 Open System


 E2701V00
THAWED PLASMA
FROZEN WITHIN 24 HOURS AFTER PHLEBOTOMY
 _____ mL from CPD Whole Blood
 Store at 1 to 6 C


SB128-40-15


 E2732V00
THAWED PLASMA
FROZEN WITHIN 24 HOURS AFTER PHLEBOTOMY
 _____ mL from CP2D Whole Blood
 Store at 1 to 6 C


SB128-40-18 Open System


 E2737V00
THAWED PLASMA
FROZEN WITHIN 24 HOURS AFTER PHLEBOTOMY
 _____ mL from CP2D Whole Blood
 Store at 1 to 6 C

SB128-40-19



 E1237V00
THAWED APHERESIS FRESH FROZEN PLASMA
 _____ mL containing approx _____ mL
 ACD-A
 Store at 1 to 6 C

SB128-41-15


 E5275V00
THAWED POOLED PLASMA
FROZEN WITHIN 24 HOURS AFTER PHLEBOTOMY
 _____ mL
 Number of units in pool _____
 From CPD Whole Blood
 Store at 1 to 6 C

SB128-40-23 Open System


PLATELETS


 E2833V00
PLATELETS
 Approx 45-65 mL
 From 450 mL CPDA-1 Whole Blood
 Store at 20 to 24 C


SB128-50-1


 E2807V00
PLATELETS
 Approx 45-65 mL
 From 450 mL CPD Whole Blood
 Store at 20 to 24 C

SB128-50-7


 E2824V00
PLATELETS
 Approx 45-65 mL
 From 500 mL CPD Whole Blood
 Store at 20 to 24 C

SB128-50-23


 E2866V00
PLATELETS
LEUKOCYTES REDUCED
 Approx 45-65 mL
 From 500 mL CP2D Whole Blood
 Store at 20 to 24 C


SB128-50-24


 E2940V00
APHERESIS PLATELETS
 _____ mL containing approx _____ mL ACD-A
 Store at 20 to 24 C

SB128-50-8


 E3046VA0
APHERESIS PLATELETS
DIVIDED IRRADIATED LEUKOCYTES REDUCED
 _____ mL containing approx _____ mL
 ACD-A
 Store at 20 to 24 C

SB128-50D-3A Part A0


 E3056VA0
APHERESIS PLATELETS
DIVIDED IRRADIATED LEUKOCYTES REDUCED
 _____ mL containing approx _____ mL
 ACD-A
 Store at 20 to 24 C
 1st Container

SB128-50D-20A Part A0
1st Container


 E3057VA0
APHERESIS PLATELETS
DIVIDED IRRADIATED LEUKOCYTES REDUCED
 _____ mL containing approx _____ mL
 ACD-A
 Store at 20 to 24 C
 2nd Container

SB128-50D-21A 2nd Container

Specialty Product Labels for ISBT 128



Blood Product Labels measure 2" x 2" and are supplied 500 labels per roll.
Shown at 70% of actual size. If a scannable proof is needed, please contact Shamrock.


PLATELETS




SB128-50D-13A Part A0




SB128-50D-25A 1st Container




SB128-50D-26A 2nd Container



SB128-50-3




SB128-50-20 1st Container




SB128-50-21 2nd Container




SB128-50-13




SB128-50-25 1st Container




SB128-50-26 2nd Container



SB128-50-9 Open System




SB128-50-12 Open System




SB128-50-34 Bacterial Monitoring


CRYOPRECIPITATED AHF




SB128-60-4



SB128-60-9



SB128-60-1




SB128-60-5 Open System

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




Shamrock Blood Bank Labels for ISBT 128 are printed in compliance with the US Industry Consensus Standard for the Uniform Labeling of Blood and Blood Components Using ISBT 128. They feature adhesive that is compositionally compliant with FDA 21 CFR 175.105.




SB128-60-2



SB128-60-7 Open System



SB128-60-8 Open System



SB128-60-6

SHAMROCK CAN PROVIDE ANY ISBT 128 LABEL

If you do not find the exact ISBT 128 compatible label needed in your Blood Bank, please contact Shamrock for:

- Base Labels • Products Labels • ABO Group Labels • Auxiliary Blood Bag Labels • Consecutively Numbered Labels • Blank Labels for On-Demand Printing •

Shamrock will work with you to craft the precise labels demanded for your procedures. Call 800-323-0249 for immediate assistance.

Specialty Blood Group Labels for ISBT 128

Blood Group Labels are placed in the upper right quadrant of the Base Label. These labels measure 2" x 2" and are supplied 500 labels per roll. Labels shown at 60% of actual size. For complete information on all stock Group Labels, or if a scan-able proof is needed, please contact Shamrock.

ABO/Bh GROUP LABELS




SB128-AP62



SB128-AN06



SB128-BP73




SB128-BN17



SB128-OP51



SB128-ON95



SB128-ABP84



SB128-ABN28

Specialty ISBT 128 Blood Group Labels



Blood Group Labels are placed in the upper right quadrant of the Base Label. These labels are available for all Blood Groups, measure 2" x 2" and are supplied 500 labels per roll. Labels shown at 55% of actual size. If a scannable proof is needed, please contact Shamrock.

POOLED PRODUCT GROUP LABELS

GROUP A POOLED Rh	GROUP B POOLED Rh	GROUP AB POOLED Rh	GROUP O POOLED Rh	POOLED ABO POOLED Rh
SB128-A000	SB128-B000	SB128-C000	SB128-D000	SB128-G000

POOLED PRODUCT GROUP LABELS

Rh POSITIVE	Rh NEGATIVE	A Rh POSITIVE	A Rh NEGATIVE	A Rh POSITIVE
FOR AUTOLOGOUS USE ONLY	FOR AUTOLOGOUS USE ONLY	FOR AUTOLOGOUS USE ONLY	FOR AUTOLOGOUS USE ONLY	BIOHAZARD FOR AUTOLOGOUS USE ONLY
SB128-OP53	SB128-ON97	SB128-AP64	SB128-AN08	SB128-AP65

POOLED PRODUCT GROUP LABELS

A Rh POSITIVE	O Rh POSITIVE	A Rh POSITIVE
FOR DESIGNATED RECIPIENT ONLY	FOR DESIGNATED RECIPIENT ONLY	BIOHAZARD FOR DESIGNATED RECIPIENT ONLY
SB128-AN02	SB128-ON91	SB128-OP49

Contact Shamrock for information on the complete listing of all ABO/Rh Blood Group Labels.

Specialty Auxiliary Labels for ISBT 128 FACILITY LABELS

Printed with your facility information, the Primary processing labels are placed below the Donation ID Label in the upper left quadrant of the Base Label. Secondary processing facility labels are placed on the bottom third portion of the lower right quadrant. Labels are supplied 500 labels per roll. Labels shown at 70% of actual size.

Accurate Blood Center Anywhere, USA FDA Registration Number _____ US License Number _____ SB128-7 5/8" x 2"	Further Processed By Second Blood Center Anywhere, USA FDA Registration Number _____ US License Number _____ SB128-5 5/8" x 2"	Accurate Blood Center Anywhere, USA FDA Registration Number _____ US License Number _____ Properly Identify Intended Recipient See Circular of Information for indications, contraindications, cautions and methods of infusion. This product may transmit infectious agents. Rx Only SB128-6 1 1/8" x 2"	Accurate Blood Center Anywhere, USA FDA Registration Number _____ US License Number _____ Properly Identify Intended Recipient See Circular of Information for indications, contraindications, cautions and methods of infusion. This product may transmit infectious agents. Rx Only VOLUNTEER DONOR SB128-6V 1 1/4" x 2"
--	--	---	--

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



INTENDED RECIPIENT LABELS

Needed for all Autologous or Directed Donation units, labels are placed on the front of the container, immediately above the Donation Identification Number and ABO/Rh Blood Group bar codes. Choose blank labels for on-demand printing or pre-printed labels, which may also be ordered as tie tags. Labels shown at 60% of actual size.

INTENDED RECIPIENT INFORMATION
 Patient Name _____
 ID Number _____
 Date of Birth _____
 Facility _____

INTENDED RECIPIENT INFORMATION
 WB ___ Irrad ___ Patient Name _____
 RBC ___ LKORD ___ ID Number _____
 FFP ___ Other ___ Facility _____
 PLT ___ Birthdate ___/___/___ Collected ___/___/___
 CRYO _____
 Blood Relative: Yes _____

SB128-5001
 1 3/8" x 3"
 500 per roll

SB128-5000
 1 3/8" x 3"
 500 per roll

THERM-466R 1" x 2 1/2"
 3,000 per roll 3" core

DONATION ID NUMBER LABELS

Applied at the time of collection, these consecutively numbered 5/8" x 2" labels cover the top third portion of the upper left quadrant. They are available in rolls or tablets. Each number may be printed as many times as needed. Minimum order 1,000 labels. Labels shown at 60% of actual size.

SB128-SCT
 Tablet Form

SB128-SCR
 Roll Form

COLLECTION DATE & TIME LABELS

For use with bags of Recovered Plasma, these 5/8" x 2" bar coded labels are placed in the middle one third of the upper left quadrant of the Base Label. Labels shown at 80% of actual size.

Collection Date
 016031
31 JAN 2016

SB128-COLDATE
 5/8" x 2"

Collection Date/Time
 0160311200
31 JAN 2016 12:00

SB128-COLDATE-T
 5/8" x 2"

EXPIRATION DATE & TIME LABELS

Measuring 5/8" x 2", these bar coded labels are placed in the top third portion of the lower right quadrant. They feature both bar code and eye readable text. Labels shown at 80% of actual size.

Expiration Date
 0160312359
31 JAN 2016

SB128-DATE
 5/8" x 2"

Expiration Date/Time
 0160311200
31 JAN 2016 12:00

SB128-DATE-T
 5/8" x 2"

SPECIAL TESTING LABELS

Shamrock Special Testing Labels are placed in the lower right quadrant of the Base Label, below the Expiration Date Label. The Emergency Release Label may also be ordered as a tie tag. Call Shamrock at 800-323-0249 for information. Labels shown at 75% of actual size, supplied 500 per roll. If a scan-able proof is needed, please contact Shamrock.

N0106
 Hemoglobin S Negative

SB128-TEST-106
 5/8" x 2"

N0107
 CMV seronegative
 Hemoglobin S Negative

SB128-TEST-107 5/8" x 2"

N0008
 Negative for antibodies to CMV

SB128-TEST-8
 5/8" x 2"

N0009
 Positive for antibodies to CMV

SB128-TEST-9 5/8" x 2"

DONOR UNTESTED

SB128-DU
 5/16" x 1 1/4"


Specialty ISBT 128 DIN Label Sets



PIGGYBACK LABELS PRE-PRINTED/BLANK FOR THERMAL TRANSFER PRINTING

Moisture & Chemical Resistant · Glove Friendly · FDA 21 CFR 175.105 Compliant


Consecutively Numbered DIN sets may be pre-printed with your specified number sequence or supplied blank for thermal transfer printing on-demand. A set of labels, each showing the same number, is “piggybacked” on an adhesive carrier sheet to be placed on the primary blood container. Individual labels are then dispensed as needed. Labels feature standard Thermal Transfer material or, for more demanding procedures, are available in special PolyFlex Transfer and Kimdura Transfer materials. Please call for ordering assistance. Labels shown at 55% of actual size.



4 1/2" x 3 3/4"

Printed: SCBAR-1285 (Thermal Transfer)
 SCBAR-1285-P (PolyFlex Transfer)
 SCBAR-1285-KT (Kimdura Transfer)
 Blank: THERM-PB1711RT (Thermal Transfer)
 THERM-PB1711RP (PolyFlex Transfer)
 THERM-PB1711RKT (Kimdura Transfer)



Bar Codes can be printed at 600 DPI upon request. Call 800-323-0249.



2" x 4"

Printed: SCBAR-1281 (Thermal Transfer)
 SCBAR-1281-P (PolyFlex Transfer)
 SCBAR-1281-KT (Kimdura Transfer)
 Blank: THERM-PB1808RT (Thermal Transfer)
 THERM-PB1808RP (PolyFlex Transfer)
 THERM-PB1808RKT (Kimdura Transfer)


Piggyback Construction



3 1/4" x 4"

Printed: SCBAR-1288 (Thermal Transfer)
 SCBAR-1288-P (PolyFlex Transfer)
 SCBAR-1288-KT (Kimdura Transfer)
 Blank: THERM-PB1816RT (Thermal Transfer)
 THERM-PB1816RP (PolyFlex Transfer)
 THERM-PB1816RKT (Kimdura Transfer)

Color coding is available.



7 3/16" x 4"


Printed: SCBAR-1283 (Thermal Transfer)
 SCBAR-1283-P (PolyFlex Transfer)
 SCBAR-1283-KT (Kimdura Transfer)
 Blank: THERM-PB1815RT (Thermal Transfer)
 THERM-PB1815RP (PolyFlex Transfer)
 THERM-PB1815RKT (Kimdura Transfer)

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




Shamrock Blood Bank Labels for ISBT 128 are printed in compliance with the US Industry Consensus Standard for the Uniform Labeling of Blood and Blood Components Using ISBT 128. They feature adhesive that is compositionally compliant with FDA 21 CFR 175.105.




2⁵/₈" x 4"
Printed: SCBAR-1282 (Thermal Transfer)
SCBAR-1282-P (PolyFlex Transfer)
SCBAR-1282-KT (Kimdura Transfer)
Blank: THERM-PB1787RT (Thermal Transfer)
THERM-PB1787RP (PolyFlex Transfer)
THERM-PB1787RKT (Kimdura Transfer)



3¹/₈" x 4"
Printed: SCBAR-1286 (Thermal Transfer)
SCBAR-1286-P (PolyFlex Transfer)
SCBAR-1286-KT (Kimdura Transfer)
Blank: THERM-PB1814RT (Thermal Transfer)
THERM-PB1814RP (PolyFlex Transfer)
THERM-PB1814RKT (Kimdura Transfer)





6" x 4"
Printed: SCBAR-1284 (Thermal Transfer)
SCBAR-1284-P (PolyFlex Transfer)
SCBAR-1284-KT (Kimdura Transfer)
Blank: THERM-PB1801RT (Thermal Transfer)
THERM-PB1801RP (PolyFlex Transfer)
THERM-PB1801RKT (Kimdura Transfer)



5³/₄" x 3³/₄"
Printed: SCBAR-1287 (Thermal Transfer)
SCBAR-1287-P (PolyFlex Transfer)
SCBAR-1287-KT (Kimdura Transfer)
Blank: THERM-PB1810RT (Thermal Transfer)
THERM-PB1810RP (PolyFlex Transfer)
THERM-PB1810RKT (Kimdura Transfer)

If you do not see the label format you are using, or if you would like to design a completely new format, contact Shamrock. Call 800-323-0249 for complete information and design assistance.




Specialty Blank Cellular Therapy Labels




Shamrock supplies labels designed to meet your cellular therapy needs. The Blank Thermal Labels illustrated below are sized to ISBT 128 standards and available for on-site printing as needed. The printed samples are for illustration only. For ISBT 128 labels pre-printed to your specifications, please contact Shamrock.

BASE LABELS




SB128-44 800/roll
SB128-44-PL (perfs thru label & liner) 1,400/roll
4" x 4" Blank



SB128-22 2" x 2" Blank 1,600/roll



SB128-23
SB128-23-PL (perf thru label & liner)
2" x 3" Blank 1,600/roll




SB128-32
SB128-32-PL (perf thru label & liner)
3" x 2" Blank 1,100/roll

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




INTENDED RECIPIENT LABELS



SB128-15375-PL (perfs thru label & liner)
1 1/2" x 3 3/4" Blank 2,100/roll

CRYO VIAL LABELS



SB128-7515 3/4" x 1 1/2"
Blank 3,900/roll

PRE-PRINTED MESSAGE LABEL




SB128-CTSM-1 2" x 2"


Specialty Blank Tissue Bank Labels

Shamrock blank thermal labels in ISBT 128 format for tissue banks are also available for on-site printing or pre-printed to your specifications. Contact Shamrock for complete information on the variety of tissue bank labels available.


CONTAINER LABELS



SB128-44 4" x 4" Blank 800/roll



THERM-487R 1/2" x 1 1/2"
Blank 5,000/roll



SB128-22 2" x 2"
Blank 1,600/roll

Blank labels are stock. Contact Shamrock for information on custom pre-printed labels.

Specialty Codabar Blood Bank Labels

SHAMROCK BLOOD BANK LABELS WITH CODABAR SYMBOLOGY ARE STILL AVAILABLE

Codabar labels feature top quality materials including glove free adhesive that is **COMPOSITIONALLY COMPLIANT WITH FDA 21 CFR 175.105**. Shamrock Codabar labels stick tightly throughout standard procedures and storage. Blood Product Labels measure 1" x 2¼" and are supplied 500 labels per roll. Shown at 70% of actual size. If a scannable proof is needed, please contact Shamrock.

CODABAR PRODUCT LABELS

<p>WHOLE BLOOD</p>	<p>CPDA-1 RED BLOOD CELLS From 500 mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>04060</p> 	<p>RED BLOOD CELLS LEUKOCYTES REDUCED WASHED From 450 mL CPD Whole Blood. Store at 1 to 6 C.</p> <p>06000</p> 	<p>FRESH FROZEN PLASMA ____ mL from 450 mL CPD Whole Blood. Store at -18 C or colder.</p> <p>18201</p> 
<p>CPDA-1 WHOLE BLOOD Approx 450 mL plus 63 mL CPDA-1. Store at 1 to 6 C.</p> <p>00160</p> 	<p>CPDA-1 RED BLOOD CELLS LOW VOLUME Approx ____ mL from ____ mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>04061</p> 	<p>CPDA-1 RED BLOOD CELLS DIVIDED Part ____ of ____ Approx ____ mL from 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>34061</p> 	<p>FRESH FROZEN PLASMA ____ mL from 500 mL CPD Whole Blood. Store at -18 C or colder.</p> <p>18201</p> 
<p>CPDA-1 WHOLE BLOOD IRRADIATED Approx 450 mL plus 63 mL CPDA-1. Store at 1 to 6 C.</p> <p>00260</p> 	<p>CPDA-1 RED BLOOD CELLS LEUKOCYTES REDUCED From 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>04360</p> 	<p>AS-1 RED BLOOD CELLS ADENINE - SALINE ADDED DIVIDED 15.4 mEq Sodium added. Part ____ of ____ Approx ____ mL from 450 mL CPD Whole Blood. Store at 1 to 6 C.</p> <p>34211</p> 	<p>FRESH FROZEN PLASMA ____ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C or colder.</p> <p>18201</p> 
<p>CPDA-1 WHOLE BLOOD LOW VOLUME Approx ____ mL plus 63 mL CPDA-1. Store at 1 to 6 C.</p> <p>00161</p> 	<p>AS-1 RED BLOOD CELLS ADENINE - SALINE ADDED LEUKOCYTES REDUCED 15.4 mEq Sodium added. From 450 mL CPD Whole Blood. Store at 1 to 6 C.</p> <p>04710</p> 	<p>AS-3 RED BLOOD CELLS ADENINE - SALINE ADDED DIVIDED 15.0 mEq Sodium added Part ____ of ____ Approx ____ mL from 450 mL CP2D Whole Blood. Store at 1 to 6 C.</p> <p>34231</p> 	<p>FRESH FROZEN PLASMA ____ mL from 450 mL CP2D Whole Blood. Store at -18 C or colder.</p> <p>18201</p> 
<p>AUTOLOGOUS CPDA-1 WHOLE BLOOD Approx 450 mL plus 63 mL CPDA-1. Store at 1 to 6 C.</p> <p>30160</p> 	<p>AS-1 RED BLOOD CELLS ADENINE - SALINE ADDED LEUKOCYTES REDUCED 16.9 mEq Sodium added. From 500 mL CPD Whole Blood. Store at 1 to 6 C.</p> <p>04710</p> 	<p>CPDA-1 RED BLOOD CELLS LEUKOCYTES REDUCED DIVIDED Part ____ of ____ from 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>34361</p> 	<p>PLASMA ____ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C or colder.</p> <p>18401</p> 
<p>RED BLOOD CELLS</p>	<p>AS-3 RED BLOOD CELLS ADENINE - SALINE ADDED LEUKOCYTES REDUCED 15.0 mEq Sodium added. From 450 mL CP2D Whole Blood. Store at 1 to 6 C.</p> <p>04730</p> 	<p>PLASMA</p>	<p>PLASMA ____ mL from 450 mL CPD Whole Blood. Store at -18 C or colder.</p> <p>18401</p> 
<p>CPDA-1 RED BLOOD CELLS From 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>04060</p> 	<p>RED BLOOD CELLS LEUKOCYTES REDUCED WASHED From 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>06000</p> 	<p>PLASMA FROZEN WITHIN 24 HOURS AFTER PHLEBOTOMY ____ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C or colder.</p> <p>18101</p> 	<p>THAWED PLASMA ____ mL from 450 mL CPD/CPDA-1 Whole Blood. Store at 1 to 6 C.</p> <p>29381</p> 

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Blood Product Labels measure 1" x 2 1/4". Shown at 70% of actual size. If a scannable proof is needed, please contact Shamrock.

THAWED FRESH FROZEN PLASMA
 ___ mL from 450 mL CPD Whole Blood.
 Store at 1 to 6 C.
 18900

SB-649

PLATELETS POOLED
 Approx ___ mL from ___ donors. Collected in CPD. Store at 20 to 24 C.
 12091

SB-661

CRYOPRECIPITATED AHF POOLED
 Contains Cryoprecipitate from ___ individual donors. From 450 mL CPD Whole Blood. Store at 20 to 24 C.
 10191

SB-673

AS-3 RED BLOOD CELLS
 ADENINE-SALINE ADDED LEUKOCYTES REDUCED IRRADIATED
 15.0 mEq Sodium added. From 450 mL CPD2 Whole Blood. Store at 1 to 6 C.
 05730

SB-944

FRESH FROZEN PLASMA DIVIDED
 Part ___ of ___ mL from 450 mL CPD Whole Blood. Store at -18 C or colder.
 48201

SB-650

PLATELETS PHERESIS LEUKOCYTES REDUCED
 ___ mL containing approx ___ mL of ACD-A Anticoagulant. Store at 20 to 24 C.
 12710

SB-664

IRRADIATED

AS-5 RED BLOOD CELLS
 ADENINE-SALINE ADDED LEUKOCYTES REDUCED IRRADIATED
 15.0 mEq Sodium added. From 450 mL CPD Whole Blood. Store at 1 to 6 C.
 05750

SB-945

FRESH FROZEN PLASMA DIVIDED
 Part ___ of ___ mL from 450 mL CPDA-1 Whole Blood. Store at -18 C or colder.
 48201

SB-651

PLATELETS PHERESIS LEUKOCYTES REDUCED
 ___ mL containing approx ___ mL of Anticoagulant. Store at 20 to 24 C.
 12710

SB-665

CPDA-1 RED BLOOD CELLS IRRADIATED
 From 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.
 05060

SB-937

RED BLOOD CELLS LEUKOCYTES REDUCED BY WASHING IRRADIATED
 From 450 mL Whole Blood. Store at 1 to 6 C.
 05800

SB-946

PLATELETS

CRYOPRECIPITATED

AS-1 RED BLOOD CELLS
 ADENINE-SALINE ADDED IRRADIATED
 15.4 mEq Sodium added. From 450 mL CPD Whole Blood. Store at 1 to 6 C.
 05210

SB-939

RED BLOOD CELLS DEGLYCEROLIZED IRRADIATED
 From 450 mL Whole Blood. Store at 1 to 6 C.
 07400

SB-949

PLATELETS
 Approx 45-65 mL from 450 mL CPD Whole Blood. Store at 20 to 24 C.
 12000

SB-908

CRYOPRECIPITATED AHF
 From 450 mL CPD Whole Blood. Store at -18 C or colder.
 10100

SB-905

CPD RED BLOOD CELLS LEUKOCYTES REDUCED IRRADIATED
 From 450 mL CPD Whole Blood. Store at 1 to 6 C.
 05350

SB-941

PLATELETS PHERESIS IRRADIATED
 ___ mL containing approx ___ mL of Anticoagulant. Store at 20 to 24 C.
 12610

SB-954

PLATELETS
 Approx 45-65 mL from 450 mL CPDA-1 Whole Blood. Store at 20 to 24 C.
 12000

SB-907

CRYOPRECIPITATED AHF POOLED
 No. Pooled Units ___ Volume ___ mL. From 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.
 10191

SB-670

CPDA-1 RED BLOOD CELLS LEUKOCYTES REDUCED IRRADIATED
 From 450 mL CPDA-1 Whole Blood. Store at 1 to 6 C.
 05360

SB-942

PLATELETS PHERESIS LEUKOCYTES REDUCED IRRADIATED
 ___ mL containing approx ___ mL of Anticoagulant. Store at 20 to 24 C.
 12810

SB-667

PLATELETS PHERESIS
 Platelet count ___ x 10⁵/µL ___ mL containing approx ___ mL of Anticoagulant. Store at 20 to 24 C.
 12011

SB-657

CRYOPRECIPITATED AHF POOLED
 No. Pooled Units ___ Volume ___ mL. From 450 mL units of 0.9% Sodium Chloride added. Store at 20 to 24 C.
 10191

SB-672

AS-1 RED BLOOD CELLS
 ADENINE-SALINE ADDED LEUKOCYTES REDUCED IRRADIATED
 15.4 mEq Sodium added. From 450 mL CPD Whole Blood. Store at 1 to 6 C.
 05710

SB-943

Contact Shamrock for more information on Codabar Labels.

Specialty Codabar Blood Bank Labels

Shamrock supplies labels for every Codabar based application. Samples of various labels are shown below. Call for information on the complete line of stock Codabar labels. Labels shown at 80% of actual size. If a scannable proof is needed please contact Shamrock.

BASE LABELS

Collection Date	Unit Number	EXPIRES
Place Product Label here after bag is filled.		Affix grouping label here after all required testing has been completed.
See circular of information for indications, contraindications, cautions and methods of infusion. VOLUNTEER DONOR This product may transmit infectious agents. R _x Only		Affix Collection/Processing I.D. Label Here
PROPERLY IDENTIFY INTENDED RECIPIENT		

SBB-703 3 1/4" x 4"
SBB-723 4" x 4"


SATELLITE LABEL

Collection Date	Unit Number	EXPIRES
Place Product Label Here		Place Grouping Label Here
See circular of information for indications, contraindications, cautions and methods of infusion. VOLUNTEER DONOR This product may transmit infectious agents. R _x Only		
PROPERLY IDENTIFY INTENDED RECIPIENT		

SBB-740 3" x 4"

FACILITY LABELS


YOUR HOSPITAL BLOOD BANK
Your City, State 01234
Phone (123) 456-7890
License #123



FDA Registration #1234567

SB-5 1 1/8" x 1 3/4"
500 labels per roll

Collected and Processed by
YOUR HOSPITAL BLOOD BANK
Your City, State 01234
Phone (123) 456-7890
License #123



FDA Registration #1234567

SB-6 1 1/8" x 1 3/4"
500 labels per roll

DONOR OPTION LABELS




USE MY BLOOD




DO NOT USE MY BLOOD

SBD-3 2 1/4" x 2 1/2" 500 per roll

MESSAGE LABELS

BIOHAZARD

REASON _____




ABO-50 1 1/2" x 1 3/4"
500 labels per roll

BIOHAZARD




ABO-45 1 7/8" x 1 3/4"
500 labels per roll

HOLD FOR FURTHER PROCESSING

ABO _____


Rh _____

STS _____

HBs _____

Ag _____

OTHER _____



ABO-44 1 1/2" x 1 3/4"
500 labels per roll

FOR AUTOLOGOUS USE ONLY
Reserved for _____

NAME _____

ABO/Rh _____

COLLECTION ON _____

HOSPITAL _____


PATIENT ID _____

ABO-40 1 1/2" x 1 3/4"
500 labels per roll

NOT FOR TRANSFUSION

REASON _____

Date _____ Initial _____



ABO-41 1 1/2" x 1 3/4"
500 labels per roll

GROUP TYPE LABELS

A

Rh POSITIVE



ABO-10 1 1/2" x 1 3/4"
500 labels per roll

O

Rh POSITIVE



ABO-14 1 1/2" x 1 3/4"
500 labels per roll

EXPIRES _____

A

Rh POSITIVE



ABO-20 1 3/4" x 1 3/4"
500 labels per roll

EXPIRES _____

O

Rh POSITIVE



ABO-24 1 3/4" x 1 3/4"
500 labels per roll

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




Specialty Consecutive Number Labels

BAR CODED NUMBER LABEL SETS

- Roll Form or Tablets • Single-Ply or Piggyback* Material
- You Specify Number Sequence, Prefix, Suffix, Bar Code Symbology.
- One or Two Line Custom Imprint on Bar Code Labels if Desired.

Custom Formats and Special Label Materials Available

*Piggyback Format



SCBAR-21 Overall: 1 1/2" x 1 1/4"




SCBAR-08000 Overall: 3/4" x 2"
2 Eye/Bar: 3/4" x 1 1/2"



SCBAR-20 Overall: 1 1/4" x 1 7/8"

Bar Codes can be printed at 600 DPI. Call 800-323-0249



SCBAR-1 Labels: 7/16" x 1 7/8"
Single-Ply, Roll Form Only with No Custom Imprint



SCBAR-23 Overall: 2 1/2" x 4 5/8"



SCBAR-22 Overall: 3 3/4" x 1 1/4"




SCBAR-01000 Overall: 2 1/2" x 2"




SCBAR-01042 Overall: 2 1/2" x 2"

2 Line Custom Imprint



SCBAR-09000 3" x 2"



SCBAR-09021 3" x 2"

1 Line Custom Imprint



SCBAR-07000 Overall: 4 1/2" x 2"

NOTE: Label width is dependent on length of bar code and may vary from dimension indicated under illustration.

Specialty Consecutive Number Labels

PIGGYBACK CONSECUTIVE NUMBER LABELS


Specialty Piggyback Label stock provides a set of individual labels with the same number printed as many times as desired, all held in place on a sheet of adhesive backed carrier paper. The pressure-sensitive carrier holds tightly to a central item; the individual number labels may then be used one at a time.



Standard Specialty Piggyback Consecutive Number Labels measure 1/2" x 1 1/4". They are supplied printed black on white label stock with "No" prefix.

Piggyback Material

SCPB-30 (Piggyback labels supplied in tablets)
SCPB-31 (Piggyback labels supplied in rolls)



SPECIALTY CUSTOM NUMBERING SYSTEMS

The labels pictured illustrate Shamrock's versatility in producing specially designed labels sequentially numbered or bar coded. Send us a sketch of the label you are using or let us work with you to design a label to fit your special needs.



	R 126	
	R 126	MEDICAL CENTER
R 126	TECH	4885 4885 4885 4885
R 126	PATIENT	4885 4885 4885 4885
R 126	CHART RECORD	4885 4885 4885 4885
R 126	1st UNIT	4885 4885 4885 4885
R 126	2nd UNIT	4885 4885 4885 4885
R 126	3rd UNIT	4885 4885 4885 4885
R 126	4th UNIT	4885 4885 4885 4885
R 126	5th UNIT	4885 4885 4885 4885
R 126	6th UNIT	4885 4885 4885 4885
R 126	7th UNIT	4885 4885 4885 4885
R 126	8th UNIT	4885 4885 4885 4885
R 126	9th UNIT	4885 4885 4885 4885
R 126	10th UNIT	4885 4885 4885 4885

B 0546	B 0546	B 0546	B 0546
SPECIMEN TUBE	CHART RECORD	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546
ST. MARY'S HOSPITAL Your Town, USA	B 0546	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546
B 0546	B 0546	B 0546	B 0546

T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000
T00000	T00000	T00000	T00000

Name	SSN#	Time	Init.
MR#			
Date			

1. Fill out card, 2. Push & slide into recipient's hand, 3. Snap off strip at dotted line.

TAPE AND PRE-CUT NUMBER LABELS


Numbers on Shamrock 1/2" wide Labeling Tape may be printed up to 16 times. Pre-Cut labels, 1/2" x 1 1/4", may be printed as many times as needed. Both are available in your choice of 19 colors. Call Shamrock for complete information.

Tape Labels



SCR-10 (500 labels/roll)
SCR-20 (2,000 labels/roll)

Pre-Cut Labels



SCS-12 (Tablet)
SCR-14 (Roll)

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM





SpecialtyTags for Blood Banking

Shamrock Blood Tags feature sturdy card stock that accepts handwritten information in pen or marker. Stock Tags, supplied in packages of 100, are illustrated. Tags are shown at 40% of actual size.


Blood Tags Customized to Fit Your Procedures

Shamrock specializes in Tags produced to your exact specifications. Order tags imprinted to meet the requirements of your blood bank procedures and/or consecutively numbered tags for greater accountability.

Contact Shamrock at 800-323-0249 for complete information.



SBT-1002 6 1/4" x 3 1/8"



SBT-1000 6 1/4" x 3 1/8"

Specialty Award Stickers

These vividly printed labels measure 1 3/4" in diameter, are backed with removable adhesive and feature perforations between each label for easy dispensing. Award Stickers are supplied 200 per roll. Call Shamrock at 800-323-0249 for complete information on stickers personalized to your specifications.




Award-31




Award-44



Award-45




Award-40



Award-17

Specialty Labels and Tapes of your own design

Shamrock can supply labels with wording you specify in tape or pre-cut label form. Labeling Tape offers a choice of 31 colors, see page BB-36. Pre-cut labels, available in the same 31 colors, may be ordered in special materials and adhesives. Contact Shamrock for design assistance and ordering information.



Specialty Blood Bank Communication Labels

Shamrock supplies labels for every Codabar based application. Samples of various labels are shown below. Call for information on the complete line of stock Codabar labels. Labels shown at 80% of actual size. If a scannable proof is needed please contact Shamrock.

Message Labels

UPCR-1000 series labels measure 5/16" x 1 1/4".

Store at Room Temperature UPCR-1091	CALL REPORT TO DOCTOR UPCR-1005 Fl. Red	ISOLATION UPCR-1248 Fl. Red	RHO GAM UPCR-1305 Fl. Orange	ASAP UPCR-1077	DO NOT REFRIGERATE UPCR-1012 Fl. Green	NAME ALERT UPCR-1020 Fl. Red
STAT UPCR-1044 Fl. Red	BLOOD BANK UPCR-1175 Fl. Red	EXTRA TUBE UPCR-1332 Fl. Chartreuse	FROZEN UPCR-1083 Fl. Red	SURGERY AM _____ PM _____ UPCR-1046 Fl. Pink	BLOOD/BODY FLUID PRECAUTIONS UPCR-1004 Fl. Pink	HEPATITIS UPCR-1244 Fl. Red
PRE-OP UPCR-1030 Fl. Red	BLOOD RELEASED UPCR-1331 Fl. Green	Hold for Possible Additional Testing UPCR-1066 Fl. Pink	BLOOD DRAWN UPCR-1055 Fl. Chartreuse	QUARANTINE UPCR-1334 Fl. Chartreuse	SICKLE CELL NEGATIVE UPCR-1351 Fl. Pink	STAT UPCR-3005 3/4" x 1 1/2" Fl. Red

PATIENT REFUSES BLOOD OR BLOOD PRODUCTS
 UPCR-6071 1" x 2 1/4" Fl. Red

QUARANTINE
 UPCR-6072 1" x 2 1/4" Fl. Red

Name Alert
 UPCR-6070 1" x 2 1/4" Fl. Red

THE BLOOD PREPARED FOR THIS PATIENT HAS NOW BEEN RELEASED
 Date _____ Time _____
 No. of units _____
 UPCR-7009 1" x 3" Fl. Chartreuse

Rh₀ POSITIVE
 UPCR-6084 1" x 2 1/4" Fl. Green

Rh NEGATIVE
 UPCR-6083 1" x 2 1/4" Fl. Red

LOW VOLUME UNIT
 _____ mL Red Blood Cells
 UPCR-7031 1" x 3"

LEUKOCYTES REDUCED
 UPCR-7030 1" x 3"

PLEASE RETURN TO BLOOD BANK WHEN EMPTY
 UPCR-7032 1" x 3" Fl. Orange

Patient Identification Labels

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Shamrock Communication Labels feature glove free adhesive that is COMPOSITIONALLY COMPLIANT with FDA 21 CFR 175.105.



Stock Communication Labels are supplied 500 per roll unless otherwise indicated.

Blood Group Labels

BLOOD GROUP
A

ABO-70 1" x 1"

BLOOD GROUP
B

ABO-72 1" x 1"

BLOOD GROUP
O

ABO-74 1" x 1"

BLOOD GROUP
AB

ABO-76 1" x 1"

BLOOD GROUP AND Rh TYPE
RECHECKED BY
TECH. _____
HOSP. _____
DATE _____

SB-70 1" x 1 3/4"
275 labels/roll

Blood Group and Rh Type
Rechecked By
Your Blood Bank Name Here
Tech. _____
Date _____

SB-70C 1" x 1 3/4" (custom imprint)
275 labels/roll

Product and Component Labels

PLASMA

UPCR-1070
5/16" x 1 1/4"

WHOLE BLOOD

UPCR-1352 5/16" x 1 1/4"

Fresh Frozen Plasma
Thawed: _____
by _____ date _____ time _____
Expires: date _____ time _____

UPCR-5015 3/4" x 2"

FRESH FROZEN PLASMA

Amt. _____
Expires _____

UPCR-6076 1" x 2 1/4"

FRESH FROZEN PLASMA
Date Thawed _____
Time Thawed _____
Expiration Date _____
Expiration Time _____

UPCR-6077 1" x 2 1/4"

POOLED PLATELETS

Pool # _____ Tech _____
Date Pooled _____ Time Pooled _____
Infuse Within _____ Hours After Pooling
Volume _____ # of Units _____
ABO & Rh _____

UPCR-6078
1" x 2 1/4"

RED BLOOD CELLS

Expiration Date _____
Time _____

UPCR-6079 1" x 2 1/4" Fl. Red

PRODUCT LABEL MATCHES LABEL AFFIXED TO BOX

INITIALS _____ DATE _____

SB-62 1" x 2 1/2"
1,000 labels/roll

FRESH FROZEN PLASMA

PLASMA THAWED: DATE _____ HR _____ (AM)(PM)
EXPIRATION DATE _____ HR _____ (AM)(PM)

UPCR-7012 1" x 3" Fl. Pink

ALIQUOT UNIT

DONOR # _____ Vol _____ cc
PRODUCT _____ TECH _____
EXPIRATION DATE _____
TIME _____
TECH: _____

UPCR-9407 2" x 3"

POOLED BLOOD COMPONENT

Blood Component Pooled _____ Tech _____
Date and Time Thawed (for Cryoprecipitate only) _____
Pooled Date and Time _____ Expiration Date and Time _____

UNIT NUMBER	ABO GROUP AND Rh TYPE	EXPIRATION DATE

UPCR-9705 2 1/2" x 4"

POOLED PLATELETS

Volume _____ Date Pooled _____ Time Pooled _____
Must be used within _____ Hrs. Donor Numbers _____

UPCR-9502 2" x 4"

PEDIATRIC/NEONATAL COMPONENT LABEL

Product _____ Amt. _____
Patient _____
ID # _____ ABO, Rh _____
Donor # _____ ABO, Rh _____
Expir. Date _____ Time _____
Tech _____ Date _____

UPCR-941 2" x 3" Fl. Orange

Specialty Blood Bank Communication Labels

Irradiation Labels

IRRADIATE
ALL PRODUCTS

UPCR-1335 5/16" x 1 1/4"
Fl. Pink

IRRADIATED

UPCR-1336 5/16" x 1 1/4"

TO BE IRRADIATED

UPCR-1337 5/16" x 1 1/4"
Fl. Pink

DO NOT
IRRADIATE

UPCR-1350 5/16" x 1 1/4"

IRRADIATED

Date _____
By _____

UPCR-6082 1" x 2 1/4"

DO NOT
IRRADIATE

UPCR-6098 1" x 2 1/4"

IRRADIATED
BLOOD PRODUCT

_____ Rads _____ Date

UPCR-6080 1" x 2 1/4"

IRRADIATED

UPCR-6081 1" x 2 1/4"

Compatibility Labels

COMPATIBLE

UPCR-1338 5/16" x 1 1/4"
Fl. Pink

DOUBLE
X-MATCHED

UPCR-1339 5/16" x 1 1/4"

TYPE & SCREEN

UPCR-1340 5/16" x 1 1/4"

PATIENT HAS ATYPICAL
ANTIBODIES. ALLOW
EXTRA TIME FOR FURTHER
CROSSMATCHES.

UPCR-6086 1" x 2 1/4" Fl. Red

UNCROSSMATCHED
BLOOD

UPCR-6088 1" x 2 1/4" Fl. Red

Name _____
CROSSMATCH INCOMPLETE

UPCR-6085 1" x 2 1/4" Fl. Pink

THIS UNIT IS:
 NOT CROSSMATCHED
 PARTIALLY X-MATCHED
Keep patient under observation during
transfusion. Blood Bank will notify Doctor
of any incompatibility as X-Matched.

UPCR-6087 1" x 2 1/4"

Patient's Name _____
_____ units of blood have been
typed and crossmatched. This blood will be
released _____, unless doctor
indicates otherwise.

UPCR-7014 1" x 3" Fl. Red

COMPATIBILITY LABEL

Pt. Name _____
Hosp. # _____ Rm. # _____
Pt. Group _____ Rh _____
Donor # _____
Donor Group _____ Rh _____
Crossmatch Compatible _____
Date _____ Tech. _____
Dr. _____

UPCR-9408 2" x 3"

Patient's Blood Contains
UNEXPECTED ANTIBODIES.
ADDITIONAL UNITS OF COMPATIBLE BLOOD
MAY NOT
BE AVAILABLE IN AN EMERGENCY.

UPCR-7039 1" x 3" Fl. Chartreuse

This donor unit No: _____ has been processed for:
PATIENT: _____
HOSPITAL NO: _____

LABORATORY STUDIES

PATIENT: ABO Group _____ Type: _____
DONOR: ABO Group _____ Type: _____

INTERPRETATION OF COMPATIBILITY TESTS:

Compatible
 Emergency release-Uncrossmatched

TECH: _____ **Date** _____

UPCR-9409 2" x 3"

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Stock Communication Labels are supplied 500 per roll unless otherwise indicated.

Autologous Labels

AUTOLOGOUS BLOOD

UPCR-1302 5/16" x 1 1/4"

ALERT

Autologous and/or designated donor units were donated for this patient. Please check availability.

UPCR-6089 1" x 2 1/4" Fl. Pink

AUTOLOGOUS BLOOD

FOR _____

UPCR-6090 1" x 2 1/4" Fl. Green

AUTOLOGOUS DONOR

UPCR-6091 1" x 2 1/4"

FOR AUTOLOGOUS USE ONLY

UPCR-6092 1" x 2 1/4"

FOR AUTOLOGOUS USE ONLY

Pt. Name _____

Pt. Sign _____

I.D. No. _____ ABO & Rh _____

Date _____ Tech. _____

UPCR-6094 1" x 2 1/4" Fl. Green

For Autologous Use Only CAUTION

Biohazardous Material Autoclave Prior to Disposal

UPCR-6093 1" x 2 1/4" Fl. Red

ATTENTION

THIS PATIENT HAS _____ AUTOLOGOUS UNITS
_____ DIRECTED UNITS
CHECK TYPE AND X-MATCH BEFORE GIVING. AVAILABLE IN BLOOD BANK

UPCR-7015 1" x 3" Fl. Green

For Autologous Use Only

Tests for Infectious Diseases and Unexpected Antibodies Not Done

UPCR-7036 1" x 3" Fl. Green

NOT HOMOLOGOUS BLOOD FOR AUTOLOGOUS TRANSFUSION ONLY

UPCR-7037 1" x 3"

Directed Donor Labels

DIRECTED DONATION

UPCR-1304 5/16" x 1 1/4" Fl. Green

DIRECTED DONATION

FOR _____

UPCR-6035 1" x 2 1/4" Fl. Red

DIRECTED DONOR

HOLD FOR _____
UNTIL (DATE) _____

UPCR-6095 1" x 2 1/4" Fl. Red

Cord Blood Labels

CORD BLOOD

UPCR-1303 5/16" x 1 1/4" Fl. Pink

DELIVERY ROOM CORD BLOOD

NAME _____

HOSP. # _____

DATE _____

UPCR-6097 1" x 2 1/4" Fl. Green

CORD BLOOD

Mother's Name: _____

Sex of Baby: Male _____ Female _____

Doctor's name: _____

Date: _____ Time: _____

UPCR-6096 1" x 2 1/4"

CORD BLOOD

UPCR-7010 1" x 3"

CORD BLOOD

Baby's Name _____ Sex M F

Baby's ID # _____ Birthdate _____ Time _____

Mother _____ Type _____ Nurse Initials _____

UPCR-7011 1" x 3"

CORD BLOOD

I.D. NO. _____ DATE ____ / ____ / ____

PATIENT NAME _____

PEDIATRICIAN _____

UPCR-7038 1" x 3" Fl. Pink

CORD BLOOD

Baby # _____

Mother # _____ Name _____

Mother Group _____ Rh _____

Rm _____ Dr. _____

Delivery Date _____ Time _____

UPCR-9410 2" x 3"

Specialty Blood Bank Communication Labels

Quality Control Labels & Tapes

DO NOT REFRIGERATE

UPCR-1012 5/16" x 1 1/4"
Fl. Green

REFRIGERATE

UPCR-1034 5/16" x 1 1/4"
Fl. Chart.

**SHORT DATED
USE FIRST**

UPCR-1090 5/16" x 1 1/4"
Fl. Red

USE FIRST

UPCR-1054 5/16" x 1 1/4"
Fl. Chart.

USE FIRST

UPCR-3022 3/4" x 1 1/2"
Fl. Chart.

**USE THIS
LOT FIRST**

UPCR-3023 3/4" x 1 1/2"
Fl. Chart.

**BLOOD OUTDATES
USE FIRST**

UPCR-3034 3/4" x 1 1/2"
Fl. Pink

DATE RECEIVED

UPCR-3035 3/4" x 1 1/2"
Fl. Pink

**NEW LOT
DO NOT USE**

UPCR-3036 3/4" x 1 1/2"
Fl. Orange

DO NOT INFUSE AFTER

AM
PM

TIME _____ DATE _____

UPCR-6100 1" x 2 1/4"

DO NOT REFRIGERATE

UPCR-6101-T 1" x 500" Fl. Green
(Tape)

DO NOT REFRIGERATE

UPCR-6101 1" x 2 1/4" Fl. Green

DO NOT REFRIGERATE

UPCR-6102 1" x 2 1/4" Fl. Orange

REFRIGERATE

UPCR-6103-T 1" x 500" Fl. Chart.
(Tape)

REFRIGERATE

UPCR-6103 1" x 2 1/4" Fl. Chart.

**STORE AT ROOM
TEMPERATURE**

UPCR-6104-T 1" x 500" (Tape)

**STORE AT ROOM
TEMPERATURE**

UPCR-6104 1" x 2 1/4"

Date Prep. _____
Exp. Date _____
Made By _____
Lot No. _____

UPCR-3013 3/4" x 1 1/2"
Fl. Green

Date Rec'd _____
Date Opened _____
Exp. Date _____

UPCR-3015 3/4" x 1 1/2"

Date Rec'd _____
Date Opened _____
Exp. Date _____
Reag. Exp. Control _____

UPCR-3016 3/4" x 1 1/2"

Date Rec'd _____ Initial _____
Date Checked _____ Initial _____
Date in Use _____
Expires _____

UPCR-6056 1" x 2 1/4"
Fl. Orange

Date Rec'd _____
Date Opened _____
Exp. Date _____
Room Temp. Refrig.

UPCR-6058 1" x 2 1/4"

LABORATORY REAGENT

PREP. DATE _____ CHECKED _____

EXPIRATION DATE _____

REAGENT _____

CONC. _____

PROCEDURE _____

MADE BY _____

Reagent _____

Prep. Date _____

Exp. Date _____

Storage _____

Room Refrig.

Initials _____

UPCR-9500 2" x 4"

UPCR-31 1 1/2" x 3"

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Shamrock Communication Labels feature glove free adhesive that is COMPOSITIONALLY COMPLIANT with FDA 21 CFR 175.105.



Stock Communication Labels are supplied 500 per roll unless otherwise indicated.

Antigen Testing Labels

ANTIBODY SCREENING POSITIVE

UPCR-1342 5/16" x 1 1/4"
Fl. Red

CMV NEGATIVE

UPCR-1344 5/16" x 1 1/4"
Fl. Green

CAUTION

Tests For Infectious Diseases Not Done

UPCR-3039 3/4" x 1 1/2"
Fl. Red

Cold Agglutinin Use a Blood Warmer

UPCR-3040 3/4" x 1 1/2"

Cold Agglutinin Warm Before Infusion

UPCR-3041 3/4" x 1 1/2"
Fl. Green

ATYPICAL ANTIBODIES

UPCR-1343 5/16" x 1 1/4"
Fl. Red

Negative Antibody Screen

UPCR-1345 5/16" x 1 1/4"
Fl. Green

ABO & Rh CONFIRMED

Date/Tech.

UPCR-5004 3/4" x 2" Fl. Pink

TYPE & Rh RETESTED

Date/Tech.

UPCR-5005 3/4" x 2" Fl. Orange

TYPE & Rh CONFIRMED

Date/Tech.

UPCR-5016 3/4" x 2" Fl. Orange

ABO GROUP CONFIRMED Rh OF NEGATIVE UNITS CONFIRMED

DATE/TECH

UPCR-6106 1" x 2 1/4"

EMERGENCY BLOOD RELEASE

THE FOLLOWING TESTS HAVE NOT BEEN COMPLETED.

X-match HTLV-I HIV HBcAb
 HBsAg Ab screen ALT _____

UPCR-6107 1" x 2 1/4"

Negative for tests for antibodies to HIV-1/2, HB_s, HCV, Chagas and HTLV-I/II and nonreactive for HB_s Ag, STS, WNV, HCV RNA, HIV-1 RNA and HBV DNA. This product may transmit infectious agents.

UPCR-6099 1" x 2 1/4"

SPECIAL ANTIGEN TYPING

Ag. type	Date	Tech.

UPCR-6110 1" x 2 1/4" Fl. Pink

UNIT NEGATIVE

When tested for: _____ Antigen.

Date: _____ Tech: _____

UPCR-6111 1" x 2 1/4"

UNIT POSITIVE

When tested for: _____ Antigen.

Date: _____ Tech: _____

UPCR-6112 1" x 2 1/4" Fl. Red

This donor unit has been antigen tested and found positive or negative for:

C _____	K _____	Jk ^a _____	A _i _____
c̄ _____	k _____	Jk ^b _____	
E _____	Fy ^a _____	S _____	
ē _____	Fy ^b _____	s̄ _____	

Incompatible with:
Date: _____ Tech: _____

UPCR-9111 1 1/2" x 3"

HB_sAg, HB_cAb, STS
ALT, Anti-HIV & HTLVI
Testing has not been done

UPCR-6108 1" x 2 1/4"

Filter Labels

FILTER

UPCR-1346 5/16" x 1 1/4"
Fl. Pink

PRE-FILTERED

UPCR-1347 5/16" x 1 1/4"
Fl. Pink

FILTERED BLOOD

Expiration
Date _____ Time _____

UPCR-6114 1" x 2 1/4"

PRE-FILTERED LEUKOCYTE-REDUCED BLOOD PRODUCT

UPCR-6116 1" x 2 1/4" Fl. Green

Use Leukocyte Reduction Filter Prior to Transfusion

UPCR-6118 1" x 2 1/4" Fl. Orange

ATTACH MICRO-PORE FILTER PRIOR TO USE

UPCR-6113 1" x 2 1/4" Fl. Pink

PRE-FILTERED LEUKOCYTE-POOR BLOOD PRODUCT

UPCR-6115 1" x 2 1/4" Fl. Orange

TRANSFUSE THROUGH _____ FILTER
(BRAND NAME)
FOR LEUKOCYTE-POOR RED CELLS

UPCR-6117 1" x 2 1/4" Fl. Orange

Specialty Blood Bank Communication Labels

Biohazard Warning Labels & Tapes



SBH-128 5/8" Dia.
Fl. Red



UPCR-9906
1" x 1" Fl. Red




SBH-48 2" x 500" Fl. Red (Tape)



UPCR-9401 2" x 3" Fl. Red



UPCR-9904 3 1/2" x 2 1/4" Fl. Red




SBH-61 2" x 4" Fl. Red Package of 100



UPCR-9903 3 1/2" x 2 1/4"
Fl. Red



UPCR-9905 4" x 4" Fl. Red Package of 10




UPCR-9503 2" x 4"


Shipping Labels




HAZ-400 2 1/2" x 2 1/2"



HAZ-403 4" x 3"



HAZ-405 4" x 6"



HAZ-319 6" x 6"

Shipping Labels are supplied 250 per package.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

Specialty Products for Color Coding



SHAMROCK LABELING TAPE

Ideal for color coding, Shamrock Labeling Tape accepts handwritten notations. Labeling Tape will **stick tightly, leaves no residue when removed** and can be transferred from surface to surface. Call for information on tape imprinted to your specifications.

- 2 Lengths: 500", 2160"
- 5 Widths: 1/2", 3/4", 1", 1 1/2", 2"
- 26 Standard Colors:

1 White	2 Yellow	3 Green	4 Red	5 Orange
6 Blue	7 Pink	8 Tan	9 Salmon	10 Copper
11 Rose	12 Violet	13 Lavender	14 Gray	15 Lime
16 Aqua	17 Chartreuse	18 Gold	19 Silver	20 Magenta
21 Dk. Red	22 Maroon	23 Dk. Brown	24 Dk. Blue	25 Dk. Green
26 Black				

- 5 Fluorescent Colors:

27 Fl. Red	28 Fl. Orange	29 Fl. Green	30 Fl. Chartreuse	31 Fl. Pink
------------	---------------	--------------	-------------------	-------------



TAPE IN RAINBOW PACKS

- 500" Tape Rolls
- Your Choice of 3 Widths: 1/2", 3/4", 1"
- 7 Standard Colors in Each Box:

1 White	2 Yellow	3 Green	4 Red
5 Orange	6 Blue	7 Pink	

Tape Width	Cat. No.	Rolls/Pack
1/2"	ST-12-RNBW	24
3/4"	ST-34-RNBW	16
1"	ST-10-RNBW	12



TAPE ROLLS


Tape Width	500" Cat. No.	2160" Cat. No.
1/2"	ST-12	ST-120
3/4"	ST-34	ST-340
1"	ST-10	ST-100
1 1/2"	ST-11	ST-110
2"	ST-20	ST-200

TAPE DISPENSERS



Easy to fill multiple dispensers offer convenient access from a central location.

HAND DISPENSERS



Durable plastic, can be reused many times.


Specialty Products for Color Coding

SHAMROCK DOT LABELS

Versatile Shamrock Dot Labels offer multiple opportunities for color coding and can be imprinted to your specifications. They are supplied on sheets or in roll form. Permanent adhesive is standard; removable adhesive is available upon request. Dots are supplied 1,000 per roll or 1,000 per package of individual sheets.



- Permanent or Removable Adhesive • On Sheets or in Rolls • Plain or Your Special Imprint
- 6 Sizes: 1/4" (SDL-25), 7/16" (SDL-43), 1/2" (SDL-50), 3/4" (SDL-75), 1"(SDL-100), 1 1/4" (SDL-125)
- 31 Distinctive Colors:



Dot labels can be imprinted to your specifications. Contact Shamrock for information.



SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

Terms and Conditions of Sale

All prices subject to change without notice.

A \$25.00 minimum order is required.

Prices are F.O.B. Bellwood, Illinois, unless otherwise specified. Foreign import duties and exchange rates, if any, are additional.

Terms: Net 30 days.

RETURN POLICY

No returns will be accepted without prior authorization. Shamrock's numbered return authorization form must accompany all merchandise returned.

Defective Merchandise: Should any Shamrock product prove to be defective, a claim for return must be placed within 60 days of shipment. Products shown to be defective and returned within this period will be replaced.

Stock Items: No return of stock items will be accepted after six months. A 15% restocking charge is imposed on return of all stock items.

Custom Items: Custom orders may not be returned.

QUANTITIES

Over runs or under runs, not to exceed 10% of quantities ordered, will constitute acceptable delivery.

PLATE CHARGES

Due to natural deterioration of printing plates, and in order to provide only top-quality products, Shamrock must remake its plates every three years. If an order is not placed within that time, a new plate charge must be made upon receipt of the next order.





shamrock
labels

SHAMROCK LABELS

34 Davis Drive, Bellwood, Illinois 60104

National Toll Free: 800-323-0249

California Toll Free: 877-461-6908

Fax Toll Free: 800-248-1907

Email: sales@shamrocklabels.com

Shop Online: shamrocklabels.com