



Tapes & Labels

FOR ANESTHESIA, SURGERY & CENTRAL SERVICE





Shamrock's Mission

Our mission is to supply quality, innovative labeling systems to the healthcare industry that help protect patient safety, fulfill modern technological needs and remain cost effective.



Environmental Policy

It is Shamrock's policy to provide a safe and healthy workplace for our employees and to be an environmentally responsible neighbor in our community. From product manufacture through shipping, in both office and production facilities, Shamrock supports comprehensive programs designed to positively impact the safety of our population and natural resources.



Joint Commission National Patient Safety Goals Sustained with Shamrock Labels

- Goal 2.** Improve the effectiveness of communication among caregivers.
- Goal 3.** Accurately and completely reconcile medications across the continuum of care.
- Goal 7.** Reduce the risk of health care-associated infections.

Anesthesia, Surgery and Central Service

TAPES & LABELS FOR ANESTHESIA	Anesthesia Drug Syringe Labels	CS-2
	Individualized Anesthesia Labels	CS-5
	2D Bar Coded Labels	CS-5
	Blank Labels for On-Demand Printing	CS-5
	Dispensers for Anesthesia Labels	CS-5

TAPES & LABELS FOR SURGERY	Drug Syringe Labels for the Sterile Field	CS-6
	Instrument Marking Tape	CS-7

PRODUCTS FOR SURGERY AND RECOVERY	Wristbands for Same Day Surgery	CS-8
	IV Tubing Control Labels	CS-8
	IV Tubing Change Labels	CS-9
	IV Auxiliary Labels	CS-9
	Communication Labels	CS-9

TAPES & LABELS FOR STERILE PROCESSING	Tapes & Labels for Steam Autoclaves	CS-10
	Blank Thermal Transfer Poly Labels	CS-11
	Labeling Guns & Gun Labels	CS-11
	Communication Labels	CS-11

LABELS FOR PATIENT CHARGE	Blank Labels for Laser Printing	CS-12
	Blank Labels for Thermal Printing	CS-13
	Blank Labels for Dot Matrix Printing	CS-13
	Blank Labels for Tabletop Printing	CS-14
	Three-Part Labels for the Copier	CS-14
	Patient Charge Cards	CS-14

PRODUCTS FOR COLOR CODING	Dot Labels	CS-15
	Labeling Tape	CS-16
	Tape Dispensers	CS-16

Specialty Anesthesia Drug Labels for Syringe Identification

MEET JOINT COMMISSION NATIONAL PATIENT SAFETY GOAL 3








Reconcile Medications Across the Continuum of Care

FOLLOW ASTM STANDARD SPECIFICATION D 4774 - 11e1

Label Size and Color Backgrounds by Drug Classification

FOLLOW ASA STATEMENT ON CONSISTENCY & CLARITY IN SYRINGE LABELING

TO ORDER THE CORRECT IMPRINT

BASIC IMPRINT		WITH Date/Time/Initial	
<p>Tape Labels (1/2" x 500" Tape Roll)</p>  <p>SA-2413 Use drug's catalog number listed.</p>	<p>Pre-Cut Labels (1/2" x 1 1/2" in 500" Roll)</p>  <p>SA-2413-PRE Add the suffix -PRE.</p>	<p>Tape Labels (1/2" x 500" Tape Roll)</p>  <p>SA-2413-DTI Add the suffix -DTI.</p>	<p>Pre-Cut Labels (1/2" x 1 1/2" in 500" Roll)</p>  <p>SA-2413-DTI-PRE Add the suffix -DTI-PRE.</p>
WITH Expiration Date/Time/Initial		PRE-CUT SHEETED LABELS IN PACKS	
<p>Tape Labels (1/2" x 500" Tape Roll)</p>  <p>SA-2413-EXP Add the suffix -EXP.</p>	<p>Pre-Cut Labels (1/2" x 1 1/2" in 500" Roll)</p>  <p>SA-2413-EXP-PRE Add the suffix -EXP-PRE.</p>	<p>Labels available in all three imprint styles, sheeted and supplied in packs. Labels measure 1/2" x 1 1/2" and are supplied 10 labels per sheet, 25 sheets per pack.</p> <p>Basic Imprint With Date/Time/Initial SA-2413-PK SA-2413-DTI-PK</p> <p>With Expiration Date/Time/Initial SA-2413-EXP-PK</p> 	

Labels Coded to ASTM Color Standards

INDUCTION AGENTS

YELLOW WITH BLACK IMPRINT

- SA-208 Brevital 10 mg/mL
- SA-206 Diprivan ___ mg/mL
- SA-209 DIPRIVAN Propofol 10 mg/mL
- SA-200 Etomidate ___ mg/mL
- SA-2006 Etomidate 2 mg/mL
- SA-2001 Ketalar ___ mg/mL
- SA-201 Ketamine ___ mg/mL
- SA-2007 Ketamine 10 mg/mL
- SA-2008 Ketamine 50 mg/mL
- SA-2009 Ketamine 100 mg/mL
- SA-2002 Methohexital 10 mg/mL
- SA-202 Methohexital 1%
- SA-2003 Pentobarbital ___ mg/mL
- SA-207 Pentothal ___ mg/mL
- SA-2004 Pentothal 25 mg/mL
- SA-205 Propofol ___ mg/mL
- SA-2005 Propofol 10 mg/mL
- SA-203 Thiамylal 2.5%
- SA-204 Thiopental 2.5%

BENZODIAZEPINES

ORANGE WITH BLACK IMPRINT

- SA-218 Ativan ___ mg/mL
- SA-210 Diazepam ___ mg/mL

BENZODIAZEPINES

ORANGE WITH BLACK IMPRINT

- SA-211 Lorazepam ___ mg/mL
- SA-212 Midazolam ___ mg/mL
- SA-217 Midazolam 1 mg/mL
- SA-219 Midazolam 5 mg/mL
- SA-214 Valium ___ mg/mL
- SA-213 Versed ___ mg/mL
- SA-215 Versed 1 mg/mL

BENZODIAZEPINE RECEPTOR ANTAGONIST

ORANGE & WHITE DIAGONAL STRIPES WITH BLACK IMPRINT

- SA-2110 Flumazenil ___ mg/mL
- SA-2111 Flumazenil 0.1 mg/mL

MUSCLE RELAXANTS

RED FLUORESCENT WITH BLACK IMPRINT

- SA-227 Anectine ___ mg/mL (rev. imprint)
- SA-2210 ATRAcurium ___ mg/cc
- SA-220 ATRAcurium ___ mg/mL
- SA-2227 ATRAcurium 10 mg/mL
- SA-2216 CISATRAcurium ___ mg/mL
- SA-2228 CISATRAcurium 2 mg/mL

MUSCLE RELAXANTS

RED FLUORESCENT WITH BLACK IMPRINT

- SA-221 Curare ___ mg/mL
- SA-2219 DOXAcurium ___ mg/mL
- SA-222 Gallamine ___ mg/mL
- SA-223 Metocurine ___ mg/mL
- SA-2212 Mivacron ___ mg/mL
- SA-2213 MIVAcurium ___ mg/mL
- SA-2217 Nimbex ___ mg/mL
- SA-2211 Norcuron ___ mg/cc
- SA-2220 Norcuron ___ mg/mL
- SA-224 PANcuronium ___ mg/mL
- SA-2229 PANcuronium 1 mg/mL
- SA-229 Pavulon ___ mg/mL
- SA-2214 ROcuronium ___ mg/mL
- SA-2223 ROcuronium 10 mg/mL
- SA-226 Succinylcholine ___ mg/mL (rev. imprint)
- SA-2225 Succinylcholine 20 mg/mL (rev. imprint)
- SA-228 TRAcrium ___ mg/cc
- SA-2221 TRAcrium ___ mg/mL
- SA-2222 Tubocurarine ___ mg/mL
- SA-225 VECuronium ___ mg/mL
- SA-2224 VECuronium 1 mg/mL
- SA-2215 Zemuron ___ mg/mL
- SA-2226 Zemuron 10 mg/mL

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Labels Coded to ASTM Color Standards

RELAXANT ANTAGONISTS

RED FLU. & WHITE DIAGONAL STRIPES WITH BLACK IMPRINT

- SA-234 Attilirium ___ mg/mL
- SA-230 Edrophonium ___ mg/mL
- SA-2310 Edrophonium 10 mg/mL
- SA-231 Neostigmine ___ mg/mL
- SA-2311 Neostigmine 1 mg/mL
- SA-237 Prostigmine ___ mg/mL
- SA-232 Pyridostigmine ___ mg/mL
- SA-235 Regonol ___ mg/mL
- SA-238 Romazicon ___ mg/mL
- SA-233 Tensilon ___ mg/mL

NARCOTICS

BLUE WITH BLACK IMPRINT

- SA-2417 Alfentanil ___ mcg/mL
- SA-2421 Alfentanil 500 mcg/mL
- SA-240 Alfentanil ___ mg/mL
- SA-248 Demerol ___ mg/mL
- SA-2414 Dilaudid ___ mg/mL
- SA-2422 Dilaudid 0.5 mg/mL
- SA-2415 Duramorph ___ mg/mL
- SA-2423 Duramorph 0.5 mg/mL
- SA-2413 Fentanyl ___ mcg/mL
- SA-2424 Fentanyl 10 mcg/mL
- SA-2420 Fentanyl 50 mcg/mL
- SA-2416 Hydromorphone ___ mg/mL
- SA-2425 Hydromorphone 0.2 mg/mL
- SA-242 Meperidine ___ mg/mL
- SA-2426 Meperidine 10 mg/mL
- SA-243 Morphine ___ mg/mL
- SA-2427 Morphine 1 mg/mL
- SA-2428 Morphine Sulfate 1 mg/mL
- SA-244 Nisentil ___ mg/mL
- SA-247 Nubain ___ mg/mL
- SA-2411 Remifentanil ___ mcg/mL
- SA-2429 Remifentanil 50 mcg/mL
- SA-2410 Stadol ___ mg/mL
- SA-246 Sublimaze ___ mcg/mL
- SA-249 Sufenta ___ mcg/mL
- SA-2419 Sufentanil ___ mcg/mL
- SA-2430 Sufentanil 50 mcg/mL
- SA-245 Sufentanil ___ mg/mL
- SA-2412 Ultiva ___ mcg/mL

NARCOTIC ANTAGONISTS

BLUE & WHITE DIAGONAL STRIPES WITH BLACK IMPRINT

- SA-250 Levallorphan ___ mg/mL
- SA-251 Naloxone ___ mg/mL
- SA-253 Naloxone 0.04 mg/mL
- SA-252 Narcan ___ mg/mL

MAJOR TRANQUILIZERS

SALMON WITH BLACK IMPRINT

- SA-261 Chlorpromazine ___ mg/mL
- SA-263 Droperidol ___ mg/mL
- SA-264 Droperidol 2.5 mg/mL
- SA-262 Inapsine ___ mg/mL
- SA-265 Milrinone ___ mg/mL

VASOPRESSORS

VIOLET WITH BLACK IMPRINT

- SA-288 Adrenalin ___ mcg/mL (rev. imprint)
- SA-287 Adrenalin ___ mg/mL (rev. imprint)
- SA-280 Dobutamine ___ mg/mL
- SA-281 Dopamine ___ mg/mL
- SA-282 Ephedrine ___ mg/mL
- SA-2805 Ephedrine 5 mg/mL
- SA-2806 Ephedrine 10 mg/mL
- SA-2804 Epinephrine ___ mcg/mL (rev. imprint)
- SA-2808 Epinephrine 10 mcg/mL (rev. imprint)
- SA-284 Epinephrine ___ mg/mL (rev. imprint)
- SA-2807 Epinephrine 1 mg/mL (rev. imprint)
- SA-2801 Isuprel ___ mg/mL
- SA-2802 Levophed ___ mg/mL
- SA-285 Neo-Synephrine ___ mcg/mL
- SA-2809 Neo-Synephrine 100 mcg/mL
- SA-289 Neo-Synephrine ___ mg/mL
- SA-2803 Norepinephrine ___ mg/mL
- SA-283 Phenylephrine ___ mg/mL

HYPOTENSIVE AGENTS

VIOLET & WHITE DIAGONAL STRIPES WITH BLACK IMPRINT

- SA-2911 Hydralazine ___ mg/mL
- SA-2912 Hydralazine 4 mg/mL
- SA-2913 Hydralazine 20 mg/mL
- SA-295 Inderal ___ mg/mL
- SA-294 Nipride ___ mcg/mL
- SA-2916 Nitroglycerin ___ mcg/mL
- SA-290 Nitroglycerin ___ mg/mL

HYPOTENSIVE AGENTS

VIOLET & WHITE DIAGONAL STRIPES WITH BLACK IMPRINT

- SA-2917 Nitroprusside ___ mcg/mL
- SA-291 Nitroprusside ___ mg/mL
- SA-292 Phentolamine ___ mg/mL
- SA-293 Trimethaphan ___ mg/mL

LOCAL ANESTHETICS

GRAY WITH BLACK IMPRINT

- SA-302 Bupivacaine ___ mg/mL
- SA-3006 Bupivacaine ___ %
- SA-3008 Bupivacaine 0.25% 2.5 mg/mL
- SA-3007 Bupivacaine 0.5%
- SA-307 Carbocaine ___ %
- SA-3016 Chlorprocaine ___ %
- SA-305 Cocaine ___ %
- SA-300 Lidocaine ___ mg/mL
- SA-3009 Lidocaine ___ %
- SA-3003 Lidocaine 1%
- SA-3010 Lidocaine 1% 10 mg/mL
- SA-3004 Lidocaine 2%
- SA-3011 Lidocaine 2% 20 mg/mL
- SA-3001 Marcaine ___ mg/mL
- SA-306 Marcaine ___ %
- SA-3005 Mepivacaine ___ mg/mL
- SA-309 Naropin ___ mg/mL
- SA-304 Nesacaine ___ %
- SA-301 Procaine ___ mg/mL
- SA-308 Ropivacaine ___ mg/mL
- SA-3012 Ropivacaine ___ %
- SA-3013 Ropivacaine 0.2%
- SA-3014 Ropivacaine 0.25%
- SA-3015 Ropivacaine 0.5%
- SA-3017 Tetracaine ___ mg/mL
- SA-3002 Xylocaine ___ mg/mL
- SA-303 Xylocaine ___ %

ANTICHOLINERGIC AGENTS

GREEN WITH BLACK IMPRINT

- SA-310 Atropine ___ mg/mL
- SA-314 Atropine 0.4 mg/mL
- SA-316 Atropine 1 mg/mL
- SA-313 Atropine Sulfate ___ mg/mL
- SA-311 Glycopyrrolate ___ mg/mL
- SA-315 Glycopyrrolate 0.2 mg/mL
- SA-312 Robinul ___ mg/mL
- SA-317 Scopolamine ___ mg/mL

Specialty Anesthesia Drug Labels for Syringe Identification

Labels Coded to ASTM Color Standards

BETA BLOCKERS

**COPPER BAR
WITH REVERSE IMPRINT**

SA-333 Esmolol ___ mcg/mL
SA-331 Esmolol ___ mg/mL
SA-332 Esmolol 10 mg/mL
SA-334 Labetalol ___ mg/mL
SA-335 Labetalol 5 mg/mL
SA-336 Metoprolol ___ mg/mL
SA-337 Metoprolol 1 mg/mL

MISCELLANEOUS

WHITE WITH BLACK IMPRINT

SA-329 Antibiotic _____
SA-3225 Antibiotic ___ mg/mL
SA-3214 Apresoline ___ mg/mL
SA-3226 Benadryl ___ mg/mL
SA-321 Calcium Chloride
SA-3227 Calcium Chloride ___ mg/mL
SA-3221 Cefazolin ___ mg/mL
SA-3222 Cefotaxime ___ mg/mL
SA-3215 Date Opened _____
SA-326 Decadron ___ mg/mL

MISCELLANEOUS

WHITE WITH BLACK IMPRINT

SA-3228 Dexamethasone ___ mg/mL
SA-3243 Dexamethasone 4 mg/mL
SA-3244 Diphenhydramine 50 mg/mL
SA-3229 Drug ___ mg/mL
SA-3231 Furosemide ___ mg/mL
SA-3245 Furosemide 10 mg/mL
SA-3232 Heparin ___ mg/mL
SA-3260 Heparin 1000 mg/mL
SA-320 Heparin ___ units/mL
SA-324 Heparinized Saline
SA-3246 Insulin 100 units/mL
SA-3210 Lasix ___ mg/mL
SA-3247 Lispro 100 units/mL
SA-3234 Metoclopramide ___ mg/mL
SA-3249 Metoclopramide 5 mg/mL
SA-3250 Metoclopramide 10 mg/mL
SA-3220 Nomodyne ___ mg/mL
SA-3235 Normal Saline ___ mg/mL
SA-3242 Normal Saline (NaCl 0.9%)
SA-3251 Ondansetron ___ mg/mL
SA-3252 Ondansetron 2 mg/mL
SA-3253 Ondansetron 4 mg/mL

MISCELLANEOUS

WHITE WITH BLACK IMPRINT

SA-323 Oxytocin ___ mg/mL
SA-3236 Oxytocin ___ units
SA-3254 Oxytocin 10 units/mL
SA-327 Papaverine
SA-322 Pitocin ___ units/mL
SA-3255 Promethazine 2.5 mg/mL
SA-3256 Promethazine 25 mg/mL
SA-328 Protamine ___ mg/mL
SA-3261 Protamine 10 mg/mL
SA-3218 Reglan ___ mg/mL
SA-3217 Saline _____
SA-3239 Sodium Bicarbonate ___ mg/mL
SA-3240 Sodium Chloride ___ mg/mL
SA-325 Solu-Medrol ___ mg/mL
SA-3258 Terbutaline 1 mg/mL
SA-3223 Toradol ___ mg/mL
SA-3259 Vasopressin ___ units/mL
SA-3211 Verapamil ___ mg/cc
SA-3241 Zofran ___ mg/mL
SA-3216 _____
SA-3219 ___ mg/mL

Drug Labels in Traditional Colors

CAT#	DRUG NAME	COLOR
SA-17	Adrenalin ___%	White
SA-1	Anectine ___%	White
SA-20	Arfonad ___%	Red
SA-2	Atropine ___%	Yellow
SA-21	Atropine Sulfate ___%	Orange
SA-22	Benadryl ___%	Copper
SA-3	Brevital ___%	Green
SA-23	Calcium Chloride ___%	Pink
SA-56	Carbocaine ___%	Gray
SA-4	Curare ___%	Red
SA-24	Demerol ___%	Tan
SA-25	Digoxin ___%	Salmon
SA-26	Dopram ___%	Blue
SA-27	Droperidol ___%	Rose
SA-28	Ephedrine ___%	Violet
SA-30	Heparin ___%	Lime

CAT#	DRUG NAME	COLOR
SA-32	Ketalar ___%	Gray
SA-33	Lanoxin ___%	White
SA-34	Lasix ___%	Yellow
SA-35	Levophed ___%	Green
SA-36	Lidocaine ___%	Red
SA-37	Mannitol ___%	Orange
SA-38	Marcaine ___%	Copper
SA-7	Methedrine ___%	Pink
SA-39	Morphine ___%	Pink
SA-58	Neo-Synephrine ___%	Violet
SA-40	Nitroprusside ___%	Tan
SA-41	Pavulon ___%	Salmon
SA-42	Penicillin ___%	Blue
SA-9	Pentobarbital ___%	Salmon
SA-10	Pentothal ___%	Blue
SA-43	Pitocin ___%	Rose

CAT#	DRUG NAME	COLOR
SA-59	Potassium ___%	White
SA-44	Potassium Chloride ___%	Violet
SA-11	Prostigmine ___%	Rose
SA-60	Protamine Sulfate ___%	White
SA-62	Robinul ___%	Green
SA-45	Sodium Bicarbonate ___%	Lavender
SA-46	Sodium Chloride ___%	Lime
SA-47	Sodium Pentothal ___%	Aqua
SA-48	Solu-Cortef ___%	Gray
SA-49	Solu-Medrol ___%	White
SA-12	Sublimaze ___%	Violet
SA-13	Succinylcholine ___%	Lavender
SA-15	Thorazine ___%	Aqua
SA-52	Valium ___%	Red
SA-55	Versed ___%	Orange

All labels available in rolls of 1/2" x 500" tape, as 1/2" x 1 1/2" pre-cut labels or in packs.
See page CS-2 for ordering information.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



Individualized Anesthesia Labels

Tape and Pre-Cut Syringe Labels tinted, imprinted and/or sized to your specifications.


DRUG NAME




SIZE



COLOR




DOSAGE




Pre-Cut Labels with 2D Bar Code

Labels feature your choice of background color and eye-readable imprint.




Blank Labels are also available for Thermal, Inkjet and Laser Printing


Call Shamrock at 800-323-0249 for complete information on the wide variety of stock blank labels available for on-site printing.



THERMAL




INKJET



LASER

Dispensers

Tape Label Dispensers



MD-500

MD-350

MD-250



SHD-1 holds 1 tape roll.
SHD-2 holds 2 tape rolls.

MD-250 holds 8 rolls of 1/2" wide tape.
MD-350 holds 12 rolls of 1/2" wide tape.
MD-500 holds 16 rolls of 1/2" wide tape.

Pre-Cut Label Dispensers



AUX-10 holds 10 boxes of pre-cut labels 1/2" wide.
AUX-20 holds 20 rolls of pre-cut labels 1/2" wide.

LD-1 holds 2 boxes of labels 1/2" wide.
LD-2 holds 3 boxes of labels 1/2" wide.
LD-3 holds 6 boxes of labels 1/2" wide.
LD-4 holds 12 boxes of labels 1/2" wide.


Specialty Products for Surgery

DRUG SYRINGE LABELS FOR THE STERILE FIELD

Shamrock Drug Syringe Labels feature latex free, glove friendly adhesive that is water and solvent resistant, will not flag or fall off the syringe. Sterile Label Packs feature stock printed labels or labels imprinted to your specifications.

Custom Label Sterile Packs

- Custom Sterile Pack with optional Ruler and Marker
- Labels imprinted to your specifications
- You choose color bar, label color or imprint color



Contact Shamrock for custom pack ordering assistance.

CUSTOM LABELS AVAILABLE IN 4 SIZES

SMALL LABELS - 5/16" x 1-1/4"

- 27 labels per sheet (18 imprinted, 9 blank)
- Up to 2 line imprint with 1 color bar or 1 line imprint with 2 color bars
- 25 characters per line maximum

MEDIUM LABELS - .56" x 1-1/4"

- 15 imprinted labels per sheet
- Up to 2 line imprint with 1 color bar or 1 line imprint with 2 color bars
- 25 characters per line maximum

EXTRA LARGE LABELS - 1" x 2"

- 6 labels per sheet (4 imprinted, 2 blank)
- Up to 3 line imprint with 1 color bar or 2 line imprint with 2 color bars
- 32 characters per line maximum
- Removable adhesive available on Extra Large Labels only for use on stainless steel

LARGE LABELS - 3/4" x 1-1/2"

- 10 imprinted labels per sheet
- Up to 3 line imprint with 1 color bar or 2 line imprint with 2 color bars
- 25 characters per line maximum

Each sheet contains labels of the same size. Sheets of different size labels may be included in one custom sterile pack.

Stock Label Sterile Packs

- Two Label Sheets for stock sterile packs
- Labels pre-printed as illustrated below
- Ruler and Marker options available



Lidocaine Strength _____ Exp. Date _____ Exp. Time _____	Lidocaine & Epinephrine Strength _____ Exp. Date _____ Exp. Time _____	Antibiotic Strength _____ Exp. Date _____ Exp. Time _____	Optiray Strength _____ Exp. Date _____ Exp. Time _____
Marcaine Strength _____ Exp. Date _____ Exp. Time _____	Marcaine & Epinephrine Strength _____ Exp. Date _____ Exp. Time _____	Conray Strength _____ Exp. Date _____ Exp. Time _____	Dye Strength _____ Exp. Date _____ Exp. Time _____
Heparin Strength _____ Exp. Date _____ Exp. Time _____	Thrombin Strength _____ Exp. Date _____ Exp. Time _____	Sterile Water Exp. Date _____ Exp. Time _____	Drug Strength _____ Exp. Date _____ Exp. Time _____
Saline Strength _____ Exp. Date _____ Exp. Time _____	Ancef Strength _____ Exp. Date _____ Exp. Time _____	Drug Strength _____ Exp. Date _____ Exp. Time _____	Medication Added Drug _____ Strength _____ Exp. Date _____ Exp. Time _____
Bacitracin Strength _____ Exp. Date _____ Exp. Time _____	Kantrex Strength _____ Exp. Date _____ Exp. Time _____	Drug Strength _____ Exp. Date _____ Exp. Time _____	Drug Strength _____ Exp. Date _____ Exp. Time _____


SOR-1

SOR-2

STOCK STERILE PACKS

Catalog No.	Pack Contains
SORS-1	2 Label Sheets & Permanent Marker
SORS-2	2 Label Sheets
SORS-3	2 Label Sheets & Ruler
SORS-4	2 Label Sheets, Permanent Marker & Ruler
SORS-5	2 Label Sheets & Dual Tipped Marker
SORS-6	2 Label Sheets, Dual Tipped Marker & Ruler

Options Available for Both Stock and Custom Packs

- Single Tip Permanent Marker 
- Dual Tipped Marker: Skin Marker at one end, Permanent Marker at the other 
- 6" Ruler 

SAVE MONEY • ELIMINATE WASTE

Stock label sheets SOR-1 and SOR-2 and sheets of custom labels may be ordered individually for on-site sterilization. Contact Shamrock for information.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



INSTRUMENT MARKING TAPE

Use Specialty Instrument Marking Tape to designate your surgical instruments by function, doctor or operating room. This durable, vinyl tape with pressure-sensitive adhesive sticks on contact and remains securely attached. It will not inhibit sterilization by steam or ethylene oxide. Resistant to acids, alkalies and oils, it is unaffected by regular cleaning and won't change the weight or feel of the instrument.

Specialty Instrument Marking Tape fits all instrument sizes and conforms to the smallest, most irregular shapes. It is available in 11 colors for maximum flexibility in color coding. Instrument Marking Tape is supplied in 200" rolls, packed one roll per handy, reusable metal tin. Minimum order is six rolls; the colors may be assorted.



SMT-38 (1/4" x 200")



Instrument Marking Tape Colors Available



To indicate color, please add the number of color desired to the catalog number SMT-38. For example: SMT-38 tape in white would be ordered as SMT-38-1, in dark green as SMT-38-25.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

Specialty Labels for Surgery and Recovery

PATIENT WRISTBANDS FOR SAME DAY SURGERY

Shamrock supplies wristbands in adult or pediatric sizes for imprinters, thermal direct printers or laser printers. Snap-on or adhesive closure wristbands, white or color coded are available. A few examples are illustrated below. Please call Shamrock at 800-323-0249 for complete information.

Tyvek® Wristbands



SWB-TV-4 Adult Wristband 1" x 10" Red




SWB-TVSTR-6 Adult Wristband 1" x 10" Blue



Tyvek® wristbands may be ordered printed to your specifications. Contact Shamrock for complete information.

Thermal Direct Wristband



SWB-ACADT-3301 Adult Wristband Length 1.545" x 13.25"




Printed Example

IV TUBING CONTROL LABELS

Labels measure 1/2" x 3" and feature a vertical perforation at the center for easy separation. Place the left half of the label on the IV bag. The right half is placed on the tubing line at the port of entry. Labels are supplied 500 per roll.

ASTM Color Standards for Drug Labels



- TCL-10-1 DATE/TIME/INITIAL (White)
- TCL-10-104 AMIODARONE (White)
- TCL-10-28 ANTIBIOTIC (White)
- TCL-10-29 ATIVAN (Orange)
- TCL-10-107 CISATRACURIUM (Flu. Red)
- TCL-10-108 DOBUTAMINE (Violet)
- TCL-10-32 DOBUTREX (Violet)
- TCL-10-12 DOPAMINE (Violet)
- TCL-10-116 EPINEPHRINE (Violet)

- TCL-10-133 ESMOLOL (Violet Stripes)
- TCL-10-34 HEPARIN (White)
- TCL-10-110 INSULIN (White)
- TCL-10-37 LASIX (White)
- TCL-10-138 LEVOPHED (Violet)
- TCL-10-14 LIDOCAINE (Gray)
- TCL-10-124 LORAZEPAM (Orange)
- TCL-10-5 MIDAZOLAM (Orange)
- TCL-10-118 MILRINONE (White)
- TCL-10-39 NARCOTIC (Blue)

- TCL-10-40 NIPRIDE (Violet Stripes)
- TCL-10-119 NITROGLYCERIN (Violet Stripes)
- TCL-10-120 NITROPRUSSIDE (Violet Stripes)
- TCL-10-122 NOREPINEPHRINE (Violet)
- TCL-10-123 PHENYLEPHRINE (Violet)
- TCL-10-41 PRIMACOR (White)
- TCL-10-2 PROPOFOL (Yellow)
- TCL-10-45 SALINE FLUSH (White)
- TCL-10-126 VASOPRESSIN (White)
- TCL-10-48 VERSED (Orange)



Alternative Colors for Drug Labels

- TCL-10-3 AMICAR (Green)
- TCL-10-4 AMIODARONE (Red)
- TCL-10-31 CARDIZEM (White)
- TCL-10-7 CISATRACURIUM (Pink)
- TCL-10-8 DOBUTAMINE (Tan)
- TCL-10-16 EPINEPHRINE (Teal)
- TCL-10-33 ESMOLOL (White)

- TCL-10-9 FENOLDOPAM (Salmon)
- TCL-10-10 INSULIN (Copper)
- TCL-10-36 INTEGRILIN (White)
- TCL-10-38 LEVOPHED (Green)
- TCL-10-24 LORAZEPAM (Dk. Blue Reverse)
- TCL-10-18 MILRINONE (Gold)
- TCL-10-19 NITROGLYCERIN (Gray)



- TCL-10-20 NITROPRUSSIDE (Magenta Reverse)
- TCL-10-22 NOREPINEPHRINE (Maroon Reverse)
- TCL-10-23 PHENYLEPHRINE (Dk. Brown Reverse)
- TCL-10-43 PROTONIX (White)
- TCL-10-46 SANDOSTATIN (White)
- TCL-10-26 VASOPRESSIN (Black Reverse)

Shamrock can supply Tubing Control labels imprinted and color tinted to match your individual procedures. Call for complete information and design assistance.


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



IV TUBING CHANGE LABELS


Labels measure 3/4" x 3", feature permanent adhesive and are supplied 1,000 labels per roll.



To use, separate liner from tab and wrap it around tubing, leaving printed label on its liner until after tubing change.

IV AUXILIARY LABELS

Labels supplied 500 per roll with a perf down the center for easy fold over.



COMMUNICATION LABELS

Labels supplied 500 per roll unless otherwise noted.

USE ME FIRST

STSA-1 Fl. Red
Supplied in 1/2" x 500" Tape Roll

DIFFICULT INTUBATION

SA-510-PK 1/2" x 3" Fl. Red
Removable Adhesive
Package of 100 Labels

MEDICATION ADDED

DRUG ADDED _____
AMOUNT _____
ADDED BY _____
DATE _____ TIME _____
BAG NO. _____

SN-2000 1 1/2" x 2 1/2" Fl. Red
Roll of 1,000 Labels

MEDICATION ADDED

DRUG _____
AMOUNT _____
ADDED BY _____
DATE _____ TIME _____
THIS LABEL MUST BE AFFIXED TO ALL INFUSION FLUIDS CONTAINING ADDITIONAL MEDICATION

SN-2001 1 1/2" x 2 1/2" Fl. Red
Roll of 1,000 Labels

BIOHAZARD
HAZARD IDENTITY:
BIOHAZARD CAT. NO. SBH-4

SBH-4 1" x 3" Fl. Red

EPIDURAL | EPIDURAL

STB-1 3/4" x 3" Fl. Pink
Perf Down Center For Easy Fold Over

INTRATHECAL | INTRATHECAL

STB-2 3/4" x 3" Fl. Chartreuse
Perf Down Center For Easy Fold Over

MEDICATION ADDED

Name _____ Rm. No. _____
Drug Added _____
Date _____ Time _____ A.M. P.M.
Added By _____
Storage _____
Return to Pharmacy After _____ Hrs.

SMA-23 2" x 3" Roll of 1,000 Labels

RETURN TO ANESTHESIA DEPT. | RETURN TO ANESTHESIA DEPT.

SA-500-PK 3/4" x 4" Yellow with Clear Overlay
Perf Down Center For Easy Fold Over
Package of 100 Labels

WARNING THIS ITEM IS NOT DISPOSABLE

UPCR-5020 3/4" x 2" Fl. Chartreuse

BIOHAZARD
BIOHAZARD CAT. NO. SBH-1

SBH-1 3" x 2 1/4" Fl. Red

Time Out Verification

Correct Patient Correct Procedure Correct Consent
Correct Site/Side/Level/Multiple Structures N/A
Correct Imaging Studies N/A
Correct Implants/Equipment Available N/A

Initials of ALL STAFF involved in timeout verification

Physician _____ Anesthesia _____
R.N. _____ Other _____
Date: _____ Time: _____

TO-1 2" x 4" Fl. Chartreuse

INTRAVENOUS SOLUTION ADDITIVES

Patient _____ Room _____
Date _____ Time _____ By _____
Time Started _____

DRUGS ADDED	STRENGTH

EXP. DATE _____
This label must be affixed to all infusion fluids containing additional medication

SMA-2125 2 1/2" x 4" Roll of 1,000 labels


Specialty Tapes and Labels for Sterile Processing

INDICATOR TAPES AND LABELS FOR STEAM AUTOCLAVES


Plain, printed and color coded tape for use with steam autoclaves allows you to select exactly what you need.

Specialty Autoclave Tape sticks tightly to paper, cloth, glass, metal and plastic. Supplied in 60 yd. rolls, tape can be used to seal the bundle and identify its contents. Mark the tape with pen, pencil or marker; order it pre-printed or with any of six color edges.

Dark, diagonal stripes appear on the tape to provide a visual indication of the bundle's status after processing in a steam autoclave. It is easily removed after use, will not damage the material or leave any residue.




CS-010 3/4" x 60 yds.
CS-020 1" x 60 yds.



CS-003-3 3/4" x 60 yds.

Color coded edges may be printed with a color bar of yellow, green, red, blue, brown or black. Colors may be assorted to achieve quantity pricing. To indicate color border, please put the color number after the Catalog Number CS-003. For example: CS-003-4 for tape with a red border.



CS-113 thru CS-286 3/4" x 60 yd tape rolls
Contact Shamrock for information on additional stock imprints.

Cat. No.	Label
CS-113	Baby Blankets
CS-114	Basin
CS-115	Basin Sets
CS-148	Distilled Water
CS-154	Dressing Tray
CS-158	Emergency Room
CS-169	Forceps
CS-184	Gown-Pack
CS-186	Gowns-O.R.
CS-203	Lap Packs
CS-223	O.B. Instruments
CS-224	O.B. Pack
CS-225	O.R. Instruments
CS-239	Prep Basins
CS-241	Prep Tray
CS-250	Scissors
CS-251	Sheets
CS-261	Suture Removal Set
CS-262	Suture Sets
CS-263	Syringe Size _____ Date _____
CS-276	Towels
CS-280	Tracheotomy Suction Tray
CS-286	Vaginal Speculum Size _____

Shamrock steam autoclave tape may be ordered with borders in colors other than standard or with custom imprints. Contact Shamrock for complete information.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



THERMAL TRANSFER POLY LABELS

Shamrock offers Thermal Transfer Labels, supplied in a white polyester material on a 3" core, designed to accept crisply printed words and barcodes, they withstand autoclaving and storage, instantly identify sterile pack contents.

- White opaque polyester label material.
- Strong adhesive holds labels securely in place throughout use.
- 4 label sizes accommodate most Central Sterile tracking software such as Impress and Rosebud.
- Shamrock can supply thermal poly labels in different sizes and configurations. Call for complete information.



LABELING GUN AND GUN LABELS FOR STERILIZATION

Paxar/Monarch Handheld 1136® Labeling Gun for sterilization prints 2 lines, with 8 characters available for each line. It features an easy to use, snap in ink roller. Choose from three standard printing band setups.

Gun labels for sterilization are supplied white with red ink. Permanent adhesive sticks securely during sterilization and shelving.



STERILIZATION COMMUNICATION LABELS

Labels are supplied in rolls of 1,000 labels unless otherwise specified.

<p>WARNING GAS STERILIZE ONLY</p> <p>SGS-1212 1" x 2 1/2" Fl. Red</p>	<p>READY FOR USE REMOVED FROM AERATOR AT ___:___ PM ON ___-___-___ BY _____</p> <p>SRU-1213 1" x 2 1/2" Fl. Green</p>	<p>RETURN TO CENTRAL SERVICE</p> <p>SRC-1214 1" x 2 1/2"</p>	<p>STERILE Unless Package Opened or Damaged Check Before Use</p> <p>SSI-1218 1" x 2 1/2" Fl. Red</p>	<p>STERILIZED Indefinite Shelf Life Unless Damaged or Opened</p> <p>SSI-1216 1" x 2 1/2"</p>
<p>WAIT! DO NOT USE FOR ___ HOURS THIS PACK WAS GAS STERILIZED ON ___-___-___ DATE AT ___ A.M. ___ P.M.</p> <p>SW-122 1" x 2 1/2"</p>	<p>STERRAD® SYSTEM STERILIZED</p> <p>SSI-1226 1" x 2 1/2" Fl. Orange</p>	<p>RETURN TO CENTRAL SUPPLY</p> <p>SSI-1227 1" x 2 1/2" Fl. Red</p>	<p>IMPLANT RUN WITH BIOLOGICAL TEST</p> <p>SSI-1224 1" x 2 1/2" Fl. Pink</p>	<p>AUTOCLAVED ON ___-___-___ EXPIRES ON ___-___-___ RETURN TO CENTRAL SERVICE FOR RE-STERILIZATION</p> <p>SAE-1215 1" x 2 1/2" Fl. Orange</p>
<p>CLEAN</p> <p>Initials _____ Date _____</p> <p>SSI-1219 2" x 3" 500 per roll</p>	<p>EXPIRATION DATE _____</p> <p>SSI-1225 1" x 2 1/2" Fl. Chartreuse</p>	<p>Consignment Tray for Dr. _____ of _____ Vendor _____ Date _____</p> <p>SSI-1223 1" x 3" Fl. Chartreuse</p>	<p>This Tray Is INCOMPLETE Missing _____ _____ Date _____</p> <p>SSI-1229 1 3/4" x 2 3/4" Fl. Orange 500 per roll</p>	<p>Dust Cover This Item</p> <p>SSI-1228 1" x 2 1/2" Fl. Pink</p>
	<p>OUTER WRAPS NOT STERILE. ONLY INNER CONTENTS STERILE.</p> <p>SSI-1230 1" x 3" Fl. Red</p>	<p>Implant Sterilized with Biological Indicator</p> <p>SSI-1222 1" x 3" Fl. Pink</p>		

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

Specialty Labels for Patient Charge


LASER PRINTER CHARGE LABELS

LASER PIGGYBACK LABELS


Shamrock offers patient charge labels in 8½" x 11" sheets for laser printing. Labels feature piggyback construction and may be tinted any of the thirteen colors listed below.




MSL-4FLAT-PB-2 5/8" x 1.7"
Yellow (carrier size 1" x 2")



MSL-4MT-PB-1 7/8" x 2" (Meditech)
White (carrier size 1" x 2½")




MSL-4FPB-TAB-11 7/8" x 1¼"
Rose (carrier size 1" x 2")




LASER SINGLE-PLY LABELS

Shamrock also offers non-piggyback labels in a wide range of sizes, supplied on 8½" x 11" sheets. Single-ply labels may be ordered in any of the thirteen colors listed below. Call 800-323-0249 for complete information.




LASER-100-1
1" x 2⅝" White

LASER-100-PERF-1 1" x 2⅝" White



MSL-100-FLAT-6 1" x 2" Blue



LASER-3-1 1" x 2⅝" White

LASER-24-1
¾" x 2⅝" White

8 STANDARD COLORS



5 FLUORESCENT COLORS



Colors available for Laser, Thermal and Pinfed Labels for Patient Charge. To order, please add the color number shown to the label's catalog number. For example: LASER-3 in orange would be LASER-3-5.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



THERMAL PRINTER CHARGE LABELS

Patient Charge Thermal Labels, available in the colors listed on page CS-12, are supplied on material for thermal direct printing. Call for information on labels for thermal transfer printing and additional sizes available. Thermal labels are supplied on a 3" core unless otherwise indicated.


PIGGYBACK THERMAL LABELS



THERM-PB548R-2
1" x 2" Yellow
(carrier size 1 1/8" x 2 1/4")
2,000/roll




THERM-PB992R-2 3/4" x 1 3/4"
Yellow (carrier size 7/8" x 1 7/8")
4,000/roll



MSL-10-5 3/4" x 1 3/8"
Orange color bar
(carrier size 7/8" x 1 7/8")
2,000/roll 4" core



THERM-PB1142R
7/8" x 1 7/8"
(carrier size 1" x 2")
2,000/roll



MSL-X-2 1 5/16" x 1 9/16"
Yellow
(carrier size 1 5/16" x 1 9/16")
3,000/roll 1 5/8" core




THERM-PB1148R-7
1" x 2" Pink
(carrier size 1 1/4" x 2 1/4")
1,500/roll

SINGLE PLY THERMAL LABELS



THERM-487R
1/2" x 1 1/2"




THERM-520R
3/4" x 2"



THERM-396R
7/8" x 2"



THERM-513R
1" x 2"



THERM-1791R-DO
1 1/4" x 2 1/4"
1,000 per roll
1" core




THERM-1416R-DO
2 1/8" x 1"
500 per roll
1" core

DOT MATRIX PRINTER CHARGE LABELS

Patient Charge Pinfed Labels may be ordered in any of the 8 standard and 5 fluorescent colors illustrated on page CS-12. Call for information on additional sizes available.


PIGGYBACK PINFED LABELS




MSL-2PF-2 3/4" x 1 1/2" Yellow
(carrier size 7/8" x 1 11/16")



MSL-3PF-7 5/8" x 1.7" Pink (carrier size 1 5/16" x 2")



MSL-4PF-3 5/8" x 1.7" Green (carrier size 1 5/16" x 2")




MSL-5PF-6 3/4" x 1 3/4" Blue (carrier size 7/8" x 2")



MSL-5PF-STRPD-6 3/4" x 1 3/4" Blue (carrier size 7/8" x 1 15/16")

NOTE: MSL-2PF thru MSL-5PF series labels are supplied on standard or flexible piggyback label material. To order labels in flexible material, place "-F-" in the catalog number before the color number. For example: The MSL-4PF-3 (green) would be ordered as MSL-4PF-F-3 (flexible material - green).



MSL-PEEL-5 2/3" x 1 7/10" Orange with pull tab (labels & carrier same size)



PCS-PB1516PF-7 1 5/16" x 1 1/2" Pink with pull tab (labels & carrier same size)



PCS-PB1710PF-2 1 5/16" x 1 7/10" Yellow with pull tab (labels & carrier same size)

SINGLE PLY PINFED LABELS



SNO-31
1 5/16" x 2 1/2"



SNO-32 1 5/16" x 2 3/4"




SNO-42 1 7/16" x 2 3/4"

Specialty Labels for Patient Charge



LABELS FOR TABLETOP IMPRINTERS

Specialty Piggyback Labels for Tabletop Imprinters produce clearly printed, smudge resistant patient charge labels. For accurate color coding, these labels are available in your choice of the seven distinct colors shown below.




MSL-2 3/4" x 1 3/8" Yellow (carrier size 7/8" x 1 5/8")




- MSL-1 White
- MSL-2 Yellow
- MSL-3 Green
- MSL-4 Red
- MSL-5 Orange
- MSL-6 Blue
- MSL-7 Pink

THREE-PART LABELS FOR THE COPIER


Specialty Three-Part Labels have two perforations through the face and backing to insure easy separation of the three sections. Space is provided on two sections of the labels to record the charge item's code number. One section also accommodates the patient's name and room number. Each Three-Part Label is available in your choice of the seven colors shown below.



PCS-11-5 1 1/2" x 2 3/4" Orange (sheet form)



PCS-10-6 1" x 2 3/4" Blue (sheet form)




PCS-20-3 2" x 2 3/4" Green (sheet form)

- 1 White
- 2 Yellow
- 3 Green
- 4 Red
- 5 Orange
- 6 Blue
- 7 Pink

To order PCS series labels in color, please add the color number to the PCS label catalog number. For example: PCS-11 in green would be ordered as PCS-11-3. To order PCS-10 in green use catalog number PCS-10-3.

PATIENT CHARGE CARDS

Your Specialty Patient Charge Label is easily transferred to these handy, 5" cards when the chargeable item is dispensed. Cards are printed to aid label placement.



PCC-4 5" x 8"

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM


Specialty Products for Color Coding

Shamrock Dot Labels

Dot Labels offer multiple opportunities for color coding shelves, bins or chargeable items.
Dot Labels are supplied in lots of 1,000 Dots of the same color.



- Permanent or Removable Adhesive • On Sheets or in Rolls • Plain or Your Special Imprint
- 6 Sizes: 1/4" (SDL-25), 7/16" (SDL-43), 1/2" (SDL-50), 3/4" (SDL-75), 1" (SDL-100), 1 1/4" (SDL-125)
- 26 Standard Colors • 5 Fluorescent Colors



Dot labels can be imprinted to your specifications. Call Shamrock for information.



SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

Specialty Products for Color Coding

Shamrock Labeling Tape

Ideal for color coding, versatile Shamrock Labeling Tape accepts handwritten notations or can be imprinted to your specifications. Tape sticks tightly, leaves no residue and can be transferred from surface to surface.

- **2 Lengths:** 500", 2160"
- **5 Widths:** 1/2", 3/4", 1", 1 1/2", 2"
- **26 Standard Colors:**

1 White	2 Yellow	3 Green	4 Red	5 Orange
6 Blue	7 Pink	8 Tan	9 Salmon	10 Copper
11 Rose	12 Violet	13 Lavender	14 Gray	15 Lime
16 Aqua	17 Chartreuse	18 Gold	19 Silver	20 Magenta
21 Dk. Red	22 Maroon	23 Dk. Brown	24 Dk. Blue	25 Dk. Green
26 Black				

- **5 Fluorescent Colors:**

27 Fl. Red	28 Fl. Orange	29 Fl. Green	30 Fl. Chartreuse	31 Fl. Pink
------------	---------------	--------------	-------------------	-------------



TAPE IN RAINBOW PACKS

- 500" Tape Rolls
- Your Choice of 3 Widths: 1/2", 3/4", 1"
- 7 Standard Colors in Each Box:

1 White	2 Yellow	3 Green	4 Red
5 Orange	6 Blue	7 Pink	

Tape Width	Cat. No.	Rolls/Pack
1/2"	ST-12-RNBW	24
3/4"	ST-34-RNBW	16
1"	ST-10-RNBW	12



TAPE ROLLS


Tape Width	500" Cat. No.	2160" Cat. No.
1/2"	ST-12	ST-120
3/4"	ST-34	ST-340
1"	ST-10	ST-100
1 1/2"	ST-11	ST-110
2"	ST-20	ST-200

TAPE DISPENSERS



Easy to fill multiple dispensers offer convenient access from a central location.

HAND DISPENSERS



Durable plastic, can be reused many times.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

Terms and Conditions of Sale

All prices subject to change without notice.

A \$25.00 minimum order is required.

Prices are F.O.B. Bellwood, Illinois, unless otherwise specified. Foreign import duties and exchange rates, if any, are additional.

Terms: Net 30 days.

RETURN POLICY

No returns will be accepted without prior authorization. Shamrock's numbered return authorization form must accompany all merchandise returned.

Defective Merchandise: Should any Shamrock product prove to be defective, a claim for return must be placed within 60 days of shipment. Products shown to be defective and returned within this period will be replaced.

Stock Items: No return of stock items will be accepted after six months. A 15% restocking charge is imposed on return of all stock items.

Custom Items: Custom orders may not be returned.

QUANTITIES

Over runs or under runs, not to exceed 10% of quantities ordered, will constitute acceptable delivery.

PLATE CHARGES

Due to natural deterioration of printing plates, and in order to provide only top-quality products, Shamrock must remake its plates every three years. If an order is not placed within that time, a new plate charge must be made upon receipt of the next order.



813



shamrock
labels

SHAMROCK LABELS

34 Davis Drive, Bellwood, Illinois 60104

National Toll Free: 800-323-0249

California Toll Free: 877-461-6908

Fax Toll Free: 800-248-1907

Email: sales@shamrocklabels.com

Shop Online: shamrocklabels.com