



---

# Tapes & Labels

**FOR THE LABORATORY**

---





## Shamrock's Mission

Our mission is to supply quality, innovative labeling systems to the healthcare industry that help protect patient safety, fulfill modern technological needs and remain cost effective.



## Environmental Policy

It is Shamrock's policy to provide a safe and healthy workplace for our employees and to be an environmentally responsible neighbor in our community. From product manufacture through shipping, in both office and production facilities, Shamrock supports comprehensive programs designed to positively impact the safety of our population and natural resources.



## Joint Commission National Patient Safety Goals For the Laboratory Sustained with Shamrock Labels

- Goal 1.** Improve the accuracy of patient identification.
- Goal 2.** Improve the effectiveness of communication among caregivers.
- Goal 7.** Reduce the risk of health care-associated infections.

# Laboratory

<b>TAPE &amp; DOTS FOR COLOR CODING</b>	Labeling Tape	L-2
	Dot Labels	L-4
<b>HI/LO TEMP &amp; INDICATOR TAPES</b>	Hi/Lo Temp Tape & Hi Temp Indicator Tape	L-5
	Steam Indicator Tape	L-5
<b>CRYOGENIC LABELS</b>	Labels for all Cryogenic Temperatures	L-6
	CryoPens & Thermal Ribbons	L-7
<b>SPECIMEN IDENTIFICATION &amp; TESTING LABELS</b>	Specimen & Pathology Labels	L-8
	Urine Collection & Glucose Tolerance Labels	L-8
	Bacterial Genus Tape	L-9
<b>RED BORDER LABELS AND CONTROL LABELS</b>	Red Border & Reagent Labels	L-9
	Calibration Labels & Lot Signal Labels	L-10
<b>TAMPER RESISTANT LABELS</b>	Tamper Resistant Labels & Tapes	L-11
<b>CONSECUTIVE NUMBER LABELS</b>	Tape & Pre-Cut Number Labels	L-12
	CNI System & Piggyback Number Labels	L-13
	Solvent Resistant Number Labels	L-13
	Pre-Printed Bar Coded Number Labels	L-14
	Blank Thermal, Laser, Dot Matrix Printer Labels	L-15
<b>MICROSCOPE SLIDE LABELS</b>	Slide Labels	L-16
	Overlay Label	L-17
<b>PHOTOGRAPHY LABELS</b>	Metric Measure Labels	L-17
<b>WATER TESTING LABELS</b>	Preservative Labels	L-17
<b>ANIMAL CARE LABELS</b>	Cage Cards, Communication Labels	L-18
	Medication & Drug Syringe Labels	L-18
<b>WARNING LABELS, TAPES &amp; BAGS</b>	Biohazard Warning Labels & Tapes	L-19
	Waste Control Labels	L-21
	Radioactive Warning & Chemical Hazard Labels	L-21
	GHS Warning Labels	L-22
	NFPA Chemical Hazard Labels	L-24
	Caution Labels & Tapes	L-25
<b>SHIPPING LABELS &amp; TRANSPORT BAGS</b>	D.O.T. Labels	L-26
	Reclosable Poly Bags	L-26
<b>COMMUNICATION LABELS</b>	Departmental Information Labels	L-27
<b>MULTI-PURPOSE BLANK LABELS</b>	Roll Form Labels & Imprinter Labels	L-33
	VAC, Satin Cloth & Tyvek® Labels	L-34
<b>TAPE &amp; LABEL DISPENSERS</b>	Tape & Label Dispensers	L-35
<b>CUSTOM LABELS</b>	Ordering Custom Tape & Pre-Cut Labels	L-36
<b>LABELS FOR ON-DEMAND PRINTING</b>	Introduction to Blank Labels	L-38
	Thermal Printer Labels	L-40
	Laser Printer Labels	L-62
	Dot Matrix Printer Labels	L-66

# Specialty Products for Color Coding

## Shamrock's Multi-Purpose Labeling Tape

Provides the Perfect Solution to a Wide Variety of Labeling Needs.

- Versatile
- Dependable
- Transferable
- Acid Resistant
- Waterproof
- Oilproof



### IDEAL FOR COLOR CODING

**Adheres** to any dry surface – glass, plastic, metal, paper and rubber.

**Sticks** throughout the most demanding procedures – is oilproof, waterproof and acid resistant – can be steam autoclaved.

**Holds** in a wide range of temperatures – transfers cleanly from one surface to another, leaving no messy residue.

**Accepts** handwritten notations with pen, pencil or marker, or order pre-printed to your specifications.

### 31 Colors

1 White	11 Rose	22 Maroon
2 Yellow	12 Violet	23 Dark Brown
3 Green	13 Lavender	24 Dark Blue
4 Red	14 Gray	25 Dark Green
5 Orange	15 Lime	26 Black
6 Blue	16 Aqua	27 Fl. Red
7 Pink	17 Chartreuse	28 Fl. Orange
8 Tan	18 Gold	29 Fl. Green
9 Salmon	19 Silver	30 Fl. Chartreuse
10 Copper	20 Magenta	31 Fl. Pink

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



**Individual Rolls** of Shamrock Labeling Tape are available in your choice of:

- 31 Colors
- 2 Lengths
- 5 Widths

Width	500" Cat. No.	2160" Cat. No.
1/2"	ST-12	ST-120
3/4"	ST-34	ST-340
1"	ST-10	ST-100
1 1/2"	ST-11	ST-110
2"	ST-20	ST-200

**Rainbow Packs** of tape rolls in white, yellow, green, red, orange, blue and pink are also available.

- 7 Colors
- 500" Tape Rolls
- Your Choice of 3 Widths

Cat. No.	Width	Rolls per Pack
ST-12-RNBW	1/2"	24
ST-34-RNBW	3/4"	16
ST-10-RNBW	1"	12




### MULTIPLE ROLL TAPE DISPENSERS

Easy to fill dispensers, available in four sizes, offer convenient access to tape from a central location.



### HAND DISPENSERS

Durable plastic, can be reused many times.




# Specialty Products for Color Coding

## Shamrock Dot Labels

Versatile Dot Labels offer multiple opportunities for color coding in the laboratory, for instance by test, physician or department. Dot Labels are supplied in lots of 1,000 Dots of the same color.



- Permanent or Removable Adhesive • On Sheets or in Rolls • Plain or Your Special Imprint
- 6 Sizes: 1/4" (SDL-25), 7/16" (SDL-43), 1/2" (SDL-50), 3/4" (SDL-75), 1" (SDL-100), 1 1/4" (SDL-125)
- 31 Distinctive Colors



Dot labels can be imprinted to your specifications. Call Shamrock for information.



SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



# Specialty Hi/Lo Temperature and Indicator Tapes



## HI/LO TEMPERATURE TAPE



Shamrock Specialty Hi/Lo Tape won't tear, must be cut with scissors.


Cloth textured tape sticks to any container, withstands acid, alkalis, oils, greases, dry ice and liquid gases, won't crack or peel in dry oven or cryogenic temperatures. Write on it with Shamrock's CryoPen-26 or ask for a price quotation on tape imprinted to your specifications. Tape, available in three widths, is supplied in 250" rolls.

SHLT-12 (1/2" Tape), SHLT-34 (3/4" Tape), SHLT-10 (1" Tape)



**CryoPen-26**  
Supplied in packs of 4, ink won't fade or smudge.

## HI TEMPERATURE INDICATOR TAPE



Before



After

Imprinted with heat-sensitive indicator ink, cloth textured tape sticks to any surface or shape, withstands acid, alkalis, oils, greases and shows the word "Sterilized" in bold letters after standard dry oven processing. Supplied in 250" rolls, the tape is available in 3 widths. May be custom imprinted at your request.

SHTI-12 (1/2" Tape), SHTI-34 (3/4" Tape), SHTI-10 (1" Tape)


## INDICATOR TAPE FOR STEAM AUTOCLAVES

Versatile tape, available plain or imprinted, indicates completion of normal steam sterilization processing. Tape accepts handwritten notations, may be used to seal and identify laboratory materials before and after processing.

### Blank Indicator Tape

2 Widths – 3/4", 1"

Length – 2160"



CS-010 (3/4" x 2160")  
CS-020 (1" x 2160")

3 Widths – 1/2", 3/4", 1"  
2 Lengths – 500", 2160"



SIT-34 (3/4" x 500"), SIT-10 (1" x 500"),  
SIT-340 (3/4" x 2160"), SIT-100 (1" x 2160")




SIT-12 (1/2" x 500"),  
SIT-120 (1/2" x 2160")

### Imprinted Indicator Tape


3 Imprint Choices

2 Widths – 3/4", 1"


2 Lengths – 500", 2160"



SIL-34-B-4 (3/4" x 500"), SIL-10-B-4 (1" x 500"),  
SIL-340-B-4 (3/4" x 2160"), SIL-100-B-4 (1" x 2160")  
Available with Red, Yellow, Green, Blue or Black Bar for Color Coding



SIL-34 (3/4" x 500"), SIL-10 (1" x 500"),  
SIL-340 (3/4" x 2160"), SIL-100 (1" x 2160")




SIL-34-1 (3/4" x 500"), SIL-10-1 (1" x 500"),  
SIL-340-1 (3/4" x 2160"), SIL-100-1 (1" x 2160")

# Specialty Labels for Cryogenic Procedures


## Laser Printer Labels on 8-1/2" x 11" Sheets

### LIQUI-STICK LABEL LASER SHEETS -196°C


- Withstand Liquid or Vapor Phase Liquid Nitrogen
- Resistant to UV and Gamma Irradiation • 7 Color Choices for Labels Starred
- Smudge Resistant Surface for Laser Printing and CryoPen Marking



LASER-CL12T  
(7 Color Choices)




LASER-CL23T  
(7 Color Choices)




LASER-CL48T  
(7 Color Choices)



LASER-CL60 (White)  
for Tubes & Tube Tops



SLS-LASER-CL32T  
(White)



LASER-TRCL12  
(Transparent)

LASER-TRCL23  
(Transparent)

Catalog Number	Label Size	Description	Labels/Sheet	Sheets/Pack
* LASER-CL12T	.5" x .94"	7 Color Choices	212	16
* LASER-CL23T	.5" x 1.28"	7 Color Choices	176	16
* LASER-CL48T	.4333" Dia.	7 Color Choices	120	16
LASER-CL60	.94" x 1.125" & .44" Dia.	White Only	84	16
SLS-LASER-CL32T	.78" x .94"	White Only	70	16
LASER-TRCL12	.5" x .94"	Transparent	84	16
LASER-TRCL23	.5" x 1.28"	Transparent	70	16


\* Seven Colors Available 1 White 2 Yellow 3 Green 4 Red 5 Orange 6 Blue 7 Pink

To indicate label color when ordering, place the color number after the catalog number above.  
For instance LASER-CL12T in yellow would be ordered as LASER-CL12T-2.

### LASER SMUDGEPROOF KIMDURA® SHEETS -54°C

- Withstand Temperatures as low as -54°C • Supplied on 8½" x 11" Sheets/5,000 Labels per Carton
- Glove friendly, smudge resistant • 13 Color Choices

Catalog Number	Label Size	Labels/Sheet	Sheets/Ctn
LASER-KDS-43	7/16" Dia.	212	24
LASER-KDS-59	1/2" x 9/10"	176	29
LASER-KDS-5125	1/2" x 1 1/4"	120	42
LASER-KDS-515	1/2" x 1 1/2"	100	50
LASER-KDS-759	3/4" x 9/10"	104	49
LASER-KDS-7512	3/4" x 1 1/4"	84	60
LASER-KDS-7515	3/4" x 1 1/2"	70	72



LASER-KDS-7512-1

LASER-KDS-7515-1

LASER-KDS-5125-3  
(Green)

#### STANDARD COLORS

1 White 2 Yellow 3 Green 4 Red 5 Orange 6 Blue 7 Pink 11 Rose

#### FLUORESCENT COLORS

27 Fl. Red 28 Fl. Orange 29 Fl. Green 30 Fl. Chartreuse 31 Fl. Pink

To order labels in color, please add the color number shown to the label's catalog number.  
For instance LASER-KDS-43 in fluorescent red would be LASER-KDS-43-27.

**Printing template information is supplied with all Laser Labels.**


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM





## CRYO-STICK LABELS IN COLORS ON SHEETS OR ROLLS -80°C



FCS-43



FCS-43-4  
(Red)



FCS-43-2  
(Yellow)




FCS-43-30  
(Fl. Chartreuse)

7/16" Dia. Dots (32 Dots per 2 1/2" x 4" Sheet)



FCS-5125-11  
(Rose)




FCS-5125-6  
(Blue)



FCS-5125-3  
(Green)

1/2" x 1 1/4" Labels (32 Labels per 5" x 4" Sheet)



FCS-5125

8 Standard Colors · 5 Fluorescent Colors

Glove Friendly · Smudge Resistant · Supplied in Lots of 1,000 Labels

### STANDARD AND FLUORESCENT COLORS



To order labels in color, please add the color number shown to the label's catalog number.  
For instance FCS-43 in fluorescent red would be FCS-43-27.

## LABELS FOR ALL CRYOGENIC TEMPERATURES ON SHEETS OR ROLLS

Glove Friendly · Smudge Resistant · Supplied in Lots of 500 or 1,000 Labels

### Cryogenic Labels in Rolls for Thermal Printers



Rolls Cat. No.	Label Size	Materials
KDS-43R FCS-43R LCS-43R	7/16" Dia.	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-59R FCS-59R LCS-59R	1/2" x 9/10"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-5125R FCS-5125R LCS-5125R	1/2" x 1 1/4"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-515R FCS-515R LCS-515R	1/2" x 1 1/2"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-759R FCS-759R LCS-759R	3/4" x 9/10"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7512R FCS-7512R LCS-7512R	3/4" x 1 1/4"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7515R FCS-7515R LCS-7515R	3/4" x 1 1/2"	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C

### Cryogenic Labels in Sheets for Hand Marking




Sheets Cat. No.	Label Size	Sheet Size	Materials
KDS-43 FCS-43 LCS-43	7/16" Dia.	2 1/4" x 4" 32 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-59 FCS-59 LCS-59	1/2" x 9/10"	4 1/2" x 4" 40 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-5125 FCS-5125 LCS-5125	1/2" x 1 1/4"	5" x 4" 32 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-515 FCS-515 LCS-515	1/2" x 1 1/2"	4" x 4 1/2" 24 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-759 FCS-759 LCS-759	3/4" x 9/10"	4 1/2" x 4 1/2" 36 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7512 FCS-7512 LCS-7512	3/4" x 1 1/4"	4 1/2" x 5" 24 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C
KDS-7515 FCS-7515 LCS-7515	3/4" x 1 1/2"	4 1/2" x 4 1/2" 18 labels/sheet	Dura-Stick -40°C Cryo-Stick -80°C Liqui-Stick -196°C

Use catalog numbers above for lots of 1,000 labels. For 500 labels add the suffix "-D" to the appropriate catalog number (KDS-43R-D, KDS-43-D).

### Wrap-Around Labels


These unique labels for sub-zero use feature a matte white area for thermal transfer printing or marking with CryoPen-26.

**Note:** The extended clear portion of each label covers and protects printed information. The ACTG-WAPC-1 label features an indexing notch between the printable portion of each label.



### CryoPen-26

Supplied in packs of 4, ink won't fade or smudge.



### Thermal Transfer Ribbons

Shamrock offers a wide variety of ribbons specifically designed to produce an imprint that remains clear and legible throughout cryogenic procedures. See page L-61 or call Shamrock for complete information.

# Specialty Labels for Specimen Testing

## PRE-CUT AND TAPE IDENTIFICATION LABELS

Smudge proof labels accept written notations in pen, pencil or marker, adhere instantly to glass, plastic or waxed surfaces. Labels are supplied on tape or pre-cut, round cornered labels in 500" rolls.

NAME \_\_\_\_\_  
 \_\_\_\_\_  
 NAME \_\_\_\_\_  
 Rm. No. \_\_\_\_\_  
 SL-1-PRE Rm. No. \_\_\_\_\_ DATE \_\_\_\_\_  
 SL-1 1" x 2 1/4" (tape)

NAME \_\_\_\_\_  
 \_\_\_\_\_  
 SPECIMEN NAME \_\_\_\_\_  
 Rm. No. \_\_\_\_\_  
 SL-2-PRE Rm. No. \_\_\_\_\_ DATE \_\_\_\_\_  
 SL-2 1" x 2 1/4" (tape)

NAME \_\_\_\_\_  
 HOSP. NO. \_\_\_\_\_  
 NAME \_\_\_\_\_  
 Rm. No. \_\_\_\_\_  
 HOSP. NO. \_\_\_\_\_  
 SL-3-PRE Rm. No. \_\_\_\_\_ DATE \_\_\_\_\_  
 SL-3 1" x 2 1/4" (tape)

**PATHOLOGY SPEC**

NAME	SEX	NAME	PATHOLOGY SPEC
DATE		DATE	SEX
ROOM		ROOM	AGE
DOCTOR		DOCTOR	NUMBER
SPECIMEN		SPECIMEN	

SL-4-PRE SL-4 1 1/2" x 3 1/2" (tape)

IMPRINT YOUR NAME HERE  
**PATHOLOGY SPEC**

NAME	IMPRINT YOUR NAME HERE
DATE	<b>PATHOLOGY SPEC</b>
ROOM	
DOCTOR	NAME
SPECIMEN	DATE
	SEX
	AGE
	NUMBER

SL-5-PRE SL-5 2" x 3 1/2" (tape)




## PRE-CUT IDENTIFICATION LABELS

Supplied 500 labels per roll. Available as pre-cut labels only.

NAME \_\_\_\_\_ Age \_\_\_\_\_  
 Number \_\_\_\_\_ Date \_\_\_\_\_  
 Doctor \_\_\_\_\_  
 UPCR-6052 1" x 2 1/4"

NAME \_\_\_\_\_  
 Hosp. No. \_\_\_\_\_ Tech. \_\_\_\_\_  
 Location \_\_\_\_\_ Date \_\_\_\_\_ Time \_\_\_\_\_  
 UPCR-6054 1" x 2 1/4"

**PATHOLOGY SPECIMEN**

Name \_\_\_\_\_ Date \_\_\_\_\_  
 Specimen \_\_\_\_\_  
 Doctor \_\_\_\_\_  
 10% FORMALIN POISON   
 UPCR-7021 1" x 3"

**PATHOLOGY SPECIMEN**

UCPR-9404 2" x 3"

## URINE COLLECTION LABELS

Clearly printed in black, these labels offer at a glance identification, instructions and warnings. Where specific information is needed, make notations in pen, pencil or marker. Supplied 500 labels per roll.

**URINALYSIS**  
 UPCR-1324 5/16" x 1 1/4"  
 Fl. Chartreuse

**URINE**  
 UPCR-5010 3/4" x 2"

**CAUTION:**  
 Caustic Preservative Added  
 UPCR-6062 1" x 2 1/4"  
 Fl. Red

**GENERAL 24-HOUR URINE COLLECTION**  
 No Preservative Added  
 Keep Refrigerated  
 UPCR-6063 1" x 2 1/4"  
 Fl. Chartreuse

**URINE FOR CULTURE**  
 UPCR-6064 1" x 2 1/4"  
 Fl. Pink

**URINE SPECIMEN NEEDED**  
 UPCR-6065 1" x 2 1/4"  
 Fl. Chartreuse

**CAUTION:**  
 THIS BOTTLE CONTAINS DANGEROUS ACID WHICH WILL BURN THE SKIN.  
**DO NOT**  
 URINATE DIRECTLY INTO THIS BOTTLE!!  
 UPCR-9703 2 1/2" x 4" Fl. Red

**24 HOUR URINE COLLECTION PRESERVATIVE**

TEST \_\_\_\_\_  
 PT. \_\_\_\_\_ RM. NO. \_\_\_\_\_  
 DATE & TIME STARTED \_\_\_\_\_  
 DATE & TIME COMPLETED \_\_\_\_\_  
 INITIALS \_\_\_\_\_  
 SB-375 1 3/4" x 3"

**24 HOUR URINE COLLECTION PRESERVATIVE**

TEST \_\_\_\_\_  
 PT. \_\_\_\_\_ RM. NO. \_\_\_\_\_  
 DATE & TIME STARTED \_\_\_\_\_  
 DATE & TIME COMPLETED \_\_\_\_\_  
 INITIALS \_\_\_\_\_  
 UPCR-9105 1 1/2" x 3"


**24 HOUR URINE COLLECTION**

TESTS REQUESTED: 1. \_\_\_\_\_  
 2. \_\_\_\_\_  
 3. \_\_\_\_\_  
 PRESERVATIVE: \_\_\_\_\_  
 ADDED BY: \_\_\_\_\_  
 COLLECTION STARTED: \_\_\_\_\_  
 COLLECTION COMPLETE: \_\_\_\_\_  
 UPCR-9106 1 1/2" x 3"  
 Fl. Chartreuse

**KEEP SPECIMEN REFRIGERATED DURING COLLECTION**  
 UPCR-9107 1 1/2" x 3"  
 Fl. Pink

## GLUCOSE TOLERANCE TESTING LABELS

Shamrock's handy, eye-catching labels identify each timed specimen, reduce chance of error.



### Tablet Form Labels

Take one sheet from the convenient tablet, attach one to the patient's requisition form and use individual labels at the required times.

Tablets hold 50 sheets of 8 labels per sheet and are sold in packs of 4 tablets each, providing enough labels to identify 1,600 specimens.

GT-8 1/2" x 1 1/4"  
 Fl. Orange

### Labels on Rolls

Labels feature bright, fluorescent color, are supplied in rolls of 500 labels each, measure 5/16" x 1 1/4".

**GLUCOSE 1 HR.**

UCPR-1311  
 Fl. Green

**GLUCOSE 2 HR.**

UCPR-1312  
 Fl. Green

**GLUCOSE 3 HR.**

UCPR-1313  
 Fl. Green

**GLUCOSE 4 HR.**

UCPR-1314  
 Fl. Green

**GLUCOSE 5 HR.**


UCPR-1315  
 Fl. Green

**2 HOUR P.P.**

UCPR-1316  
 Fl. Green

**GLUCOSE TOLERANCE**

UCPR-1317  
 Fl. Pink



SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



## LABELING TAPE WITH BACTERIAL GENUS NAME

Black ink on yellow tape provides instant identification. Information may be added in pen, pencil or marker. Labels are supplied in 1/2" x 500" rolls.

STGS-1 *Aspergillus* \_\_\_\_\_

STGS-2 *Bacillus* \_\_\_\_\_

STGS-6 *Enterobacter aerogenes*

STGS-7 *Enterococcus faecalis*

STGS-8 *Escherichia coli*

STGS-9 *Micrococcus luteus*

STGS-13 *Penicillium* \_\_\_\_\_

STGS-14 *Proteus* \_\_\_\_\_

**Aspergillus** \_\_\_\_\_

STGS-1

STGS-16 *Pseudomonas aeruginosa*

STGS-18 *S. aureus*

STGS-19 *S. epidermidis*


STGS-20 *Saccharomyces*

## Specialty Red Border Labels and Tape

### RED BORDER TAPES AND LABELS


Your choice of four sizes in labeling tape or as pre-cut labels. May be marked with pen, pencil, or marker.

Tape labels transfer from one surface to another, leaving no messy residue. They are supplied in 500" rolls. Pre-Cut labels feature permanent adhesive, are supplied 500 labels per roll. Please specify tape or pre-cut when ordering.



SBL-20 2" x 3"  
SBL-11 1 1/2" x 2 1/2"  
SBL-10 1" x 1 5/8"  
SBL-34 3/4" x 1"

PRE-CUT LABELS



SBLP-20 2" x 3"  
SBLP-11 1 1/2" x 2 1/2"  
SBLP-10 1" x 1 5/8"  
SBLP-34 3/4" x 1"

## Specialty Control Labels

### REAGENT LABELS

Won't flag or fall off any reagent container, may be written on with pencil, pen or marker. REA-1 and REA-2 are supplied as individual labels in packages of 100. All others are supplied in rolls of 500 labels.

_____	REAGENT	_____
_____	ANALYSIS	_____
DATE	EXP. DATE	INIT.

REA-2 2" x 4 1/4" Red Border

LABORATORY REAGENT	
PREP. DATE	CHECKED
EXPIRATION DATE	_____
REAGENT	_____
CONC.	_____
PROCEDURE	_____
MADE BY	_____

UPCR-9500 2" x 4"

LABORATORY REAGENT	
DATE PREP.	MADE BY
DATE C'KED.	C'KED. BY
EXPIRATION DATE	_____

UPCR-9501 2" x 4"

Date Rec'd	_____
Date Opened	_____
Exp. Date	_____

UPCR-3015 3/4" x 1 1/2"

Date Prep.	_____
Exp. Date	_____
Made By	_____
Lot No.	_____

UPCR-3013  
Fl. Green 3/4" x 1 1/2"

REAGENT	
Date prep / Lot No.	_____
Date checked	Tech _____
Date in use	_____
Rec'd	Exp _____

UPCR-6048 1" x 2 1/4"

Date Rec'd	Initial
Date Checked	Initial
Date in Use	_____
Expires	_____

UPCR-6056 Fl. Orange 1" x 2 1/4"

Solution	_____
Concentration	_____
Procedure	_____
Preparation Date	_____
Expiration Date	_____
Prepared By	_____
Storage:	Room <input type="checkbox"/> Refrig. <input type="checkbox"/>

UPCR-7020 1" x 3"

DATE RECEIVED	_____
DATE OPENED	_____
DATE EXPIRES	_____


REA-3 3/4" x 1 1/2"

REAGENT	_____
ANALYSIS	_____
DATE	INIT.

REA-1 2" x 1"  
Red Border

REAGENT:	_____
TEST:	_____
DATE MADE:	DATE OUTDATED:
TESTED:	STORAGE:
COMMENTS:	_____

UPCR-9102 1 1/2" x 3" Fl. Chartreuse

	Reagent _____
	Prep. Date _____
	Exp. Date _____
	Storage _____
	Room <input type="checkbox"/> Refrig. <input type="checkbox"/>
	Initials _____

UPCR-31 1 1/2" x 3"

# Specialty Control Labels



## CALIBRATION LABELS

These colorful, eye-catching labels stick securely to difficult surfaces, draw immediate attention to important information, accept notations in pen, pencil or marker. Labels are supplied in rolls of 500.

**CALIBRATED**  
DATE \_\_\_\_\_

UPCR-1306  $\frac{5}{16}$ " x  $1\frac{1}{4}$ "

**NEW LOT  
DO NOT USE**

UPCR-3007  $\frac{3}{4}$ " x  $1\frac{1}{2}$ "  
Fl. Orange

**NEW LOT  
Record # in  
Log before using**

UPCR-3008  $\frac{3}{4}$ " x  $1\frac{1}{2}$ "

**This lot has been calibrated: No  
changes on ACA necessary.**

Before using these packs, please run  
and record controls, indicating new lot  
in notebook.

Date \_\_\_\_\_ Tech \_\_\_\_\_

UPCR-6039  $1$ " x  $2\frac{1}{4}$ "

This lot number has been  
calibrated. Changes must  
be made before use.  
Calib. Factor \_\_\_\_\_  
Reagent Rate \_\_\_\_\_  
Date \_\_\_\_\_

UPCR-3009  $\frac{3}{4}$ " x  $1\frac{1}{2}$ "

### CALIBRATION DATA

SETTING	ACTUAL

Date \_\_\_\_\_ By \_\_\_\_\_

UPCR-9051  $1\frac{1}{2}$ " x  $1\frac{1}{2}$ "

**VOLUMETRIC DEVICE  
CALIBRATED**

DATE \_\_\_\_\_

% DEVIATION \_\_\_\_\_ TECH \_\_\_\_\_

UPCR-5006  $\frac{3}{4}$ " x  $2$ "  
Fl. Chartreuse

**NOT  
CALIBRATED**  
**DO NOT  
USE**

Date  
Rec: \_\_\_\_\_

Init: \_\_\_\_\_

UPCR-9200  $1\frac{3}{4}$ " x  $2\frac{1}{2}$ " Fl. Red

**CALIBRATED AND  
PERFORMANCE TESTED BY**

Name \_\_\_\_\_ Date \_\_\_\_\_

UPCR-6036  $1$ " x  $2\frac{1}{4}$ " Fl. Green

**CALIBRATED  
NOT IN USE**

Before use  
Check SP and SF  
in Calibration Log

SP/C0 \_\_\_\_\_ SP/C1 \_\_\_\_\_

Date \_\_\_\_\_ Tech \_\_\_\_\_

UPCR-9050  $1\frac{1}{2}$ " x  $2\frac{1}{2}$ " Fl. Chartreuse

### CENTRIFUGE PERFORMANCE TEST

Setting	Cent. Tach.	Photo Tach.
		RPM
		RPM
		RPM
		RPM
		RPM
		RPM
		RPM
		RPM
		RPM
		RPM

MAX \_\_\_\_\_  
Timer \_\_\_\_\_ Stopwatch \_\_\_\_\_  
Date \_\_\_\_\_ By \_\_\_\_\_  
Due \_\_\_\_\_ Control # \_\_\_\_\_

UPCR-9101  $1\frac{1}{2}$ " x  $3$ "

**NOT  
CALIBRATED**  
**DO NOT  
USE**

Date  
Rec: \_\_\_\_\_

Init: \_\_\_\_\_

UPCR-6037  $1$ " x  $2\frac{1}{4}$ " Fl. Red

Test \_\_\_\_\_  
There are \_\_\_\_\_ bundles of  
lot # \_\_\_\_\_  
No. \_\_\_\_\_  
calibrated and ready to restore  
calibration. Please follow restore  
CAL procedure, green sticker and  
relocate all \_\_\_\_\_ bundles.  
No. \_\_\_\_\_

UPCR-9301  $1\frac{3}{4}$ " x  $3$ " Fl. Chartreuse

**READY  
FOR USE**

Carton opened \_\_\_\_\_ Date \_\_\_\_\_

Lot put in use \_\_\_\_\_ Date \_\_\_\_\_

Tech. \_\_\_\_\_

UPCR-9201  $1\frac{3}{4}$ " x  $2\frac{1}{2}$ " Fl. Green

**DEFECTIVE  
DO NOT USE**

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

DATE: \_\_\_\_\_ BY: \_\_\_\_\_

ME-17  $2\frac{1}{2}$ " x  $2\frac{1}{2}$ " Fl. Red Only

## LOT SIGNAL LABELS

Reduce chance of error in your laboratory with highly visible Shamrock Lot Signal Labels. Write on them with pen, pencil or rubber stamp. Their permanent adhesive sticks to any dry surface. Labels are supplied in rolls of 500 labels.

DATE RECEIVED \_\_\_\_\_  
**DO NOT USE  
THIS LOT  
NUMBER**

LL-142  $1\frac{1}{2}$ " x  $3$ " Fl. Red

**THIS LOT IS  
READY  
FOR USE**

DATE \_\_\_\_\_ INITIALS \_\_\_\_\_

LL-143  $1\frac{1}{2}$ " x  $3$ " Fl. Green

**THIS LOT HAS BEEN CALIBRATED  
BUT CHANGES MUST BE MADE ON  
THE INSTRUMENT BEFORE USE.**

Date \_\_\_\_\_ Initials \_\_\_\_\_

LL-144-A  $1\frac{1}{2}$ " x  $3$ "

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Tamper Resistant Labels and Tapes

Shamrock tamper resistant labels and tapes are ideal for use in special screening programs or whenever a specimen's validity might become suspect. Use them with specimen containers, envelopes, bags and boxes.

## SPECIMEN CONTAINER LABELS

Pressure-sensitive adhesive holds the labels permanently in place. Once a specimen container has been closed and sealed it cannot be reopened without tearing the label. STR-50 and STR-60 are supplied in packs of 500 individual labels, minimum order 2 packs.

**SEAL** Place SEAL over top of container before affixing label

LABEL Wrap around container, overlapping ends of seal strip

PATIENT NAME \_\_\_\_\_ ACCESSION # \_\_\_\_\_

DR/CLINIC NAME \_\_\_\_\_ DATE COLLECTED \_\_\_\_\_

COLLECTED BY \_\_\_\_\_ Signature \_\_\_\_\_

INSTRUCTION - Complete by person witness by collection

1. Place Accession # sticker on container.
2. Place SEAL strip over top of container & over one edge of Access. # sticker, leaving # signs.
3. Place LABEL around container, overlapping ends of SEAL. Don't cover Access. #

STR-50 1½" x 8"

DO NOT USE IF PROTECTIVE SEAL IS BROKEN.

DO NOT USE IF PROTECTIVE SEAL IS BROKEN.

STR-51 ½" x 4½"  
200 labels per roll



**SEAL** PLACE SEAL STRIP OVER TOP OF TUBE BEFORE AFFIXING LABEL

OVERLAP SEAL STRIP

MAILING CONTAINER LABEL  
BLOOD ALCOHOL SPECIMEN

AFFIX FIRST CLASS POSTAGE

NAME OF SUBJECT

**YOUR LABORATORY NAME AND ADDRESS**

WRAP AROUND TUBE

STR-60 3" x 6"

## DESTRUCTIBLE MATERIAL LABELS

Labels, featuring permanent adhesive, cannot be removed in one piece; portions will remain behind to signal tampering.

**INDIVIDUAL LABELS** are supplied in packages of 100. Minimum order 6 packages.

**SECURITY SEAL** DO NOT TAMPER

DATE \_\_\_\_\_

INITIALS \_\_\_\_\_

EL-11 1½" x 6½" Red

**LAB SEAL** DO NOT TAMPER

DATE \_\_\_\_\_

INITIALS \_\_\_\_\_

EL-12 1½" x 6½" Red

**EVIDENCE**

EL-13 1½" x 6½" Red



**LABELS IN CONTINUOUS ROLLS** measure 1-3/8" x 6-5/8" and may be ordered with a stock or custom imprint. Supplied 192 imprints per roll, minimum order 3 rolls of any one imprint.

**SECURITY SEAL** DO NOT TAMPER

DATE \_\_\_\_\_

INITIALS \_\_\_\_\_

EL-110-4 (Stock Imprint) 1½" x 108' Red

**YOUR MESSAGE HERE**

EL-100-6 (Custom Imprint) EL-100-4 1½" x 108' Red  
EL-100-6 1½" x 108' Blue

## "OPENED" LABELS

Labels feature your message clearly printed in black ink. When removal is attempted, labels leave a strong red or blue residue with the word "opened" repeated for unassailable evidence of tampering. Supplied in packages of 100 individual labels. Minimum order 6 packages.

**MESSAGE HERE**


EL-14-4 1½" x 6½" Red

**MESSAGE HERE**

EL-14-6 1½" x 6½" Blue

## NOTCHED, PRE-CUT LABELS

Bearing your imprint in red or black ink and backed with permanent adhesive, labels tear into pieces when removed, leaving fragments as an immediate sign of tampering. Supplied in rolls of 125 labels in dispenser boxes. Minimum order 12 rolls of each imprint.



# Specialty Consecutive Number Labels and Tapes



## TAPE LABELS

Consecutive Numbers are printed on versatile Shamrock Labeling Tape that sticks to any surface, will stand up to a wide range of temperatures. Unique adhesive allows tape to be transferred from one surface to another. Order Shamrock Tape Consecutive Number Labels printed black with your choice of prefix and/or suffix, on any one of 19 tape colors. Each number can be printed one to eight times, or 10, 12, 14 or 16 times as you specify.

Tape Consecutive Number Labels measure 1/2" x 1 1/4" and are supplied in rolls of either 500 or 2,000 labels. Unless otherwise specified, numbers will be printed black on white tape with "No." prefix.

SCR-10 (500 per roll)  
SCR-20 (2,000 per roll)

## PRE-CUT LABELS IN TABLETS OR ROLLS

Shamrock Consecutive Numbers are also available as round corner, pre-cut labels supplied in roll form or in tablet form. Each number may be printed as many times as needed. The 1/2" x 1 1/4" labels feature a unique, glove friendly, pressure-sensitive adhesive that has a special affinity for all silicone treated surfaces and hard to label shapes, holds the label securely in place without flaging.

Labels are available in 19 colors, printed black, with your choice of prefix and/or suffix. Labels are supplied printed black on white stock with "No." prefix unless otherwise specified. A minimum of 1,000 labels is required.

SCS-12 (Tablet)  
SCR-14 (Roll)

### 19 Specialty Colors

1 White	2 Yellow	3 Green	4 Red	5 Orange
6 Blue	7 Pink	8 Tan	9 Salmon	10 Copper
11 Rose	12 Violet	13 Lavender	14 Gray	15 Lime
16 Aqua	17 Chartreuse	18 Gold	19 Silver	

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM






## CNI SYSTEM FOR PATHOLOGY SPECIMENS

Shamrock's CNI (Consecutive Number Identification) System, provides three different labels, each imprinted in identical number sequence, each designed to function in a different step of tissue examination. Numbers may be printed once or repeated as many times as needed, with or without a prefix and/or suffix. You may order the complete CNI System or any single part. A 5,000 label minimum order for each system part is required.




### STEP ONE Numbering the Cassette

Cassette numbers are printed in permanent, black ink on sturdy card stock with no adhesive backing. The 1/4" x 1 1/2" label fits beside the tissue inside the cassette. The number remains unaltered throughout use in such solvents as Xylene, Benzene and alcohol, will not interfere with the tissue.



### STEP TWO Numbering the Paraffin Block and/or Specimen Storage Bottle

Labels measure 1/4" x 1" to fit the paraffin block exactly, or use them on specimen storage bottles, wherever you need to correlate the whole specimen with its sections. Pressure-sensitive adhesive holds them securely in place making any attempt at removal immediately evident.



### STEP THREE Numbering the Microscope Slide

The 15/16" square, pathology thick numbered label carries permanent, pressure-sensitive adhesive that holds it firmly to the slide throughout examination and storage. Your hospital name and department can be printed on each slide label. Blank, pathology-thick end labels are also available.

Mylar Overlay labels available. See page L-17.

## PIGGYBACK CONSECUTIVE NUMBER LABELS

Shamrock Piggyback Labels provide a set of individual pressure-sensitive labels imprinted with the same number, on a single sheet of adhesive-backed carrier paper. Labels measure 1/2" x 1 1/4", printed in black with a "No." prefix.

SCPB-30 (Piggyback labels supplied in tablets)  
SCPB-31 (Piggyback labels supplied in rolls)

### Piggyback Material



## SOLVENT RESISTANT CONSECUTIVE NUMBER LABELS

Labels hold fast during staining and storage; won't fade or wash away in Xylene, Shandon, Americlear, alcohol or acid solutions. Each number may be printed as often as needed.

- Labels measure 5/8" x 7/8".
- Xylene resistant numbers are available on clear or opaque white label stock. Please specify your choice of black, blue, green or red ink when ordering.
- Shandon and Americlear resistant numbers are supplied on shiny silver material, printed black.




SCR-15-SX Clear SCR-15-SX-OP Opaque SCR-15-AS Silver

# Specialty Consecutive Number Labels

## PRE-PRINTED BAR CODED CONSECUTIVE NUMBER SETS

Available in a variety of stock configurations that are supplied in your choice of piggyback or single label construction, in tablets or roll form. Labels can feature eye readable and bar code numbers alone or in combination with your facility name. You specify the number sequence, bar code language and any additional imprinting.




SCBAR-1 7/16" x 1 1/8"




SCBAR-08000 3/4" x 2"



SCBAR-20 1 1/4" x 1 7/8"



Roll Form



Tablet Form

Bar Codes can be printed at 600 DPI upon request. Call 800-323-0249.



SCBAR-01021 2 1/2" x 2"



SCBAR-01000 2 1/2" x 2"

SCBAR-01022 2 1/2" x 2"



SCBAR-40 2 1/2" x 4 5/8"




SCBAR-23 2 1/2" x 4 5/8"


SCBAR-09021 3" x 2 3/4"



SCBAR-22 3 1/4" x 1 1/4"



SCBAR-07082 4 1/2" x 2"



SCBAR-06021 4 1/2" x 2"

Shamrock can supply any stock or custom bar coded label with piggyback construction. Call 800-323-0249 for complete information.


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## BLANK LABELS FOR ON-DEMAND CONSECUTIVE NUMBERS

**Thermal Printer Labels** Catalog numbers designate Thermal Direct Material. For Thermal Transfer material add the suffix "-T" (THERM-504R-T). For premium Kimdura Transfer material add "-KT" (THERM-504R-KT). See pages L-40 thru L-59 for more Shamrock Thermal Labels.




**Laser Printer Labels** Shamrock can customize your label sheet. Call for complete information. See pages L-62 thru L-65 for more Shamrock Laser Labels.


Cat. No.	Size	Labels/Sheet	Sheets/Ctn
LASER-20	1/2" x 1"	160	32
LASER-21	1/2" x 1 1/4"	126	40
LASER-22	1/2" x 1 3/4"	80	63
LASER-1	1/2" x 2"	80	63
LASER-23	1 5/16" x 1 5/16"	88	57
MSL-100-FLAT	1" x 2"	44	114




LASER-20 1/2" x 1"




LASER-21 1/2" x 1 1/4"




LASER-22 1/2" x 1 3/4"



LASER-1 1/2" x 2"



LASER-23 1 5/16" x 1 5/16"



MSL-100-FLAT 1" x 2"

## Dot Matrix Printer Labels

See pages L-66 and L-67 for more Shamrock Pinfed Labels.



BDP-3 1/4" x 1"



SNO-137 7/16" x 1 1/2"



SNO-443 1 5/16" x 1 5/16"

Shamrock can Add Color or Color Highlights to Any Stock Blank Label.

### Most Popular Colors



Call Shamrock for complete information on color options available.


# Specialty Microscope Slide Labels

## STANDARD MICROSCOPE SLIDE LABELS AND END LABELS


Available in standard or pathology thickness, may be imprinted to your specifications. Call for information.

### Label Sheets 5 1/4" x 8"

#### Slide Labels 1,000 Labels/Pack




SLS-15 15/16" x 15/16" (Standard)  
SLSP-15 15/16" x 15/16" (Pathology)




SLS-15RC 15/16" x 15/16" (Standard)  
SLSP-15RC 15/16" x 15/16" (Pathology)

#### End Labels 2,000 Labels/Pack




SES-16 15/32" x 15/16" (Standard)  
SESP-16 15/32" x 15/16" (Pathology)




SES-16RC 15/32" x 15/16" (Standard)  
SESP-16RC 15/32" x 15/16" (Pathology)

### Label Rolls

#### Slide Labels 1,000 Labels/Roll




SLR-15 15/16" x 15/16" (Standard)  
SLRP-15 15/16" x 15/16" (Pathology)




SLR-15RC 15/16" x 15/16" (Standard)  
SLRP-15RC 15/16" x 15/16" (Pathology)

#### End Labels 1,000 Labels/Roll



SER-16 15/32" x 15/16" (Standard)  
SERP-16 15/32" x 15/16" (Pathology)



SER-16RC 15/32" x 15/16" (Standard)  
SERP-16RC 15/32" x 15/16" (Pathology)

Standard  
Thickness



Pathology  
Thickness




## LABELS FOR COMPUTER PRINTING

On-demand labels are available in standard thickness.

### Laser Label Sheets 8 1/2" x 11"


#### Slide Labels



SLS-LASER  
15/16" x 15/16"  
88 labels/sheet  
57 sheets/ctn

SLS-LASER-RC 15/16" x 15/16"  
80 labels/sheet  
63 sheets/ctn

#### End Labels



SES-LASER 15/32" x 15/16"  
168 labels/sheet  
30 sheets/ctn

### Thermal Printer Labels



7/8" x 7/8" 3" Core 12,800/roll  
THERM-1318R Thermal Direct  
THERM-1318R-T Thermal Transfer  
THERM-1318OPR Thermal Transfer  
(Opaque Polyester)




15/16" x 15/16" 3" Core 12,800/roll  
THERM-1213R Thermal Direct  
THERM-1213R-T Thermal Transfer  
THERM-1213OPR Thermal Transfer  
(Opaque Polyester)



THERM-1840R-DO  
1" x 1"  
1" core 750/roll



THERM-1839RSQDO  
15/16" x 15/16"  
750 per roll



THERM-1839R-DO  
15/16" x 15/16"  
1" core 750/roll

### Dot Matrix Printer Labels



BDP-1 15/16" x 15/16" Fan Folded 30,000/carton



BDP-5 15/16" x 15/16" Fan Folded 25,000/carton

SLIDE LABELS ARE AVAILABLE IN COLOR • CALL FOR COMPLETE INFORMATION

White

Yellow

Green

Red

Orange

Blue

Rose


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



### MYLAR OVERLAY FOR MICROSCOPE SLIDE LABELS






Protect handwritten or printer generated information from repeated cleanings and harsh solvents. Overlay labels are supplied 500 per roll.



## Specialty Labels for Photography

### METRIC MEASURE LABELS

Labels identify and accurately size specimens, feature removable adhesive that holds firmly throughout use then is easily removed, leaving no sticky residue. Supplied 1,000 per roll.




## Specialty Labels for Water Testing

### PRESERVATIVE IDENTIFICATION LABELS


#### Circle Labels

Identify container contents on cap and/or sides with eye-catching 1" dot labels. Permanent adhesive holds securely throughout use. Labels are supplied 1,000 per roll.



#### Tamper Evident Labels

Place identifying center circle of label on cap and smooth ends down each side of the container. Unique, clear material won't obstruct view of contents; features serrated edges that shred to indicate tampering if removal is attempted. Labels measure 1 3/8" x 9" and are supplied 100 labels per package.



Shamrock can also supply labels imprinted to meet your specific needs. Call for complete information.

# Specialty Animal Care Labels



## CAGE CARDS

Available in sturdy card stock or pressure-sensitive label material. Choose from stock imprints or customize card color, imprint and size. Contact Shamrock for ordering assistance.

Animal Name:	<input type="checkbox"/> Female <input type="checkbox"/> Male	Description:	Date in:
Tag No.			Date out:
Owner:	Phone No:	Feeding Instructions:	
Notes:			

SCC-35-1 3" x 5"  
SCC-46-1 4" x 6"

## UNDER VETERINARY OBSERVATION

Animal ID \_\_\_\_\_  
Investigator \_\_\_\_\_  
Protocol # \_\_\_\_\_  
From: \_\_\_\_\_ To: \_\_\_\_\_  
Date: \_\_\_\_\_ Date: \_\_\_\_\_  
Comments: \_\_\_\_\_

SCC-53-1 5" x 3"  
SCC-64-1 6" x 4"

## MEDICATION LABELS

Pre-cut labels feature permanent adhesive, may be marked with pen, pencil or maker.

### VIAL OPENED

EXPIRES  
DATE \_\_\_\_\_ INITIALS \_\_\_\_\_

UPCR-3402 3/4" x 1 1/2"  
Fl. Red 500/roll

### MEDICATION ADDED DRUG

AMOUNT ADDED BY  
DATE \_\_\_\_\_ TIME \_\_\_\_\_

SN-2001 1 1/2" x 2 1/2"  
1,000/roll

### MEDICATION ADDED

DRUG	ADDITIONAL	AMOUNT

ADDED BY \_\_\_\_\_  
DATE \_\_\_\_\_ TIME \_\_\_\_\_ A.M.  
P.M.  
PATIENT \_\_\_\_\_  
ROOM NO. \_\_\_\_\_  
SPECIAL INFORMATION \_\_\_\_\_

SMA-231 3" x 2"  
1,000/roll

## TAPE & PRE-CUT COMMUNICATION LABELS

**Tape Labels**, supplied in 500" rolls, stick to plastic, metal, glass, wood; resist water, animal excretions, formaldehyde and oily substances; won't leave messy residue when removed.


All stock Tape labels are also available as **Pre-Cut Labels** measuring 1" x 2 1/2" and supplied in 500" rolls. To order, simply add the suffix -PRE to the appropriate catalog number as illustrated.

### TAPE Label



SAC-1 3/4" x 500"

### PRE-CUT Label



SAC-1-PRE 1" x 2 1/2"



SAC-1 3/4" x 500"



SAC-2 3/4" x 500"



SAC-3 3/4" x 500"



SAC-4 1" x 500"




SAC-9 3/4" x 500"




SAC-10 3/4" x 500"



SAC-11 3/4" x 500"



SAC-17 3/4" x 500"



SAC-15 3/4" x 500"



SAC-16 3/4" x 500"




SAC-20 3/4" x 500"

## DRUG SYRINGE LABELS

Labels feature permanent adhesive that holds them securely in place throughout use. They are supplied on continuous tape in 1/2" x 500" rolls; as 1/2" x 1 1/2" pre-cut labels in rolls of 323 labels; or as 1/2" x 2" labels in packs of 250 labels.



SA-201 Yellow



SA-210 Orange



SA-2413 Blue




SA-324 White



SA-205 Yellow




SA-212 Orange



SA-243 Blue




SA-3221 White




SA-2416 Blue



SAV-900 Yellow



SAV-910 Orange



SAV-911 Orange




SAV-912 Orange



SAV-930 White/Red



SAV-940 Blue




SAV-941 Blue



SAV-950 White/Blue



SAV-1000 Red



SAV-1001 Red



SAV-1002 Red

Contact Shamrock for the complete list of stock drug syringe labels or information on labels designed and printed to your specifications.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



# Specialty Warning Labels and Tapes



## BIOHAZARD WARNING LABELS

All Shamrock Biohazard Labels comply with OSHA labeling standards for bloodborne pathogens as outlined in Section 1910.1030 of Title 29 of the Code of Federal Regulations. They feature brilliant fluorescent color as mandated by OSHA and permanent, pressure-sensitive adhesive to keep them securely in place as long as needed.



SBH-1 3" x 2 1/4" Roll of 500  
CAT. NO. SBH-1




SBH-4 1" x 3" Roll of 500  
CAT. NO. SBH-4



SBH-8 1" x 1" Roll of 500




SBH-8-D 1" Diam. Roll of 500




SBH-8-L 5/16" x 1 1/4" Roll of 500



SBH-5 3 1/2" x 2 1/4" Package of 100  
CAT. NO. SBH-5




SBH-3 8" x 10" Pressure-sensitive  
SBH-2 8" x 10" Non-adhesive card stock  
Package of 10  
CAT. NO. SBH-3




SBH-103 8" x 10" Pressure-sensitive  
SBH-104 8" x 10" Non-adhesive card stock  
Package of 10  
CAT. NO. SBH-103




SBH-101 8" x 10" Pressure-sensitive  
SBH-102 8" x 10" Non-adhesive card stock  
Package of 10  
CAT. NO. SBH-101



SBH-62 1 1/2" x 3" Roll of 500  
CAT. NO. SBH-62



SBH-60 2" x 4" Package of 100  
CAT. NO. SBH-60



SBH-61 2" x 4" Package of 100  
CAT. NO. SBH-61

**BIOHAZARD WARNING LABELS**



SBH-9 3" x 4" Package of 100




SBH-70 1" x 2 1/4" Roll of 500




CW-10 2 1/2" x 4" Roll of 500



SBH-73 1" x 2 1/4" Roll of 500



SBH-74 3" x 2 1/4" Roll of 500



UPCR-9800 3" x 2" Roll of 500



UPCR-9401 2" x 3" Roll of 500



UPCR-6027 1" x 2 1/4" Roll of 500



UPCR-9550 2" x 5" Roll of 500

**BIOHAZARD WARNING TAPE**

Pressure-sensitive tape sticks to all surfaces yet is easily removed from walls, doors and cabinets, leaving no messy adhesive residue to clean up.



SBH-48



SBH-48C



SBH-48T

Each roll of 500" tape is 2" wide. Label imprint area measures 2" x 6".

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Warning Labels and Tapes



## WASTE CONTROL LABELS

Labels feature aggressive permanent adhesive for superior hold.



**UNIVERSAL WASTE**

CONTENTS \_\_\_\_\_

ACCUMULATION START DATE \_\_\_\_\_

SHIPPER \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY, STATE, ZIP \_\_\_\_\_

WC-100 6" x 6" Roll of 250



**Infectious Waste**

Waste Generator

Name \_\_\_\_\_

Address \_\_\_\_\_

Phone (\_\_\_\_) \_\_\_\_\_

Water Transporter/Recipient

Name \_\_\_\_\_

Address \_\_\_\_\_

Phone (\_\_\_\_) \_\_\_\_\_

Total Wt. \_\_\_\_\_

UPCR-9901 3" x 5" Roll of 500



**INFECTIOUS WASTE**

BIOHAZARD

UPCR-9551 2" x 5" Roll of 500



**NON-HAZARDOUS WASTE**

GENERATOR INFORMATION: (optional)

SHIPPER \_\_\_\_\_


ADDRESS \_\_\_\_\_

CITY, STATE, ZIP \_\_\_\_\_

CONTENTS \_\_\_\_\_

SOLID WASTE EXCLUDED FROM REGULATORY PAPERS 49 CFR 261.4(b)

WC-101 6" x 6" Roll of 250



**HAZARDOUS MATERIALS**

**CHEMICAL WASTE**

HANDLE WITH CARE

DATE \_\_\_\_\_ DEPT. \_\_\_\_\_ BY \_\_\_\_\_

UPCR-9403 2" x 3" Roll of 500



**FOR INCINERATION ONLY**

UPCR-6066 1" x 2 1/4" Roll of 500



Clear, protective mylar overlay labels are also available. Contact Shamrock for complete information.

## WASTE LABELS OF YOUR OWN DESIGN

Labels specific to your regulatory and procedural needs. You choose the label material, size, colors and imprint.



**HAZARDOUS WASTE**

STATE AND FEDERAL LAWS PROHIBIT IMPROPER DISPOSAL. IF FOUND, CONTACT THE NEAREST POLICE OR PUBLIC SAFETY AUTHORITY, THE U.S. ENVIRONMENTAL PROTECTION AGENCY OR THE CALIFORNIA DEPARTMENT OF TOXIC SUBSTANCE CONTROL.

GENERATOR INFORMATION:

NAME \_\_\_\_\_ Community Medical Center

ADDRESS \_\_\_\_\_ 2468 Main Street PHONE (543) 864-5000

CITY \_\_\_\_\_ Home Town STATE \_\_\_\_\_ CA ZIP \_\_\_\_\_ 90182

EPA ID NO. \_\_\_\_\_ MANIFEST DOCUMENT NO. \_\_\_\_\_

EPA WASTE NO. \_\_\_\_\_ ACCUMULATION START DATE \_\_\_\_\_

CONTENTS, COMPOSITION: \_\_\_\_\_

PHYSICAL STATE:  SOLID  LIQUID | HAZARDOUS PROPERTIES:  FLAMMABLE  TOXIC  CORROSIVE  REACTIVE  OTHER


D.O.T. PROPER SHIPPING NAME AND UN OR NA NO. WITH PREFIX

**HANDLE WITH CARE!**

Call Shamrock for complete design and ordering assistance.

## RADIOACTIVE WARNING LABELS


Printed with vivid magenta ink on bright yellow background, labels can be used on any surface where visual warning is required.




**CAUTION RADIOACTIVE MATERIAL**

CAT. NO. SRA-12

SRA-12 4" x 6" Package of 10 labels



SRA-11 1/4" Diameter Roll of 1000 labels



UPCR-1301 5/16" x 1 1/4" Roll of 500 labels



UPCR-6033 1" x 2 1/4" Roll of 500 labels

## RADIOACTIVE WARNING TAPE

Supplied on tape in 500" rolls, specific information may be added to these labels in pen, pencil, felt marker or rubber stamp.




SRA-4 1" Tape




SRA-2 1" Tape



SRA-10 1/2" Tape



SRA-1 1" Tape



**CAUTION RADIOACTIVE MATERIAL**

ISOTOPE \_\_\_\_\_

SOURCE STRENGTH \_\_\_\_\_

SERIAL NO. \_\_\_\_\_

ASSAY DATE \_\_\_\_\_

INVESTIGATOR \_\_\_\_\_

SRA-6 2" Tape

# Specialty Labels for GHS Compliance



GHS is the UN Globally Harmonized System of identifying hazardous chemicals. Aligned with the OSHA Hazard Communication Standard 29 CFR 1910.1200, it defines specific labeling guidelines to facilitate easy to understand, international communication of health and safety information.

## PRE-PRINTED GHS HAZARD LABELS

Pre-Printed GHS Labels define, in words and pictograms, dangers related to a specific chemical and provide guidelines for prevention of and response to adverse events. Labels feature a tough, chemical and solvent resistant, opaque white material with permanent adhesive that sticks securely to plastic, glass and metal. Supplied 250 labels per roll.

### Xylene



#### DANGER

Flammable liquid and vapor. Causes skin irritation. May be fatal if swallowed and enters airways. May cause drowsiness or dizziness.

#### PREVENTION

Keep away from heat, sparks, and open flames. — No smoking. Keep container tightly closed.

Avoid breathing vapors or mist. Wash hands and any other contaminated skin thoroughly after handling. Wear protective gloves and eye protection. Use only outdoors or in a well-ventilated area.

#### RESPONSE

If swallowed: immediately call a poison center or doctor. Do NOT induce vomiting.

If inhaled: Remove person to fresh air and keep comfortable for breathing. Call a poison center or doctor if you feel unwell.

If on skin: Take off all contaminated clothing. Wash with plenty of soap and water or shower. Wash contaminated clothing before reuse. If skin irritation occurs: Get medical attention.

In case of fire: Use foam, water spray or fog. Dry chemical, carbon dioxide or sand may be used for small fires only. Do NOT use water in a jet.

#### STORAGE

Keep container tightly closed. Store locked up, in a cool, well-ventilated place.

#### DISPOSAL

Dispose of contents to an EPA permitted incinerator.

#### Vertical

GHS-1-V-32 3" x 2"

GHS-1-V-43 4" x 3"

GHS-1-V-64 6" x 4"

### Xylene

#### DANGER



Flammable liquid and vapor. Causes skin irritation. May be fatal if swallowed and enters airways. May cause drowsiness or dizziness.

#### PREVENTION



Keep away from heat, sparks, and open flames. — No smoking. Keep container tightly closed.



Avoid breathing vapors or mist. Wash hands and any other contaminated skin thoroughly after handling. Wear protective gloves and eye protection. Use only outdoors or in a well-ventilated area.

#### RESPONSE

If swallowed: immediately call a poison center or doctor. Do NOT induce vomiting. If inhaled: Remove person to fresh air and keep comfortable for breathing. Call a poison center or doctor if you feel unwell. If on skin: Take off all contaminated clothing. Wash with plenty of soap and water or shower. Wash contaminated clothing before reuse. If skin irritation occurs: Get medical attention.

In case of fire: Use foam, water spray or fog. Dry chemical, carbon dioxide or sand may be used for small fires only. Do NOT use water in a jet.

#### Horizontal

GHS-1-H-23 2" x 3"

GHS-1-H-34 3" x 4"

GHS-1-H-46 4" x 6"

Labels are available in 3 sizes for both horizontal and vertical formats.



## ADDITIONAL PRE-PRINTED GHS HAZARD LABELS

Shamrock offers a wide selection of pre-printed GHS Hazard Labels, in each of the sizes shown above, for both vertical and horizontal formats. Supplied 250 per roll, additional labels available from Shamrock include, but are not limited to:

Acetone

Concentrated Hydrochloric Acid

Ethanol

(10%) Formalin

Isopropanol (propan-2-ol)

Methanol

Sodium Hydroxide

Sulfuric Acid

Contact Shamrock for complete information and ordering assistance.


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## LABELS FOR ON-DEMAND THERMAL TRANSFER PRINTING

Shamrock supplies top quality, chemical resistant, opaque white labels in each of the sizes needed for on-site thermal transfer printing of GHS hazard warning labels. See pages L-40 thru L-59 for other Shamrock stock Thermal Transfer Labels. Labels in custom sizes and materials are available at your request.




THERM-1090R-KT 3" x 2" (3" core)  
THERM-1090RKT1C 3" x 2" (1" core)



THERM-2014R-KT 4" x 3" (3" core)  
THERM-2014RKT1C 4" x 3" (1" core)




THERM-405R-KT 6" x 4" (3" core)  
THERM-405RKT1C 6" x 4" (1" core)



THERM-1714R-KT 4" x 6" (3" core)  
THERM-1714RKT1C 4" x 6" (1" core)



THERM-506R-KT 2" x 3" (3" core)  
THERM-506RKT1C 2" x 3" (1" core)



THERM-320R-KT 3" x 4" (3" core)  
THERM-320RKT1C 3" x 4" (1" core)


### RIBBONS FOR THERMAL TRANSFER PRINTING

Shamrock offers chemical-resistant, black and red ribbons, ideal for printing our thermal transfer labels. To find the correct match for your printer, contact Shamrock with the printer make and model number, and the width of the ribbon required.



## PICTOGRAM LABEL SHEETS

Pictogram Labels are required on each warning label as an instant alert to the chemical hazard present. Shamrock combines the nine pictograms illustrated on 8½" x 11" sheets, hole-punched for handy binder storage. Supplied 25 sheets per package.



GHS-PIC-1 (1" x 1" labels)  
1 row of each pictogram

GHS-PIC-34 (¾" x ¾" labels)  
2 rows each of first 3 pictograms  
1 row each of remaining 6 pictograms

GHS-PIC-12 (½" x ½" labels)  
2 rows of each pictogram

**Call Shamrock for information on Pictogram Labels supplied in other formats.**

# Specialty Chemical Hazard Labels



## NFPA HAZARD LABELS

Warn emergency, rescue and fire personnel of potential hazards within a specified area.  
Labels are supplied 500 per roll unless otherwise indicated.

Reagent \_\_\_\_\_  
 Prep. Date \_\_\_\_\_  
 Exp. Date \_\_\_\_\_  
 Storage \_\_\_\_\_  
 Room  Refrig.   
 Initials \_\_\_\_\_

UPCR-31 1 1/2" x 3"

**HEALTH HAZARD**  
 4 Deadly  
 3 Extreme danger  
 2 Hazardous  
 1 Slightly hazardous  
 0 Normal material  
**SPECIFIC HAZARD**  
 ACID-Acid  
 ALK-Alkali  
 COR-Corrosive  
 OXY-Oxidizer  
 P-Polymerization  
 \*Radioactive  
 \*\*Use No Water  
**CHEMICAL NAME** \_\_\_\_\_  
**SDS #** \_\_\_\_\_

**FIRE HAZARD**  
 Flash Points  
 4 Below 73°F  
 3 Below 100°F  
 2 Above 100°F not exceeding 200°F  
 0 Will not burn

**REACTIVITY**  
 4 May detonate  
 3 Shock & heat  
 2 Violent chemical change  
 1 Unstable if heated  
 0 Stable

UPCR-36 2" x 3"

**H** HEALTH  OWNER  
**F** FLAMMABILITY  LOT NO.  
**R** REACTIVITY  REC'D.  
 PERSONAL PROTECTION  EXP.

UPCR-30 2" x 2"

**ETHYL ALCOHOL**


UPCR-43 1" x 3"

CAT. NO. UPCR-7002

UPCR-7002 1" x 3"

Chemical \_\_\_\_\_  
 Owner \_\_\_\_\_ Date \_\_\_\_\_  
 CAT. NO. UPCR-7003

UPCR-7003 1" x 3"



SBH-32 1" x 1"  
Roll of 1,000 labels

**CAUTION. CONTAINS FORMALDEHYDE**  
 Toxic by inhalation and if swallowed, irritating to the eyes, respiratory system, and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. May cause cancer. Repeated or prolonged exposure increases the risk.

UPCR-32 1" x 3"

**CAUTION**  
**CONTAINS FORMALDEHYDE**  
 TOXIC BY INHALATION AND IF SWALLOWED, IRRITATING TO THE EYES, RESPIRATORY SYSTEM AND SKIN. MAY CAUSE SENSITIZATION BY INHALATION OR SKIN CONTACT. RISK OF SERIOUS DAMAGE TO EYES. MAY CAUSE CANCER. REPEATED OR PROLONGED EXPOSURE INCREASES THE RISK.

UPCR-33 1" x 2 1/4"

**2.5% GLUTARALDEHYDE EM FIXATIVE**  
 MADE \_\_\_\_\_  
 EXPIRE \_\_\_\_\_  
 TECH \_\_\_\_\_

UPCR-34 1" x 2 1/4"

**CHEMICAL NAME** \_\_\_\_\_  
**MSDS#** \_\_\_\_\_ **DATE** \_\_\_\_\_

**HEALTH HAZARD**   
*(Possible injury)*

**FLAMMABILITY**   
*(Susceptibility to burning)*

**REACTIVITY**   
*(Susceptibility to release energy)*

**PERSONAL PROTECTION**  
 (Check all Protective Equipment that apply)

Safety Glasses  Apron  
 Face Shield  Full Suit  
 Splash Goggles  Boots  
 Vapor Respirator  Dust Respirator  
 Gloves  Other \_\_\_\_\_

UPCR-37 2 1/4" x 2 1/4"

UPCR-38 4" x 4"

**CHEMICAL NAME** \_\_\_\_\_  
**HAZARD**

**HEALTH HAZARD**  **PERSONAL PROTECTION**   
 4 Deadly Flash Points  Safety Glasses  Apron  
 3 Extreme danger 4 Below 73°F  Face Shield  Full Suit  
 2 Hazardous 3 Below 100°F  Splash Goggles  Boots  
 1 Slightly hazardous 2 Above 100°F not exceeding 200°F  Vapor Respirator  Dust Respirator  
 0 Normal material 1 Above 200°F  Gloves  Other  
 0 Will not burn

**ROUTE OF ENTRY**  
 Inhalation  Skin Absorption  
 Ingestion  Skin or Eye Contact

**TARGET ORGAN AFFECTED**  
 Liver  Teratogen  
 Kidneys  Skin  
 Blood  Eyes  
 Lungs  Central Nervous System  
 Mutagen  Mucous Membrane  
 Other \_\_\_\_\_

**SPECIFIC HAZARD** **REACTIVITY**  
 ACID-Acid 4 May detonate  
 ALK-Alkali 3 Shock & heat  
 COR-Corrosive 2 Violent chemical change  
 OXY-Oxidizer 1 Unstable if heated  
 P-Polymerization 0 Stable  
 \*Radioactive  
 \*\*Use No Water

UPCR-39 4" x 4"

Chemical Name \_\_\_\_\_  
 Common Name \_\_\_\_\_  
 Manufacturer \_\_\_\_\_  
 SDS# \_\_\_\_\_ Date \_\_\_\_\_

UPCR-40 1" x 3"

**ACETIC ACID**

UPCR-41 1" x 3"

**HYDROCHLORIC ACID**

UPCR-44 1" x 3"

**ISOPROPYL ALCOHOL**

UPCR-45 1" x 3"

**METHYL ALCOHOL**

UPCR-46 1" x 3"


**XYLENE**

UPCR-50 1" x 3"

**EMERGENCY GUIDE—HAZARD SIGNALS**

<b>HEALTH HAZARD*</b> 4—DANGER: May be fatal on short exposure. Specialized protective equipment required. 3—WARNING: Corrosive or toxic. Avoid skin contact or inhalation. 2—WARNING: May be harmful if inhaled or absorbed. 1—CAUTION: May cause irritation. 0—No unusual hazard. *Health hazard describes short term contact or inhalation hazard only.	<b>FIRE HAZARD</b> 4—DANGER: Flammable gas or extremely flammable liquid. 3—WARNING: Flammable liquid Flash point below 100. 2—CAUTION: Combustible liquid Flash point of 100°F to 200°F. 1—Combustible if heated. 0—Not combustible.	<b>INSTABILITY HAZARD</b> 4—DANGER: Explosive material at room temperature. 3—DANGER: May be explosive if shocked, heated under confinement, or mixed with water. 2—WARNING: Unstable, or may react if mixed with water. 1—CAUTION: May react if heated or mixed with water. 0—Stable. Not reactive when mixed with water.
--	--	---

SBH-30 2 1/2" x 4" Roll of 1,000 labels



SBH-31 2" x 2"  
Roll of 1,000 labels

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM





## CAUTION LABELS

Crisp, bold printing makes these labels difficult to overlook. They stick tightly to plastic, glass, metal or wood and stay in place as long as needed. Labels are supplied 500 per roll.

**10% FORMALIN**

CAT. NO. SB-374H

SB-374H 1" x 2 1/4"


**10% FORMALIN**  
CAUTION CONTAINS FORMALDEHYDE

TOXIC BY INHALATION AND IF SWALLOWED, IRRITATING TO THE EYES, RESPIRATORY SYSTEM, AND SKIN. MAY CAUSE SENSITIZATION BY INHALATION OR SKIN CONTACT. RISK OF SERIOUS DAMAGE TO EYES. MAY CAUSE CANCER. REPEATED OR PROLONGED EXPOSURE INCREASES THE RISK.

CAT. NO. SB-381H

SB-381H 1" x 2"  
SB-380H 1 3/4" x 2 1/2"

**10% FORMALIN POISON**



SB-382H 1" x 2 1/2"

**CAUTION:** Contains FORMALDEHYDE. Toxic by inhalation and if swallowed. Irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. Potential cancer hazard. Repeated or prolonged exposure increases the risk.

UPCR-5008 3/4" x 2"

**CAUTION**  
CONTAINS FORMALDEHYDE

TOXIC BY INHALATION AND IF SWALLOWED, IRRITATING TO THE EYES, RESPIRATORY SYSTEM, AND SKIN. MAY CAUSE SENSITIZATION BY INHALATION OR SKIN CONTACT. RISK OF SERIOUS DAMAGE TO EYES. MAY CAUSE CANCER. REPEATED OR PROLONGED EXPOSURE INCREASES THE RISK.

CAT. NO. SB-385H

SB-385H 1" x 2"  
SB-386H 1 3/4" x 2 1/2"

**CAUTION**  
RELEASES FORMALDEHYDE


TOXIC BY INHALATION. IRRITATING TO THE EYES, RESPIRATORY SYSTEM AND SKIN. MAY CAUSE SENSITIZATION BY INHALATION OR SKIN CONTACT. MAY CAUSE CANCER. REPEATED OR PROLONGED EXPOSURE INCREASES THE RISK.

CAT. NO. SB-387H

**DO NOT DRINK**

UPCR-5009 3/4" x 2" Fl. Red

**CAUTION: CONTAINS FORMALDEHYDE.**  
Toxic by inhalation and if swallowed. Irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. May cause cancer. Repeated or prolonged exposure increases the risk.



BIOHAZARD


UPCR-6042 1" x 2 1/4" White with Fl. Red

**10% FORMALIN**  
CAUTION: CONTAINS FORMALDEHYDE

Toxic by inhalation and if swallowed. Irritating to the eyes, respiratory system and skin. May cause sensitization by inhalation or skin contact. Risk of serious damage to eyes. Potential of Repeated or prolonged exposure increases

SB-387H 1 3/4" x 2 1/2"

**10% NEUTRAL BUFFERED FORMALIN**  
CONTAINS FORMALDEHYDE

DANGER: May be fatal if swallowed. Harmful if inhaled or absorbed through skin. Causes irritation to skin, eyes and respiratory tract. Sensitizer. Potential carcinogen.  **POISON**


SB-500H-T 3/4" x 60 yd tape roll

**POISON**

UPCR-1300 5/16" x 1 1/4"  
Fl. Red

UPCR-6046 1" x 2 1/4" Fl. Red

**POISON**




UPCR-6034 1" x 2 1/4" Fl. Red

**CAUTION**  
**POISON**

CAT. NO. SB-373H

SB-373H 1" x 2 1/4" Fl. Red

**10% FORMALIN**

UPCR-1318 5/16" x 1 1/4"

**10% FORMALIN POISON**

UPCR-1319 5/16" x 1 1/4"  
Fl. Red

**XYLENE**  
WARNING: FLAMMABLE LIQUID

CAUTION: CONTAINS ETHYL BENZENE, TOXIC BY INHALATION AND IF SWALLOWED, IRRITATING TO EYES, RESPIRATORY SYSTEM, AND SKIN. WASH THOROUGHLY AFTER HANDLING.

CAT. NO. SB-388H

SB-388H 1" x 2"

**ALCOHOL**  
WARNING: FLAMMABLE LIQUID

CAUTION: CONTAINS ETHYL, METHYL AND/OR ISO-PROPYL ALCOHOL. TOXIC BY INHALATION AND IF SWALLOWED. WASH HANDS THOROUGHLY AFTER HANDLING.

CAT. NO. SB-389H

SB-389H 1" x 2"

**95% ALCOHOL FOR PAP SMEARS**

UPCR-5007 3/4" x 2"

**GLUTARALDEHYDE**

UPCR-7017 1" x 3" Yellow

**CYTOLOGY FIXATIVE**  
**POISON**





UPCR-6043 1" x 2 1/4" Fl. Red

**FIXATIVE**




UPCR-6044 1" x 2 1/4" Fl. Red

**FIXATIVE FOR SPUTUM CYTOLOGY**  
**POISON**

UPCR-6045 1" x 2 1/4" Fl. Red

**PAP CYTOLOGY FIXATIVE**

PATIENT \_\_\_\_\_  
DOCTOR \_\_\_\_\_  
SOURCE \_\_\_\_\_  
DATE \_\_\_\_\_ HOSP. # \_\_\_\_\_

UPCR-6047 1" x 2 1/4"

# Specialty Caution Labels and Tapes

 CAT. NO. SB-370H SB-370H 1" x 2¼" Fl. Green	 UPCR-6029 1" x 2¼" Fl. Orange	 UPCR-6028 1" x 2¼" Fl. Orange	 CAT. NO. SB-372H SB-372H 1" x 2¼" Fl. Red
 UPCR-6030 1" x 2¼" Fl. Red	 CAT. NO. SB-371H SB-371H 1" x 2¼" Fl. Red	 CAT. NO. SB-375H SB-375H 1" x 2¼" Fl. Red	 UPCR-6032 1" x 2¼" Fl. Red
 SL-7 1" x 500" tape roll	 SL-8 1" x 500" tape roll	 UPCR-7008 1" x 3" Fl. Red	 UPCR-7007 1" x 3"

# Specialty Shipping Labels and Transport Bags

## D.O.T./IATA LABELS

Compliant with DOT rules for shipping and IATA Packing Instruction 650. Supplied 500 per roll unless noted.



HAZ-120 7/8" x 7/8"  
 HAZ-220 4" x 4" 250/pkg




HAZ-121 7/8" x 7/8"  
 HAZ-221 4" x 4" 250/pkg



HAZ-400 2½" x 2½"  
 250/pkg



HAZ-403 4" x 3"  
 250/pkg



HAZ-319 6" x 6"  
 250/pkg



UPCR-6132 1" x 2¼"  
 (QTY)




UPCR-9401 2" x 3" Fl. Red

## RECLOSABLE POLY BAGS


Available in a variety of sizes featuring 2 mil or 4 mil material and 2 wall construction or 3 wall construction to provide a handy back pouch for paperwork.

### Ziplock Transport Bags




3 Wall/1,000 per ctn  
 SPPB-B22 (6" x 9") Printed  
 SPPB-B18 (6" x 9") Clear

### Write on Bags




Write-On Block  
 2 Wall/1,000 per ctn  
 SPPB-B12 (3" x 4")

### Absorb Bag



Absorb Pad  
 3 Wall/500 per ctn  
 SPPB-B100 (6" x 9")

## Adhesive Closure Transport Bags



3 Wall/500 per ctn  
 SPPB-G120 (15" x 17")

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Laboratory Communication Labels



These top quality labels stay where you put them, improve communications among hospital personnel and help prevent potentially dangerous and costly errors. Colorful labels draw attention to important information or instructions. **Labels will not stick to latex gloves** yet adhere tightly to siliconized or other hard to label surfaces. Shamrock offers the variety of Communication Laboratory Labels illustrated on the following pages and many more. Please call for information. Labels are supplied 500 per roll unless otherwise specified.

## Clerical Labels

<b>ABNORMAL</b>	<b>DISCARD</b> _____	<b>NAME ALERT</b>	<b>PRELIMINARY REPORT</b>	<b>SHORT DATED USE FIRST</b>
UPCR-1000 Fl. Red	UPCR-1011 Blue	UPCR-1020 Fl. Red	UPCR-1033 Yellow	UPCR-1090 Fl. Red
<b>Abnormal Lab Results Inform Physician</b>	<b>DO NOT REFRIGERATE</b>	<b>NOT DRAWN BY LAB</b>	<b>PRIORITY</b>	<b>SPECIAL HANDLING</b>
UPCR-1001 Fl. Red	UPCR-1012 Fl. Green	UPCR-1021 Fl. Pink	UPCR-1086 Fl. Red	UPCR-1043 Fl. Pink
<b>ABNORMAL Results verified by repeat determination.</b>	<b>DR. PHONED DATE BY</b>	<b>O.R. STAT</b>	<b>PRIORITY 1</b>	<b>STAT</b>
UPCR-1075 Fl. Red	UPCR-1013 Fl. Red	UPCR-1022 Fl. Red	UPCR-1087 Fl. Green	UPCR-1044 Fl. Red
<b>ABNORMAL VALUE</b>	<b>DRAW STAT</b>	<b>OUT PATIENT</b>	<b>PRIORITY 2</b>	<b>STAT Do _____ Draw _____</b>
UPCR-1076 Fl. Red	UPCR-1014 Fl. Pink	UPCR-1023 Fl. Red	UPCR-1088 Fl. Pink	UPCR-1045 Fl. Orange
<b>AM ADMIT PRE-OP</b>	<b>EMERGENCY</b>	<b>PANIC VALUES EXCEEDED</b>	<b>REFRIGERATE</b>	<b>Store at Room Temperature</b>
UPCR-1002 Fl. Orange	UPCR-1080 Fl. Red	UPCR-1024 Fl. Red	UPCR-1034 Fl. Chartreuse	UPCR-1091 Blue
<b>ASAP</b>	<b>E.R. — STAT</b>	<b>PATHOLOGIST REVIEW NEEDED</b>	<b>REORDERED Date _____ By _____</b>	<b>SURGERY _____ AM _____ PM</b>
UPCR-1077 Green	UPCR-1015 Fl. Red	UPCR-1320 Fl. Orange	UPCR-1035 Blue	UPCR-1046 Fl. Pink
<b>ATTENTION</b>	<b>ER STAT CALL REPORT</b>	<b>PATIENT TO PICK UP RESULTS</b>	<b>REPORT FAXED</b>	<b>TIMED SPECIMEN</b>
UPCR-1003 Fl. Orange	UPCR-1016 Fl. Red	UPCR-1025 Green	UPCR-1089 Fl. Pink	UPCR-1047 Fl. Red
<b>BLOOD/BODY FLUID PRECAUTIONS</b>	<b>EXTRA TUBE</b>	<b>PEAK LEVELS</b>	<b>REPORT PHONED</b>	<b>TIME TO REORDER</b>
UPCR-1004 Fl. Pink	UPCR-1081 Fl. Chartreuse	UPCR-1026 Gray	UPCR-1036 Fl. Pink	UPCR-1048 Fl. Chartreuse
<b>CALL REPORT TO DOCTOR</b>	<b>FAX RESULTS</b>	<b>PHONE REPORT</b>	<b>RESULTS VERIFIED</b>	<b>TO BE REPEATED</b>
UPCR-1005 Fl. Red	UPCR-1082 Fl. Chartreuse	UPCR-1027 Yellow	UPCR-1037 Blue	UPCR-1049 Fl. Green
<b>CALL RESULTS TO O.R.</b>	<b>FINAL REPORT</b>	<b>PHONE REPORT TO FLOOR</b>	<b>RUSH</b>	<b>TO FOLLOW</b>
UPCR-1006 Fl. Red	UPCR-1017 Fl. Red	UPCR-1028 Gray	UPCR-1038 Green	UPCR-1050 Blue
<b>Collecting Time Not Indicated</b>	<b>FROZEN</b>	<b>PRE ADMISSION LAB WORK</b>	<b>SAVE</b>	<b>TODAY</b>
UPCR-1007 Fl. Red	UPCR-1083 Fl. Red	UPCR-1029 Fl. Chartreuse	UPCR-1039 Fl. Pink	UPCR-1051 Green
<b>CONTAINER NOT LABELED</b>	<b>H/A PRECAUTIONS</b>	<b>PRE-OP</b>	<b>SENT TO REFERRAL LAB</b>	<b>TROUGH LEVELS</b>
UPCR-1008 Fl. Red	UPCR-1018 Fl. Pink	UPCR-1030 Fl. Red	UPCR-1040 Green	UPCR-1052 Tan
<b>CORRECTED REPORT</b>	<b>HOLD</b>	<b>PRE-SURGICAL LAB WORK</b>	<b>SERUM</b>	<b>URGENT</b>
UPCR-1009 Yellow	UPCR-1084 Fl. Orange	UPCR-1031 Fl. Pink	UPCR-1041 Fl. Red	UPCR-1053 Fl. Red
<b>DAILY LAB</b>	<b>LABORATORY REPORTS</b>	<b>PRECAUTION</b>	<b>SHARED SPECIMEN</b>	<b>USE FIRST</b>
UPCR-1010 White	UPCR-1019 Fl. Chartreuse	UPCR-1032 Fl. Chartreuse	UPCR-1042 Fl. Pink	UPCR-1054 Fl. Chartreuse

UPCR-1000 series labels are 5/16" x 1 1/4"

# Specialty Laboratory Communication Labels

Clerical Labels

**ABNORMAL LAB TEST**

**Call Results To  
OR/RR/ER/ \_\_\_\_\_**

**STAT**

UPCR-3000 3/4" x 1 1/2" Fl. Red

UPCR-3001 3/4" x 1 1/2" Fl. Red

**REPORT PHONED**

**RESULTS SIGNIFICANT  
INFORM PHYSICIAN  
AT ONCE**

UPCR-3002 3/4" x 1 1/2" Fl. Chartreuse

UPCR-3003 3/4" x 1 1/2" Fl. Red

**PANIC VALUE - CALL  
RESULTS TO DR.**

UPCR-3004 3/4" x 1 1/2" Fl. Red

UPCR-3005 3/4" x 1 1/2" Fl. Red

**ADMISSION  
LABORATORY TESTS  
HAVE BEEN COLLECTED**

**CALLED**

**DOCTOR**

UPCR-6000 1" x 2 1/4" Fl. Orange

UPCR-6001 1" x 2 1/4" Blue

UPCR-6002 1" x 2 1/4" Fl. Chartreuse

**DO NOT THROW AWAY  
HOLD SPECIMEN  
UNTIL \_\_\_\_\_**

**PATIENT REFUSES**

**PRESURGICAL  
LAB WORK COMPLETED**

UPCR-6003 1" x 2 1/4" Fl. Orange

UPCR-6004 1" x 2 1/4" Fl. Pink

UPCR-6005 1" x 2 1/4" Fl. Red

**SPECIMEN IMPROPERLY  
COLLECTED.  
FOLLOW INSTRUCTIONS,  
RE-SUBMIT.**

**Specimen Unacceptable  
Tube/Container  
Not Labeled**

**CONFIDENTIAL  
INFORMATION**

UPCR-6006 1" x 2 1/4" Fl. Red

UPCR-6007 1" x 2 1/4" Fl. Orange

UPCR-6023 1" x 2 1/4" Fl. Red

**CRITICAL VALUE**

**PRIORITY  
HANDLING**

**STAT**

UPCR-6024 1" x 2 1/4" Fl. Red

UPCR-6025 1" x 2 1/4" Fl. Red

UPCR-6026 1" x 2 1/4" Fl. Red

**CRITICAL TEST**

**CAUTION  
FOOD MUST NOT  
BE STORED IN  
THIS REFRIGERATOR**

**PRIORITY  
SPECIMEN**

UPCR-7045 1" x 3" Fl. Red

UPCR-9100 1 1/2" x 3" Fl. Pink

UPCR-5002 3/4" x 2" Fl. Green

**PLEASE RECOLLECT:**  
 Specimen improperly labeled  
 QNS for analysis  
 Specimen improperly collected  
 Remarks: \_\_\_\_\_  
 DATE \_\_\_\_\_ BY \_\_\_\_\_

**DIFFICULT VENIPUNCTURE  
BLOOD SAMPLE NOT COLLECTED**

**STORE AT ROOM  
TEMPERATURE**

UPCR-7004 1" x 3" Fl. Chartreuse

UPCR-9300 1 3/4" x 3" Fl. Chartreuse

UPCR-5003 3/4" x 2" Blue

**REFRIGERATE  
ON ARRIVAL**

**THE \_\_\_\_\_ TEST HAS  
BEEN SENT TO A REFERENCE  
LAB ON \_\_\_\_\_  
IF RESULTS ARE NOT ON CHART  
IN \_\_\_\_\_ DAY(S), PLEASE NOTIFY  
THE LABORATORY.**

UPCR-7005 1" x 3" Fl. Chartreuse

UPCR-6008 1" x 2 1/4" Fl. Chartreuse

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



## Department Designation Labels

<b>BLOOD BANK</b> UPCR-1175 Fl. Red	<b>CHEMISTRY</b> UPCR-1176 Green	<b>HEMATOLOGY</b> UPCR-1178 Fl. Red	<b>MICROBIOLOGY</b> UPCR-1180 Fl. Orange	<b>STAT LAB</b> UPCR-1182 Fl. Red
	<b>CYTOLOGY</b> UPCR-1177 Fl. Pink	<b>IMMUNOLOGY</b> UPCR-1179 Fl. Chartreuse	<b>SEROLOGY</b> UPCR-1181 Gray	UPCR-1000 series labels are 5/16" x 1 1/4"

## Chemistry Labels

<b>ACID PHOSPHATASE</b> UPCR-1183 Fl. Pink	<b>CPK</b> UPCR-1193 Pink	<b>HDL</b> UPCR-1200 Fl. Chartreuse	<b>Protein Electrophoresis</b> UPCR-1207 Yellow	<b>T3</b> UPCR-1310 Blue
<b>BUN</b> UPCR-1184 Green	<b>CPK ISOENZYMES</b> UPCR-1194 Pink	<b>HGB AIC</b> UPCR-1201 Fl. Red	<b>RPR</b> UPCR-1208 Fl. Pink	<b>T4</b> UPCR-1216 Blue
<b>CARDIAC PANEL</b> UPCR-1185 Pink	<b>CREATININE</b> UPCR-1195 Green	<b>IRON IBC</b> UPCR-1202 Gray	<b>SAVE FOR ISOENZYMES</b> UPCR-1209 Fl. Pink	<b>THEOPHYLLINE</b> UPCR-1217 Fl. Green
<b>CARDIAC PROFILE</b> UPCR-1186 Pink	<b>DIGOXIN</b> UPCR-1196 Blue	<b>LDH ISOENZYMES</b> UPCR-1203 Pink	<b>SODIUM/POTASSIUM</b> UPCR-1210 Blue	<b>THYROID PANEL</b> UPCR-1218 Blue
<b>CEA</b> UPCR-1188 White	<b>DRUG SCREEN</b> UPCR-1197 Yellow	<b>LIPID PROFILE</b> UPCR-1204 Tan	<b>SMA 23</b> UPCR-1213 Green	<b>TSH</b> UPCR-1219 Fl. Orange
<b>CO<sub>2</sub></b> UPCR-1191 Green	<b>ELECTROLYTES</b> UPCR-1198 Gray	<b>PREGNANCY TEST</b> UPCR-1205 Green	<b>SPEC. HEMOLYZED</b> UPCR-1214 Fl. Red	<b>URINE CHEMISTRIES</b> UPCR-1220 Blue
<b>CORTISOL</b> UPCR-1192 Fl. Orange	<b>GLUCOSE</b> UPCR-1199 Fl. Red	<b>PRENATAL SCREEN</b> UPCR-1206 Blue	<b>T3 &amp; T4</b> UPCR-1215 Blue	UPCR-1000 series labels are 5/16" x 1 1/4"

## Hematology Labels

<b>CBC</b> UPCR-1221 Fl. Red	<b>HGB.—HCT.</b> UPCR-1224 Fl. Red	<b>PLATELETS</b> UPCR-1226 Gray	<b>PSA</b> UPCR-1309 Pink	<b>SED. RATE</b> UPCR-1229 Fl. Red
<b>COAGULATION</b> UPCR-1222 Fl. Orange	<b>PTT</b> UPCR-1225 Fl. Chartreuse	<b>PROTHROMBIN TIME</b> UPCR-1227 Blue	<b>RETICULOCYTES</b> UPCR-1228 White	<b>SICKLE CELLS</b> UPCR-1230 Pink
<b>ESR</b> UPCR-1223 Fl. Red	UPCR-1000 series labels are 5/16" x 1 1/4"			

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Laboratory Communication Labels

## Microbiology/Serology Labels

<b>ANAEROBIC</b> UPCR-1231 Fl. Orange	<b>CULTURE</b> UPCR-1232 Blue	<b>CULTURE &amp; SENSITIVITY</b> UPCR-1233 Fl. Red	<b>POSITIVE CULTURE</b> UPCR-1234 Fl. Red	<b>RETURN SPECIMEN TO URINALYSIS</b> UPCR-1235 Yellow	<b>SAVE FOR POSSIBLE CULTURE</b> UPCR-1236 Fl. Pink
--	----------------------------------	---	--	--	--

UPCR-1000 series labels are 5/16" x 1 1/4"

<b>BLOOD CULTURE</b> NAME _____ Hosp. No. _____ Room No. _____ Date _____ Time _____ Drawn By _____ UPCR-6016 1" x 2 1/4" White	<b>SAVE FOR POSSIBLE CULTURE</b> UPCR-6018 1" x 2 1/4" Pink	<b>HOLD FOR CULTURE &amp; SENSITIVITY</b> UPCR-6051 1" x 2 1/4" Fl. Red
<b>CULTURE &amp; SENSITIVITY TO FOLLOW</b> UPCR-6017 1" x 2 1/4" Fl. Pink	<b>DOCTOR</b> NOTE: IF SENSITIVITY IS DESIRED ON THIS CULTURE PLEASE NOTIFY LAB WITHIN 24 HOURS. UPCR-6050 1" x 2 1/4" Green	<b>URINE FOR ROUTINE ANALYSIS AND CULTURE</b> UPCR-6053 1" x 2 1/4" Fl. Chartreuse

## Phlebotomy/Specimen Receiving Labels

<b>BLOOD DRAWN</b> UPCR-1055 Fl. Chartreuse	<b>DRAW BLOOD</b> AT _____ AM _____ PM _____ UPCR-1059 Fl. Red	<b>FREEZE SPECIMEN</b> UPCR-1063 Blue	<b>KEEP FROZEN</b> UPCR-1067 Fl. Red	<b>PRIORITY SPECIMEN</b> UPCR-1071 Fl. Chartreuse
<b>Blood Drawn From Art Line</b> UPCR-1056 Fl. Red	<b>EDTA</b> UPCR-1060 Pink	<b>FROZEN SPECIMEN</b> UPCR-1064 Fl. Red	<b>LINE DRAW</b> UPCR-1068 Fl. Pink	<b>SHORT SAMPLE</b> UPCR-1072 Fl. Red
<b>BLOOD DRAWN FROM BELOW IV</b> UPCR-1057 Fl. Pink	<b>FINGER STICK ONLY</b> UPCR-1061 Fl. Red	<b>HEPARIN</b> UPCR-1065 Green	<b>PEDIATRIC SPECIMEN</b> UPCR-1069 Fl. Chartreuse	<b>Spec. Icteric</b> UPCR-1073 Fl. Chartreuse
<b>CAPILLARY SAMPLE</b> UPCR-1058 Blue	<b>FREEZE SERUM</b> UPCR-1062 Fl. Green	<b>Hold for Possible Additional Testing</b> UPCR-1066 Fl. Pink	<b>PLASMA</b> UPCR-1070 Fl. Chartreuse	<b>Spec. Lipemic</b> UPCR-1074 Pink

UPCR-1000 series labels are 5/16" x 1 1/4"

<b>FROZEN SPECIMEN</b> <b>DELIVER IMMEDIATELY</b> UPCR-6009 1" x 2 1/4" Fl. Red	<b>SPECIMEN IS:</b> <input type="checkbox"/> HEMOLYZED <input type="checkbox"/> LIPEMIC <input type="checkbox"/> ICTERIC <b>THIS MAY CAUSE FALSE ELEVATIONS OF SOME RESULTS.</b> UPCR-6011 1" x 2 1/4" Fl. Chartreuse	<b>DIFFICULT VENIPUNCTURE</b> It is extremely difficult to obtain blood specimens from this patient. Please make every effort to coordinate laboratory work to minimize the number of venipunctures. Patient's Name _____ Room # _____ UPCR-7000 1" x 3" White
<b>SHARED SPECIMEN</b> TEST _____ TEST _____ TEST _____ UPCR-6010 1" x 2 1/4" Fl. Red	<b>SPECIMEN UNACCEPTABLE</b> Identification on specimen label incomplete. Please resubmit. UPCR-6012 1" x 2 1/4" Fl. Pink	<b>WRAP AROUND TUBING</b>   <b>PEAK LEVEL ORDERED</b> Lab must be called when infusion is completed UPCR-7001 1" x 3" Fl. Pink
<b>DO NOT REFRIGERATE</b> UPCR-5000 3/4" x 2" Fl. Green	<b>REFRIGERATE</b> UPCR-5001 3/4" x 2" Fl. Chartreuse	

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM





## Quality Control Labels

Date Opened \_\_\_\_\_  
Exp. Date \_\_\_\_\_  
Exp. Date Control

UPCR-3012 Blue

Date Prep. \_\_\_\_\_  
Exp. Date \_\_\_\_\_  
Made By \_\_\_\_\_  
Lot No. \_\_\_\_\_

UPCR-3013 Fl. Green

**DATE RECEIVED**

UPCR-3014 Fl. Pink

Date Rec'd. \_\_\_\_\_  
Date Opened \_\_\_\_\_  
Exp. Date \_\_\_\_\_

UPCR-3015 Yellow

Date Rec'd \_\_\_\_\_  
Date Opened \_\_\_\_\_  
Exp. Date \_\_\_\_\_  
Reag. Exp. Control \_\_\_\_\_  
LOT NO.

UPCR-3016 Green

DATE RECEIVED \_\_\_\_\_  
**DO NOT USE THIS  
LOT NUMBER**

UPCR-3017 Fl. Red

**LOT #** \_\_\_\_\_  
**PARALLEL CHECKED**  
**DATE** \_\_\_\_\_

UPCR-3018 Tan

Rec'd \_\_\_\_\_ Init. \_\_\_\_\_  
Opened \_\_\_\_\_ Init. \_\_\_\_\_

UPCR-3019 Fl. Pink

THERM # \_\_\_\_\_  
TEMPC \_\_\_\_\_  
REF THERM # \_\_\_\_\_  
TEMPC \_\_\_\_\_  
DATE: \_\_\_\_\_

UPCR-3020 White

**EXPIRATION  
DATE** \_\_\_\_\_  
UPCR-1323 5/16" x 1/4"  
Yellow

**THIS LOT IS  
READY  
FOR USE**  
DATE \_\_\_\_\_ INIT \_\_\_\_\_

UPCR-3021 Green

**USE FIRST**

UPCR-3022 Fl. Chartreuse

**USE THIS  
LOT FIRST**

UPCR-3023 Fl. Chartreuse

**REVIEWED**  
DATE \_\_\_\_\_  
SUPERVISOR \_\_\_\_\_

UPCR-3025 Fl. Orange

**DATE  
OPENED**

UPCR-9900  
3/4" diameter  
Fl. Chartreuse

UPCR-3000 series labels are 3/4" x 1 1/2"

**EMPTY  
GAS  
CYLINDER**

UPCR-9103 1 1/2" x 3" Fl. Chartreuse

**LAST BOX IN THIS LOT  
Notify Supervisor  
When Opened**

UPCR-9252 1 3/4" x 2 3/4" Fl. Red

**RELEASED**  
DATE \_\_\_\_\_ TIME \_\_\_\_\_

UPCR-9104 1 1/2" x 3" Fl. Red

DATE PREP \_\_\_\_\_  
EXP. DATE \_\_\_\_\_  
MADE BY \_\_\_\_\_  
LOT NO.

UPCR-6055 1" x 2 1/4" Tan

Rec'd \_\_\_\_\_ Lot \_\_\_\_\_  
Exp. Date \_\_\_\_\_  
Date Opened \_\_\_\_\_  
Concurrent testing of this kit or  
reagent is complete & acceptable  
Concurrent test  
on this kit  
or reagent  
has **NOT** been  
**DONE**

UPCR-7022 Fl. Green 1" x 3"

Date Rec'd \_\_\_\_\_ Initial \_\_\_\_\_  
Date Checked \_\_\_\_\_ Initial \_\_\_\_\_  
Date in Use \_\_\_\_\_  
Expires \_\_\_\_\_

UPCR-6056 1" x 2 1/4" Fl. Orange

Date Rec'd \_\_\_\_\_  
Date Opened \_\_\_\_\_  
Exp. Date \_\_\_\_\_  
Refrig. Req. \_\_\_\_\_

UPCR-6057 1" x 2 1/4" Fl. Green

Date Rec'd \_\_\_\_\_  
Date Opened \_\_\_\_\_  
Exp. Date \_\_\_\_\_  
Room Temp. \_\_\_\_\_  
LOT NO.

UPCR-6058 1" x 2 1/4" White

**Expiration  
Date** \_\_\_\_\_

UPCR-6059 1" x 2 1/4" Fl. Chartreuse

**DISCARD**  
DATE \_\_\_\_\_  
AFTER TIME \_\_\_\_\_

UPCR-6067 1" x 2 1/4" Fl. Red

**LABORATORY TEST RESULTS**  
REVIEWED BY \_\_\_\_\_  
DATE \_\_\_\_\_ SHIFT \_\_\_\_\_

UPCR-6060 1" x 2 1/4" Blue

**QUALITY CONTROL  
KITS AND REAGENTS  
PARALLEL CHECKED**  
LOT # \_\_\_\_\_  
Date \_\_\_\_\_ Tech \_\_\_\_\_

UPCR-6061 1" x 2 1/4" Yellow

# Specialty Laboratory Communication Labels

## Warning and Infection Control Labels

**Blood/Body Fluid Precaution**

UPCR-7018 1" x 3" Fl. Pink

**STOP**



WASH HANDS  
BEFORE AND AFTER PATIENT CARE

UPCR-9650 2½" x 3¼"  
White



UPCR-9902 6" x 6" White/Red  
Package of 50

**STOP** the Spread of Germs . . .  
**WASH YOUR HANDS**

UPCR-7019 1" x 3" Fl. Red

**LATEX ALLERGY**

UPCR-3010 ¾" x 1½" Pink

**LATEX PRECAUTIONS**

UPCR-3011 ¾" x 1½" Blue

**AIDS PRECAUTIONS**

UPCR-1321 Fl. Red

**CAUTION**

UPCR-1240 Fl. Red

**HEPATITIS**

UPCR-1244 Fl. Red

**INFECTIOUS MATERIAL**

UPCR-1247 Fl. Chartreuse

**SPECIMEN PRECAUTION**

UPCR-1250 White/Red

**BIOHAZARD**

UPCR-1237 Fl. Red

**CAUTION: Possible Infectious Agent**

UPCR-1241 Fl. Red

**HEPATITIS PRECAUTION**

UPCR-1245 Fl. Red

**ISOLATION**

UPCR-1248 Fl. Red

**UNIVERSAL PRECAUTIONS**

UPCR-1251 White/Red

**BLEEDING PRECAUTIONS**

UPCR-1238 Fl. Red

**CAUTION TOXIC**

UPCR-1242 Green

**HIV Positive**

UPCR-1246 White/Red

**POSSIBLE HEPATITIS**

UPCR-1249 Fl. Red

**ISOLATION PATIENT**

UPCR-1322 Fl. Red


**BLOOD PRECAUTIONS**

UPCR-1239 White/Red


**CONTAMINATED**

UPCR-1243 Green

UPCR-1000 series labels are 5/16" x 1¼"



UPCR-7002 1" x 3"



UPCR-7003 1" x 3"

**AIDS PRECAUTIONS**

UPCR-6015 1" x 2¼" Fl. Red

**BLOOD PRECAUTIONS**

UPCR-6013 1" x 2¼" White/Red

**PRECAUTIONS**

TYPE \_\_\_\_\_

UPCR-6014 1" x 2¼" Fl. Chartreuse

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Multi-Purpose Blank Labels



Shamrock offers pressure-sensitive blank labels in a wide variety of sizes, materials and adhesives to answer all labeling needs. Stock labels are white and, unless noted, supplied in rolls of 1,000 labels.



## ROLL FORM BLANK LABELS

Designed for economy and efficiency, these basic, white labels feature perforations between each label for flexibility of use, permanent adhesive that holds them securely in place and a matte finish to accept information from printer, pen, pencil or marker. Stock labels are white. Labels may be tinted in any of the color choices shown below. Call for information.

Cat. No.	Size	Cat. No.	Size
SSP-1	1" x 2" or less*	SSP-7	1½" x 2"
SSP-2	1" x 2¼"	SSP-8	1½" x 2¼"
SSP-3	1" x 2½"	SSP-9	1½" x 2½"
SSP-4	1¼" x 2"	SSP-10	1¾" x 2¾"
SSP-5	1¼" x 2¼"	SSP-11	2" x 3"
SSP-6	1¼" x 2½"		

\*Please specify size desired.



## ADDRESSOGRAPH/IMPRINTER LABELS

Smudge proof, labels provide a sharp, legible imprint. They are available in your choice of standard removable adhesive for use on any clean, flat surface or with Super-Stick removable adhesive that is adaptable to plastic containers and irregular surfaces. Stock labels are white. Color choices are shown below.

Sizes	Litho Stock	Super-Stick Stock
1" x 3"	ADS-1	ADSF-1
1" x 3½"	ADS-2	ADSF-2
1¼" x 3"	ADS-3	ADSF-3
1¼" x 3½"	ADS-4	ADSF-4
1½" x 3"	ADS-5	ADSF-5
1½" x 3½"	ADS-6	ADSF-6
2" x 3"	ADS-7	ADSF-7
2½" x 4"	ADS-8	ADSF-8

## 19 Colors for Shamrock Multi-Purpose, Addressograph and VAC Labels

Call Shamrock for complete information on ordering tinted labels.

1 White	2 Yellow	3 Green	4 Red	5 Orange
6 Blue	7 Pink	8 Tan	9 Salmon	10 Copper
11 Rose	12 Violet	13 Lavender	14 Gray	15 Lime
16 Aqua	17 Chartreuse	18 Gold	19 Silver	

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Multi-Purpose Blank Labels



## VAC LABELS

Their glove proof, permanent adhesive has a special affinity for all silicone treated surfaces, stays securely in place on Vacutainer\* tubes. Write on them, use a rubber stamp or imprinter. Information stays clear, cannot be wiped off. From lab to patient and back to the lab, these labels offer a complete system of positive identification. Stock labels are white. Color choices are illustrated on page L-33.

Cat. No.	Size	Cat. No.	Size
SVL-11	1" x 2" or less**	SVL-20	1¼" x 3½"
SVL-12	1" x 2¼"	SVL-21	1½" x 2"
SVL-13	1" x 2½"	SVL-22	1½" x 2¼"
SVL-14	1" x 3"	SVL-23	1½" x 2½"
SVL-15	1" x 3½"	SVL-24	1½" x 3"
SVL-16	1¼" x 2"	SVL-25	1½" x 3½"
SVL-17	1¼" x 2¼"	SVL-26	1¾" x 2¾"
SVL-18	1¼" x 2½"	SVL-27	2" x 3"
SVL-19	1¼" x 3"		


\*Vacutainer is a registered trademark of Becton, Dickinson & Co.  
\*\*Please specify size desired.



## SATIN CLOTH LABELS

Ideal for cryogenic procedures, these cloth textured labels feature permanent adhesive, stick securely to most surfaces. They may be marked with pen, pencil or felt marker and are designed to withstand dry ice, liquid gases, deep freezes, acids, alkalis, oils, grease, salt water and soap. Stock labels are available in a variety of sizes. Labels 5/8" wide are supplied on sheets in packages of 1,000 labels, all other sizes in rolls of 1,000 labels.

Cat. No.	Size	Cat. No.	Size
SC-1	5/8" x 1" (5-up sheet form)	SC-8	1¼" x 2½"
SC-2	5/8" x 2" (5-up sheet form)	SC-9	1½" x 2"
SC-3	1" x 2"	SC-10	1½" x 2¼"
SC-4	1" x 2¼"	SC-11	1½" x 2½"
SC-5	1" x 2½"	SC-12	1¾" x 2¾"
SC-6	1¼" x 2"	SC-13	2" x 3"
SC-7	1¼" x 2¼"		



**CryoPen-26**  
Supplied in packs of 4, ink won't fade or smudge.



## TYVEK LABELS

These polyethylene coated labels offer extra strength, won't discolor, crack or tear. Permanent adhesive holds them securely on almost any surface. Write on them with pen, pencil or marker. Available in five sizes to answer a wide range of tough labeling requirements.

Cat. No.	Size	Cat. No.	Size	Cat. No.	Size
STV-1	7/8" x 7/8"	STV-3	1" x 2"	STV-5	2" x 4"
STV-2	¾" x 1½"	STV-4	1½" x 3"		

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Tape and Label Dispensers




## Tape Dispensers

### **MULTI-TAPE DISPENSERS**

Easy dispensing from one convenient location. Unique design allows fast addition or replacement of tape rolls.


- MD-250 (for 500" rolls)
- MD-350 (for 500" rolls)
- MD-500 (for 500" rolls)
- MD-2160 (for 2160" rolls)



### **HAND DISPENSERS**

Hold 500" tape rolls. Safety cutting edge protects fingers.

- SHD-1 (for tape widths 1/2", 3/4", 1")
- SHD-2 (for tape widths 1 1/2", 2")



### **STATIONARY DISPENSERS**

For use with 60 yard rolls of tape. Available in two sizes.

- SD-2 (for one or more rolls, up to 2" wide)
- SD-6 (for multiple rolls up to 6" wide)

## Label Dispensers

### **MULTIPLE ROLL LABEL DISPENSERS**

Sturdy metal for counter top use. Dispenses labels from boxes or rolls.


- LD-1 (for widths up to 1 1/2")
- LD-2 (for widths up to 2 1/2")
- LD-3 (for widths up to 4 1/4")
- LD-4 (for widths up to 8 3/4")





# Labels of Your Own Design – A Shamrock Specialty

## Custom Printed Tape



Tape Length **500"**      Perf at the end of each imprint **YES**  
 Tape Color **GREEN**      Quantity required **24 rolls**  
 Imprint Color **BLACK**

**All the unique features and flexibility of blank tape** are available in custom printed Specialty Labeling Tape.

You specify the imprint desired and choose from:

**Tape Widths** of 1/2", 3/4", 1", 1 1/2" and 2"

**Tape Lengths** of 500" and 2160"

**Imprint** in any of Shamrock's standard colors 1 thru 26.

**Tape Colors** including 26 standard and 5 fluorescent colors.

If imprint and tape colors are not specified, imprint will be supplied in black ink on white tape.

**Send us a Sample or Sketch** of the imprinted labeling tape you need.

A minimum order of 24 rolls of each imprint is required for 500" custom printed tape, 6 rolls for 2160" printed tape.

### Choose from a variety of materials available for specially designed, pre-cut labels.

Shamrock will work with you to create a label custom crafted to meet your unique requirements. Whether you need a pre-printed label or a label for on site computer printing, Shamrock can supply the special material, adhesive, size and imprint in a label formatted to fit your procedures.

#### Kromekote

Design a label that reflects the vitality of your institution. This high gloss, white stock will enhance the colors of any design or logo. Available with permanent or removable adhesive.

#### Super-Stick

Perfect for labeling small containers, irregular surfaces. Flexible stock holds tightly where ordinary labels might not work. Features smudge proof matte finish, available in permanent or removable adhesive.

#### VAC

Features a permanent, glove proof adhesive that sticks tightly to siliconized test tubes and other hard to label surfaces without flagging.

#### Tyvek®

Ideal for tough labeling jobs, this polyethylene coated material won't crack or tear under the most trying conditions. Permanent adhesive sticks tightly to almost any surface.


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM





## Custom Pre-Cut Labels



Label Color **BLUE** Adhesive **PERMANENT**  
 Imprint Color **BLACK** Perf between each label **YES**  
 Material **SUPER-STICK** Quantity required **5000**

**Call on Shamrock's Customer Service Department or your local Shamrock Sales Representative** for help in choosing the perfect label size, material and adhesive to exactly meet your needs or **send a sample or sketch of the label you have in mind**. If colors are not specified, imprint will be supplied in black ink on white label.

**Shamrock will be happy to send a price quotation by return mail, fax or email.**

## Tape and Label Colors

Choose from a variety of colors to match your color coding system, alert personnel to important information, identify department or procedure.

1 White	16 Aqua
2 Yellow	17 Chartreuse
3 Green	18 Gold
4 Red	19 Silver
5 Orange	20 Magenta
6 Blue	21 Dark Red
7 Pink	22 Maroon
8 Tan	23 Dark Brown
9 Salmon	24 Dark Blue
10 Copper	25 Dark Green
11 Rose	26 Black
12 Violet	27 Fl. Red
13 Lavender	28 Fl. Orange
14 Gray	29 Fl. Green
15 Lime	30 Fl. Chartreuse
	31 Fl. Pink

### Cold-Stick

Perfect for labeling items to be stored in temperatures below freezing. Available with permanent, glove free adhesive in either matte or glossy finish.

### Satin Cloth

Ideal for extremely cold temperatures, this cloth textured material may be marked with pen, pencil or marker, features permanent adhesive.

### Dura-Stick

A premium material for cryogenic procedures. Stands up to moisture and chemicals. The label accepts a clear imprint that won't smear or flake. Its permanent adhesive won't stick to gloves. See page L-6 for other label materials for cryogenic temperatures.

### Acetate

Destructible material, supplied clear or tinted, that provides instant proof of tampering. Permanent adhesive.

### Mylar

Clear, polyethylene material will not tear. Order imprinted or use without imprint as a protective, overlay label.

### Dupli-Copy

Providing a permanent record, material copies everything written on the label onto the label's backing paper. Pre-printing will show only on the label unless "crash printing" is specified. Label has permanent adhesive, matte finish.

### Twin Label

Get two labels at once. All information marked on the top label is copied on the underlying, second label. Pre-printing appears only on the top label unless "crash printing" is specified. Material is matte finish with permanent adhesive.

### Thermal Materials

See page L-60 or call for information on special label materials and adhesives for both direct and transfer thermal printers.



**SHAMROCK SUPPLIES BLANK LABELS COMPATIBLE WITH ALL SOFTWARE PROGRAMS INCLUDING:**

Cerner Millennium® PathNet®

Cerner Citation

Cerner DHT

CHC

CPSI

Health Science Systems

HMS

Elekta/Impac

Epic Beaker

McKesson Horizon

Meditech

MedLab

Orchard

SCC Softlab

SMS

Sunquest/Misys

VA Hospital



**SHAMROCK BLANK LABELS FOR ON-DEMAND PRINTING**

Stock blank labels are offered in a wide variety of sizes and formats, feature top quality materials and adhesives, are glove friendly and latex free. Labels are available for all laboratory software, whether used in stand-alone operations or interfaced with your facility's IT system.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Blank Labels

## FOR ON-DEMAND PRINTING

### **SHAMROCK STOCK BLANK LABELS FOR THERMAL, LASER AND DOT MATRIX PRINTERS**

#### **Thermal Direct Labels for Printers with a Heated Printhead**

#### **Thermal Transfer Labels for Printers using a Ribbon for Printing**

- Shamrock Thermal Labels are available for standard or IR scanning and may be supplied with indexing bars, holes or notches.
- Should your printer require a core size and/or number of labels per roll different from our stock, just let us know. Shamrock can supply labels to perfectly meet your requirements.

#### **Laser Labels**

- Shamrock Laser Labels, offered in a range of sizes to answer the labeling needs of every procedure in your laboratory, run trouble free through your printer.
- Stock label formats are supplied on 8 1/2" x 11" sheets, custom configurations are available upon request.

#### **Pinfed Labels for Dot Matrix Printers**

- Shamrock Pinfed Labels won't jam in your dot matrix printer, offer a clear, easy-to-read imprint to facilitate communication.
- Special label materials are available.

### **SHAMROCK SUPPLIES BLANK LABELS FOR UNIQUE LABELING REQUIREMENTS**

- Special materials that stand up to harsh temperature extremes, water baths and solvents.
- Labels tinted in whole or with color bars to highlight information. You choose from 31 stock colors shown on page L-37.

# Specialty Thermal Labels for Laboratory Software

Shamrock stock Thermal Labels are supplied on 3" cores unless otherwise noted. Catalog numbers shown indicate Thermal Direct Material unless otherwise noted. To order stock labels in Thermal Transfer Material, please add the suffix "-T" to the appropriate catalog number, for instance: THERM-396R-T. Shamrock can also supply Thermal Labels produced to your individual specifications. Call for complete information and design assistance.

**Cerner Millennium® PathNet®**

**THERM-396R**  $\frac{7}{8}'' \times 2''$

**THERM-1213R**  $1\frac{5}{16}'' \times 1\frac{5}{16}''$  (4-up)  
**THERM-1213OPR**  $1\frac{5}{16}'' \times 1\frac{5}{16}''$  (4-up) Opaque Polyester  
 Material for use with Premium Wax/Resin Transfer Printer Ribbon

**THERM-1221R**  $1\frac{1}{8}'' \times 2''$

**THERM-1222R**  $1\frac{1}{4}'' \times 2''$

**THERM-760R**  $1\frac{1}{4}'' \times 2''$

**THERM-550R**  $1\frac{3}{8}'' \times 8\frac{1}{2}''$

**THERM-552R**  $1\frac{3}{4}'' \times 8\frac{1}{2}''$  ( $1\frac{27}{32}''$  repeat)

**THERM-553R**  $2\frac{1}{8}'' \times 8\frac{1}{2}''$

**THERM-1295R**  $2\frac{3}{16}'' \times 2''$  White  
**THERM-1295R-2**  $2\frac{3}{16}'' \times 2''$  Yellow (Die Cuts)

**THERM-1294R**  $2\frac{3}{16}'' \times 2''$  White  
**THERM-1294R-2**  $2\frac{3}{16}'' \times 2''$  Yellow (Perforations)





Stock and Special Thermal Label Materials detailed on page L-60.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## Cerner Millennium® PathNet®




Most stock labels are white. Color may be added to your specifications.


THERM-1298R 2<sup>3</sup>/<sub>16</sub>" x 7<sup>1</sup>/<sub>2</sub>"  
THERM-1299R 2<sup>3</sup>/<sub>16</sub>" x 7<sup>1</sup>/<sub>2</sub>" Horizontal Slit in Liner  
(Die Cuts)




THERM-554R 2<sup>1</sup>/<sub>4</sub>" x 8<sup>1</sup>/<sub>2</sub>"




THERM-371R 2<sup>1</sup>/<sub>2</sub>" x 2<sup>1</sup>/<sub>8</sub>" (2<sup>3</sup>/<sub>4</sub>" repeat)



THERM-7825R 2<sup>1</sup>/<sub>2</sub>" x 2<sup>1</sup>/<sub>4</sub>"



THERM-1300R 2<sup>1</sup>/<sub>2</sub>" x 2<sup>1</sup>/<sub>4</sub>"



THERM-1301R 2<sup>1</sup>/<sub>2</sub>" x 2<sup>1</sup>/<sub>4</sub>"  
THERM-1301R-IR 2<sup>1</sup>/<sub>2</sub>" x 2<sup>1</sup>/<sub>4</sub>" (IR Readable)

# Specialty Thermal Labels for Laboratory Software

**Cerner Millennium® PathNet®**

**THERM-782R** 2½" x 2¼" (Die Cuts)

**THERM-489R** 2½" x 2¼" (Perforations)

**THERM-785R** 2½" x 4"

**THERM-783R** 2½" x 2½"

**THERM-343R** 2½" x 4"  
**THERM-343R-1C** 2½" x 4" (1" core)

**THERM-786R** 2½" x 5¼" (Die Cuts)

**THERM-787R** 2½" x 5¼" (Perforations)

**THERM-327R** 2½" x 8½" (Die Cuts)

**THERM-344R** 2½" x 8½" (Perforations)

Stock and Special Thermal Label Materials detailed on page L-60.


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM







## Cerner Millennium® PathNet®




THERM-556R 2 1/2" x 8 1/2"




THERM-1303R 3" x 5" White  
THERM-1303R-2 3" x 5" Yellow




THERM-623R 3 9/16" x 4"




THERM-571R 4 1/2" x 4"  
THERM-1306R 4 1/2" x 4" Vertical perforation thru  
liner 1/16" right of center




THERM-1307R 4 1/2" x 4"  
THERM-1307R-Z 4 1/2" x 4" (1 1/2" Core for use with  
Zebra printer)




THERM-1305R 4 1/2" x 4"



4 1/2" x 4"  
THERM-1305R-2 Yellow  
THERM-1305R-4 Red  
THERM-1305R-6 Blue  
(With Highlight Bars)



4 1/2" x 4"  
THERM-788R-2 Yellow  
THERM-788R-4 Red  
THERM-788R-6 Blue  
(With Highlight Bars)




THERM-788R 4 1/2" x 4"


Most stock labels are white. Color may be added to your specifications.

# Specialty Thermal Labels for Laboratory Software


## Cerner Citation




THERM-513R 1" x 2"  
 THERM-513R-1C 1" x 2" (1" core)  
 THERM-513R-OP 1" x 2" (Opaque Polyester)



THERM-775R 1" x 3 1/8"




THERM-1457R 1" x 3 1/8"




THERM-1225R 1 1/4" x 3 1/8"

## Cerner DHT




THERM-1453R 1 3/4" x 8 15/32"


Stock and Special Thermal Label Materials detailed on page L-60.



THERM-1201R 4 7/8" x 3 9/32"



THERM-1455R 5 15/16" x 3 3/10"




THERM-1456R 7" x 2 3/8"

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



## McKesson Horizon




THERM-1212R  $7/8" \times 3^{3/8}"$




THERM-304R  $1^{3/8} \times 2^{15/16}"$


Most stock labels are white. Color may be added to your specifications.




THERM-325R  $1^{3/16} \times 3^{15/16}"$



THERM-1451R  $1^{7/16} \times 4"$




THERM-1451R-B5  $1^{7/16} \times 4"$  Orange Border



THERM-1452R  $6" \times 4"$


# Specialty Thermal Labels for Laboratory Software

Sunquest/Misys




1" x 2 1/4"  
THERM-1319R  
THERM-1319R-KT

(1/8" gap between labels)  
Thermal Transfer DuraStick  
Material for Cold Temperatures




7/8" x 7/8" (4-up)  
THERM-1318R  
THERM-1318R-T  
THERM-1318OPR

Thermal Direct  
Thermal Transfer  
Thermal Transfer Opaque  
Polyester, Xylene Resistant




1 1/4" x 2"  
THERM-760R




1 3/16" x 4 3/32"  
(Non-IR Readable)  
THERM-451R White  
THERM-451R-2 Yellow

(IR Readable)  
THERM-3R White  
THERM-3R-2 Yellow




1 3/16" x 4 3/32"  
(Non-IR Readable)  
THERM-451R-B2 Yellow Border  
THERM-451R-B3 Green Border  
THERM-451R-B4 Red Border  
THERM-451R-B6 Blue Border

(IR Readable)  
THERM-3R-B2 Yellow Border  
THERM-3R-B3 Green Border  
THERM-3R-B4 Red Border  
THERM-3R-B6 Blue Border




1 3/16" x 4 3/32" Core Size 1 1/2"  
(Non-IR Readable)  
THERM-1308R White  
THERM-1308R-2 Yellow  
THERM-1308R-SB (Sensor Bar)

(IR Readable)  
THERM-15R White  
THERM-15R-2 Yellow  
THERM-15R-SB (Sensor Bar)




1 3/16" x 4 3/32" Core Size 1 1/2"  
(Non-IR Readable)  
THERM-1308R-B2 Yellow Border  
THERM-1308R-B3 Green Border  
THERM-1308R-B4 Red Border  
THERM-1308R-B6 Blue Border

(IR Readable)  
THERM-15R-B2 Yellow Border  
THERM-15R-B3 Green Border  
THERM-15R-B4 Red Border  
THERM-15R-B6 Blue Border



1 3/16" x 4 3/32"  
THERM-1357R



1 3/16" x 4 3/32"  
THERM-1357R-L


Stock and Special Thermal Label Materials detailed on page L-60.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## Sunquest/Misys




1 7/16" x 3 5/8" (Die Cuts)  
THERM-374R White THERM-374R-2 Yellow




1 7/16" x 3 5/8" (Die Cuts)  
THERM-374R-B2 Yellow Border THERM-374R-B4 Red Border  
THERM-374R-B3 Green Border THERM-374R-B6 Blue Border




1 7/16" x 3 5/8" (Perforations)  
THERM-1309R White THERM-1309R-2 Yellow




1 7/16" x 3 5/8" (Perforations)  
THERM-1309R-B2 Yellow Border THERM-1309R-B4 Red Border  
THERM-1309R-B3 Green Border THERM-1309R-B6 Blue Border




1 7/16" x 3 1/16" (Die Cuts)  
THERM-541R White THERM-541R-2 Yellow



1 7/16" x 3 1/16" (Die Cuts)  
THERM-541R-B2 Yellow Border THERM-541R-B4 Red Border  
THERM-541R-B3 Green Border THERM-541R-B6 Blue Border




1 7/16" x 3 1/16" (Perforations)  
THERM-1310R White THERM-1310R-2 Yellow



1 7/16" x 3 1/16" (Perforations)  
THERM-1310R-B2 Yellow Border THERM-1310R-B4 Red Border  
THERM-1310R-B3 Green Border THERM-1310R-B6 Blue Border


# Specialty Thermal Labels for Laboratory Software

## Sunquest/Misys




THERM-1314R 5 1/4" x 2 5/8"  
(1/4" gap between labels)  
Core Size 1 1/2" or 3" – please specify

Stock and Special Thermal Label Materials detailed on page L-60.




THERM-1315R 6 5/8" x 3 3/4"  
(1/4" gap between labels)  
Core Size 1 1/2"




THERM-1316R 7 1/8" x 4"


## Meditech




THERM-478R 1 1/4" x 4"  
THERM-4R 1 1/4" x 4" (IR Readable)




THERM-478R-B4 1 1/4" x 4" Red Border  
THERM-4R-B4 1 1/4" x 4" (IR Readable) Red Border  
Labels also available with Yellow (B2), Green (B3) or Blue (B6) Border



THERM-358R-P 2" x 4"



THERM-1927R 2 3/4" x 4"




THERM-1794R 5 1/4" x 3 3/8"







## SCC Softlab




THERM-466R 1" x 2½"  
 THERM-466R-1C 1" x 2½" (1" core)  
 THERM-2R 1" x 2½" (IR Readable)



THERM-1458R 1 3/8" x 3"




THERM-1458R-B4 1 3/8" x 3" Red Border




THERM-1458R-B2 1 3/8" x 3" Yellow Border


## SMS




THERM-792R 6" x 4¼" (6½" repeat)  
 THERM-7R 6" x 4¼" (6½" repeat) (IR Readable)



THERM-1459R 3" x 8 15/32"



THERM-792R-4 6" x 4¼" (6½" repeat)  
 THERM-7R-4 6" x 4¼" (6½" repeat)  
 (IR Readable)



THERM-792R-3 6" x 4¼" (6½" repeat)  
 THERM-7R-3 6" x 4¼" (6½" repeat)  
 (IR Readable)

# Specialty Thermal Labels for Laboratory Software

## VA Hospital

Most stock labels are white. Color may be added to your specifications.

THERM-329R 4" x 2<sup>3</sup>/<sub>8</sub>"  
THERM-M-10R 4" x 2<sup>3</sup>/<sub>8</sub>" (IR Readable)

THERM-1788R 4" x 3<sup>3</sup>/<sub>8</sub>"

THERM-1807R 5<sup>1</sup>/<sub>2</sub>" x 3<sup>1</sup>/<sub>4</sub>"

THERM-1868R 5<sup>1</sup>/<sub>2</sub>" x 3<sup>1</sup>/<sub>4</sub>"

## Health Science Systems

Stock and Special Thermal Label Materials detailed on page L-60.

THERM-1454R 3" x 4<sup>1</sup>/<sub>8</sub>"

## HMS

THERM-1450R 5<sup>1</sup>/<sub>2</sub>" x 4"

## CPSI


THERM-757R 1<sup>5</sup>/<sub>16</sub>" x 3<sup>1</sup>/<sub>2</sub>"

THERM-1214R-L 1<sup>5</sup>/<sub>16</sub>" x 3<sup>1</sup>/<sub>2</sub>"

THERM-1214R 1<sup>5</sup>/<sub>16</sub>" x 3<sup>1</sup>/<sub>2</sub>"



## Additional Thermal Labels for Laboratory Software



THERM-100R  $\frac{3}{8}$ " x  $1\frac{1}{4}$ "



THERM-366R  
 $\frac{1}{2}$ " x  $\frac{3}{8}$ "




THERM-487R  $\frac{1}{2}$ " x  $1\frac{1}{2}$ "



THERM-1987R  $\frac{1}{2}$ " x 2"  
THERM-1987R-OP  $\frac{1}{2}$ " x 2" (Opaque Polyester)




THERM-504R  
 $\frac{5}{8}$ " x  $\frac{7}{8}$ "



THERM-617R  $\frac{5}{8}$ " x 1"


Stock and Special  
Thermal Label  
Materials detailed  
on page L-60.




THERM-754R  $\frac{5}{8}$ " x  $3\frac{1}{4}$ "



THERM-643R  $\frac{5}{8}$ " x  $2\frac{1}{2}$ "



THERM-425R  $\frac{5}{8}$ " x  $2\frac{5}{8}$ "




THERM-547R  $\frac{5}{8}$ " x  $1\frac{1}{16}$ "




THERM-528R  $\frac{5}{8}$ " x 4"




THERM-756R  $1\frac{1}{16}$ " x  $3\frac{1}{2}$ "




THERM-305R  $\frac{3}{4}$ " x  $\frac{3}{4}$ "




THERM-566R  $\frac{3}{4}$ " x  $1\frac{1}{2}$ "




THERM-520R  $\frac{3}{4}$ " x 2"  
THERM-520R-1C  $\frac{3}{4}$ " x 2" (1" core)




THERM-364R  $\frac{3}{4}$ " x  $2\frac{1}{4}$ "



THERM-486R  $\frac{3}{4}$ " x  $2\frac{1}{2}$ "



THERM-755R  $\frac{7}{8}$ " x  $1\frac{1}{4}$ "



THERM-628R  $\frac{7}{8}$ " x  $1\frac{1}{2}$ "

# Specialty Thermal Labels for Laboratory Software

## Additional Thermal Labels for Laboratory Software

**THERM-491R**  $\frac{7}{8}'' \times 1\frac{27}{32}''$   
**THERM-108R**  $\frac{7}{8}'' \times 4''$   
**THERM-123R**  $1'' \times 1\frac{15}{16}''$   
**THERM-570R**  $1'' \times 1\frac{3}{4}''$   
**THERM-125R**  $1'' \times 1\frac{1}{2}''$  (Non-IR Readable)  
**THERM-5R**  $1'' \times 1\frac{1}{2}''$  (IR Readable)  
**THERM-263R**  $1'' \times 2\frac{1}{16}''$  (1 $\frac{1}{4}''$  repeat)  
**THERM-632R**  $1'' \times 2\frac{5}{16}''$  (1 $\frac{3}{16}''$  repeat)  
**THERM-466R**  $1'' \times 2\frac{1}{2}''$   
**THERM-466R-1C**  $1'' \times 2\frac{1}{2}''$  (1" core)  
**THERM-492R**  $1'' \times 3''$   
**THERM-492R-1C**  $1'' \times 3''$  (1" core)  
**THERM-CEQ**  $1'' \times 3''$  (Transfer Kimdura Material)  
**THERM-CEQ-OP**  $1'' \times 3''$  (White Opaque Polyester Material)  
**THERM-330R**  $1'' \times 3\frac{1}{4}''$   
**THERM-467R**  $1'' \times 3''$  (1 $\frac{1}{4}''$  repeat)  
**THERM-315R**  $1'' \times 3\frac{1}{2}''$   
**THERM-126R**  $1'' \times 3\frac{1}{2}''$   
**THERM-126R-1C**  $1'' \times 3\frac{1}{2}''$  (1" core)  
**THERM-1894R**  $1'' \times 4''$   
**THERM-1894R-1C**  $1'' \times 4''$  (1" core)  
**THERM-524R**  $1\frac{1}{16}'' \times 1\frac{11}{16}''$   
**THERM-456R**  $1\frac{1}{8}'' \times 1\frac{7}{8}''$   
**THERM-759R**  $1\frac{1}{8}'' \times 2\frac{1}{4}''$   
**THERM-2017R**  $1'' \times 1''$   
**THERM-2017R-1C**  $1'' \times 1''$  (1" core)  
**THERM-2017R-OP**  $1'' \times 1''$  (Opaque Polyester)  
**THERM-2017ROP1C**  $1'' \times 1''$  (Opaque Polyester) (1" core)

Most stock labels are white. Color may be added to your specifications.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## Additional Thermal Labels for Laboratory Software

Stock and Special Thermal Label Materials detailed on page L-60.



THERM-128R 1 1/8" x 3 1/2"




THERM-545R 1 3/16" x 2"



THERM-558R  
1 1/4" x 2 1/4"  
THERM-558R-1C  
1 1/4" x 2 1/4" (1" core)


THERM-751R 1 1/4" x 2 1/2"



THERM-776R 1 1/4" x 4"

5/8" x 1"  
5/8" x 1"

1 1/4" x 2"



THERM-314R 1 3/8" x 1 3/8"


7/8" x 1 3/8"

1/2" x 1 3/8"

THERM-567R 1 3/8" x 2"



THERM-132R 1 3/8" x 3 3/4"




THERM-777R 1 3/8" x 4"


1 3/8" x 2 1/2"

1 3/8" x 5/8"


1 1/16" x 7/8"




THERM-133R 1 3/8" x 4 1/4"



THERM-523R 1 1/2" x 1"



THERM-761R 1 1/2" x 2"



THERM-1862R 1 1/2" x 4"  
THERM-1862R-1C 1 1/2" x 4" (1" core)  
THERM-1862R-OP 1 1/2" x 4" (Opaque Polyester)

3/4" x 3"


THERM-427R 1 1/2" x 1 1/2"

THERM-1R 1 1/2" x 1 1/2" (IR Readable)


THERM-779R 1 1/2" x 3"

# Specialty Thermal Labels for Laboratory Software

## Additional Thermal Labels for Laboratory Software




THERM-613R 1 $\frac{5}{8}$ " x 1 $\frac{5}{8}$ "




THERM-750R 1 $\frac{3}{4}$ " x 1 $\frac{1}{4}$ "


Most stock labels are white. Color may be added to your specifications.




THERM-150R 1 $\frac{3}{4}$ " x 2 $\frac{7}{8}$ "




THERM-563R 1 $\frac{3}{4}$ " x 3"




THERM-780R 1 $\frac{3}{4}$ " x 3 $\frac{11}{16}$ "




THERM-551R 1 $\frac{3}{4}$ " x 8 $\frac{7}{16}$ "




THERM-608R 1 $\frac{31}{32}$ " x 2 $\frac{1}{2}$ "




THERM-7554R 1 $\frac{7}{8}$ " x 1 $\frac{7}{8}$ "



THERM-1863R 2" x 1"  
THERM-1863R-OP 2" x 1"  
(Opaque Polyester)



THERM-534R 2" x 2"



THERM-506R 2" x 3"  
THERM-506R-1C 2" x 3" (1"core)


SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM







## Additional Thermal Labels for Laboratory Software




THERM-331R 2" x 3¼"




THERM-358R 2" x 4"  
THERM-358R-1C 2" x 4" (1"core)




THERM-781R 2⅛" x 4⅞"



THERM-161R 2⅝" x 3¼"



THERM-164R 2⅝" x 3¾"




THERM-535R 2⅝" x 4"  
THERM-535R-1C 2⅝" x 4" (1"core)




THERM-1312R 2½" x 1"


Stock and Special Thermal Label Materials detailed on page L-60.




THERM-165R 2⅝" x 5⅞"




THERM-479R  
2½" x 1⅞"




THERM-177R 2½" x 4" (2¼" repeat)




THERM-784R 2½" x 4"



THERM-644R 2½" x 4"



THERM-362R 2¾" x 2"



THERM-555R 2½" x 7½"

# Specialty Thermal Labels for Laboratory Software

**Additional Thermal Labels for Laboratory Software**

**Label Specifications:**

- THERM-557R** 2<sup>3</sup>/<sub>4</sub>" x 8<sup>1</sup>/<sub>2</sub>"
- THERM-611R** 2<sup>7</sup>/<sub>8</sub>" x 4<sup>1</sup>/<sub>8</sub>"
- THERM-316R** 3" x 3"
- THERM-316R-1C** 3" x 3" (1"core)
- THERM-320R** 3" x 4"
- THERM-320R-1C** 3" x 4" (1"core)
- THERM-324R** 3<sup>1</sup>/<sub>8</sub>" x 3<sup>3</sup>/<sub>4</sub>"
- THERM-345R** 3<sup>1</sup>/<sub>8</sub>" x 2<sup>1</sup>/<sub>4</sub>"
- THERM-321R** 3<sup>1</sup>/<sub>4</sub>" x 1<sup>15</sup>/<sub>16</sub>" (3<sup>1</sup>/<sub>2</sub>" repeat)
- THERM-341R** 3" x 1"
- THERM-306R** 3" x 1<sup>1</sup>/<sub>8</sub>"
- THERM-356R** 3" x 1<sup>1</sup>/<sub>2</sub>"
- THERM-310R** 3<sup>3</sup>/<sub>8</sub>" x 2<sup>1</sup>/<sub>8</sub>"
- THERM-529R** 3<sup>3</sup>/<sub>8</sub>" x 5<sup>1</sup>/<sub>4</sub>"
- THERM-642R** 2<sup>7</sup>/<sub>8</sub>" x 3<sup>1</sup>/<sub>4</sub>"
- THERM-398R** 2<sup>7</sup>/<sub>8</sub>" x 4"
- THERM-752R** 3" x 1<sup>3</sup>/<sub>4</sub>"
- THERM-500R** 3" x 5" White
- THERM-500R-2** 3" x 5" Yellow



Most stock labels are white. Color may be added to your specifications.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## Additional Thermal Labels for Laboratory Software



## Specialty Thermal Labels for Bedside Printing

Shamrock Blank Labels for handheld printers provide clear imprints and feature glove friendly adhesive. They will not fall off specimen tubes during transport or testing.




Shamrock can provide labels to suit any handheld imprinter, belt or software. Call for complete information.

# Specialty Thermal Labels for Desktop Printers


Shamrock offers thermal labels for on-demand use with your small desktop printer. Labels are precisely designed for trouble free use in these specialized printers, producing clear, easy to read text and/or bar codes. Stock desktop printer labels feature direct thermal material and are supplied on a 1" core.




THERM-1839R-DO  
1<sup>5</sup>/<sub>16</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
750 per roll



THERM-1839RSQDO  
1<sup>5</sup>/<sub>16</sub>" x 1<sup>5</sup>/<sub>16</sub>"  
750 per roll




THERM-1840R-DO  
1" x 1"  
750 per roll



THERM-2018R-DO  
1" x 1"  
750 per roll




THERM-1791R-DO  
1<sup>1</sup>/<sub>4</sub>" x 2<sup>1</sup>/<sub>4</sub>"  
1,000 per roll




THERM-1416R-DO  
2<sup>1</sup>/<sub>8</sub>" x 1"  
500 per roll



THERM-1924R-DO  
2<sup>1</sup>/<sub>8</sub>" x 1"  
500 per roll




THERM-1917R-DO  
2<sup>1</sup>/<sub>4</sub>" x 1<sup>1</sup>/<sub>16</sub>"  
500 per roll




THERM-1918R-DO  
2<sup>3</sup>/<sub>4</sub>" x 2<sup>1</sup>/<sub>8</sub>"  
320 per roll



THERM-1919R-DO  
3<sup>1</sup>/<sub>16</sub>" x 1<sup>13</sup>/<sub>16</sub>"  
300 per roll



THERM-1953R-DO  
3<sup>1</sup>/<sub>16</sub>" x 2<sup>1</sup>/<sub>16</sub>" (Die Cut)  
300 per roll




THERM-1920R-DO  
3<sup>7</sup>/<sub>16</sub>" x 1<sup>1</sup>/<sub>8</sub>"  
350 per roll




THERM-1921R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>16</sub>"  
300 per roll




THERM-1823R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>3</sup>/<sub>32</sub>"  
350 per roll




THERM-1989R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>3</sup>/<sub>32</sub>" (2-up)  
700 per roll



THERM-TL27R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>10</sub>"  
350 per roll



THERM-1099R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>8</sub>"  
350 per roll



THERM-TL55R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>10</sub>"  
350 per roll



THERM-1922R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>13</sup>/<sub>32</sub>"  
350 per roll




THERM-1805R-DO  
3<sup>1</sup>/<sub>2</sub>" x 1<sup>15</sup>/<sub>16</sub>"  
350 per roll



THERM-1710R-DO  
3<sup>1</sup>/<sub>2</sub>" x 2"  
300 per roll

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



**STANDARD AND FLUORESCENT COLORS AVAILABLE**

Stock labels for small desktop printers are available in your choice of 13 colors. To order, simply add the appropriate color number listed in the color box below to the label's catalog number. For example:

**THERM-1791R-DO** in yellow would be ordered as **THERM-1791RDO-2**.

**8 STANDARD COLORS**



**5 FLUORESCENT COLORS**



**CUSTOM OPTIONS**

Shamrock can print desktop labels with static information, using color to enhance your facility address or call attention to important information. Call for complete information.



# Material Options for Specialty Thermal Labels

## STANDARD THERMAL MATERIALS

**Thermal Direct** Permanent adhesive won't stick to gloves. May be applied at 32°F (0°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C).

**Thermal Transfer** Bright white in color and smudge resistant. Permanent adhesive won't stick to gloves. May be applied at 32°F (0°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C). Sticks to most surfaces including plastic and curved.

**Stock Thermal Labels are supplied on 3" cores for Thermal Direct printing unless otherwise specified. To order stock labels in Thermal Transfer material, please add the suffix "-T" to the catalog number.**



## PREMIUM THERMAL MATERIALS

**Thermal Direct Kimdura®** Solvent resistant, imprint won't smear or flake. Glove free permanent adhesive. May be applied at 10°F (-12°C), holds in temperatures ranging from -65°F to 300°F (-54°C to 148°C).

**Transfer Smudgeproof Kimdura®** Glove free, solvent resistant, aggressive permanent adhesive. May be applied at 10°F (-12°C), holds in temperatures ranging from -40°F to 300°F (-40°C to 149°C).

**To order stock labels in Kimdura Direct add the suffix "-KD", for Kimdura Transfer add the suffix "-KT".**

## THERMAL MATERIALS FOR SPECIAL APPLICATIONS

**Thermal Direct IR Readable** Smudge resistant, perfect for curved and plastic surfaces. Glove free permanent adhesive. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).

**Thermal Direct Piggyback** Holds tightly to a wide variety of surfaces, won't stick to gloves. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).

**Thermal Direct Tight Release** Smudge resistant, ideal for serpentine feed printers. Glove free permanent adhesive. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).

**Thermal Direct Removable** Smudge resistant material. Glove free removable adhesive leaves no residue. May be applied at 32°F (0°C), holds in temperatures ranging from -40°F to 200°F (-40°C to 93°C).

**Thermal Transfer Piggyback** Bright white, smudge resistant material won't stick to gloves, holds tightly to most surfaces. May be applied at 32°F (0°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).

**Transfer PolyFlex Piggyback** High performance, clear carrier sheet, excellent flexibility. Glove free permanent adhesive. May be applied at 40°F (5°C), holds in temperatures ranging from -40°F to 302°F (-40°C to 150°C).

**Transfer Smudgeproof Kimdura® Piggyback** Resists moisture and chemicals. Glove free permanent adhesive. May be applied at -20°F (-29°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C).

**Thermal Transfer Matte Cold-Stick** Glove free, permanent adhesive. May be applied at -10°F (-23°C), holds in temperatures ranging from -65°F to 200°F (-53°C to 93°C). Good initial adhesion that becomes permanent after 24 hours.

**Thermal Transfer Cryo-Stick** Glove free, permanent adhesive for cryogenic testing and storage applications. May be applied at 50°F (10°C), holds in temperatures ranging from -112°F to 176°F (-80°C to 80°C).

**Thermal Transfer Liqui-Stick** Glove free, permanent adhesive to withstand liquid nitrogen and dry ice. May be applied at -20°F (-29°C), holds in temperatures ranging from -320°F to 190°F (-196°C to 90°C).

**Thermal Continuous Rolls** Shamrock's top quality standard and special application thermal materials may be ordered as continuous rolls for thermal direct and transfer printers with automatic cutters.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

# Specialty Products for Thermal Transfer Printers

## RIBBONS FOR THERMAL PRINTERS

Shamrock offers ribbons in a variety of materials, widths and colors to best serve the label's application and run trouble free through your thermal transfer printer.

### For all Thermal Transfer Printers including:

Zebra • Intermec • Blaster • Blazer • Datamax • Sato • Dymo • Fargo  
Brady • Epson • Printronix • Citizen • Weber • Cognitive • And More.



### Wax-Based Ribbons

Premium wax-based film, manufactured with strict quality control. Superior back-coating provides effective ribbon flow, reduces printhead wear and ensures high performance on a wide variety of facestocks.

- Bar Coded Labels
- Patient Charge Labels

### Wax/Resin Ribbons

Wax/resin formulation ribbons are an effective alternative to resin when matched with correct materials and applications. These ribbons provide superior durability when matched with paper facestocks and provide smudge resistance when matched with softer synthetics.

- Bar Coded Labels
- Extended Use Labels

### Resin Ribbons

Resin ribbons are excellent for printing on polyester, vinyl and other synthetic materials. When matched with the proper facestocks, these ribbons provide superior smudge and smear resistance.

- Extreme Environments
- Asset Labels

**To find a ribbon with exact properties you require, simply give us the following information.**

- Make/model of printer (for instance, Zebra TLP 2844)
- Dimensions of the label to be printed
- Ribbon color desired
- Label usage requirements (temperature extremes, chemical contact, etc.)

**Shamrock can provide the Ribbon ideally suited to your thermal transfer printing needs.**

## THERMAL PRINTER CLEANING AIDS

Easy to use products offer simple solutions for maintaining thermal printer reliability and consistent quality in printing.

### CLEAN CARDS

Saturated with EZ Clean 01 Solution. These 4" x 6" cards safely and effectively remove dirt and ink from thermal print heads, paper guides and paper paths.

Supplied 20 cards per package.

Catalog # CLEAN-CARD

### SWABS

Ideal for cleaning hard to reach areas inside thermal printers. Swab stem is snapped to release Isopropyl Alcohol cleaning solution to sponge head. 4" long swabs are supplied 10 per package.

Catalog # SWABS-1

### WIPES

Cleaning wipes for thermal print heads, platen rollers. 5" x 5½" wipes are saturated with EZ Clean 01 Solution. Supplied 10 wipes per package.

Catalog # WIPES-1






# Specialty Labels for Laser Printers

Shamrock's top quality labels run trouble free through your laser printer. Labels are supplied on 8 1/2" x 11" sheets. Shamrock can customize any 8 1/2" x 11" sheet to contain the exact size label you require. Call for complete information.

## GENERAL PURPOSE LASER LABELS


Cat. No.	Label Size	Labels/Sheet	Sheets/Carton
MSL-100-FLAT	1" x 2"	44	114
MSL-175-FLAT	1 3/4" x 2 3/4"	18	139
MSL-275-FLAT	2 3/4" x 4 1/4"	8	157
LASER-1	1/2" x 2"	80	63
LASER-2	7/8" x 2 5/8"	33	152
LASER-3	1" x 2 1/8"	40	125
LASER-4	1" x 2 5/8"	30	167
LASER-5	1" x 4"	20	125
LASER-6	1 1/16" x 4"	14	179
LASER-7	2" x 4"	10	250
LASER-8	2" x 4 1/4"	10	250
LASER-9	2 1/2" x 4"	8	313
LASER-10	2 3/4" x 4 1/4"	8	313
LASER-11	3 1/3" x 4"	6	84
LASER-12	3 2/3" x 8 1/2"	3	167
LASER-13	5" x 4"	4	125
LASER-14	5 1/2" x 8 1/2"	2	250
LASER-15	11" x 8 1/2"	1	250
LASER-20	1/2" x 1"	160	32
LASER-21	1/2" x 1 1/4"	126	40
LASER-22	1/2" x 1 3/4"	80	63
LASER-23	15/16" x 1 5/16"	88	57
LASER-24	3/4" x 2 5/8"	42	119
LASER-100	1" x 2 5/8"	30	167
LASER-100-PERF	1" x 2 5/8"	30	167
LASER-101	1/2" x 2"	80	63
LASER-104	1/2" x 1 1/2"	100	50
LASER-102	1 1/2" x 2"	28	179
LASER-103	1 3/4" x 2 3/4"	18	139
LASER-105RC	2" x 3"	10	250
LASER-DOT-1	1" Dia.	63	80




MSL-100-FLAT-6  
1" x 2" Blue




MSL-175-FLAT-1  
1 3/4" x 2 3/4" White




MSL-275-FLAT-1  
2 3/4" x 4 1/4" White




LASER-1-1  
1/2" x 2" White




LASER-2-1  
7/8" x 2 5/8" White




LASER-3-30  
1" x 2 1/8" Fl. Chartreuse



LASER-4-2  
1" x 2 5/8" Yellow



LASER-5-7  
1" x 4" Pink



LASER-6-11  
1 1/16" x 4" Rose




## COLOR CHOICES FOR LASER LABELS ON PAGES L-62 THRU L-64

### STANDARD COLORS




### FLUORESCENT COLORS




To order Blank Laser Labels in color, please add the color number shown to the label's catalog number. For example: LASER-5 in orange would be ordered as LASER-5-5, in fluorescent orange as LASER-5-28.

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104


TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM



LASER-7-1  
2" x 4" White




LASER-8-1  
2" x 4 1/4" White



LASER-9-2  
2 1/2" x 4" Yellow



LASER-10-1  
2 3/4" x 4 1/4" White




LASER-11-3  
3 1/3" x 4" Green



LASER-12-1  
3 2/3" x 8 1/2" White

Printing Templates  
are available for  
all Laser Labers.




LASER-13-1  
5" x 4" White




LASER-14-1  
5 1/2" x 8 1/2" White




LASER-15-28  
11" x 8 1/2" Fl. Orange




LASER-20-1  
1/2" x 1" White




LASER-21-31  
1/2" x 1 1/4" Fl. Pink




LASER-22-1  
1/2" x 1 3/4" White




LASER-23-1  
1 5/16" x 1 5/16" White




LASER-24-31  
3/4" x 2 5/8" Fl. Pink




LASER-100-1  
1" x 2 5/8" White




LASER-100-PERF-1  
1" x 2 5/8" White




LASER-101-7  
1/2" x 2" Pink




LASER-104-1  
1/2" x 1 1/2" White




LASER-102-6  
1 1/2" x 2" Blue




LASER-103-1  
1 3/4" x 2 3/4" White




LASER-105RC-2  
2" x 3" Yellow



LASER-DOT-1-1  
1" Dia. White



LASER-DOT-11-1  
1 1/2" Dia. White




LASER-DOT-12-30  
1/2" Dia. Fl. Chartreuse


# Specialty Labels for Laser Printers

## PIGGYBACK LASER LABELS


Cat. No.	Label Size	Labels/Sheet	Sheets/Carton
MSL-4FLAT-PB	5/8" x 1.7"	44	114
MSL-4FPB-TAB	7/8" x 1 1/16"	40	125
MSL-4MT-PB	7/8" x 2"	40	125



MSL-4FLAT-PB-2  
5/8" x 1.7" Yellow  
(carrier size 1" x 2")




MSL-4FPB-TAB-1  
7/8" x 1 1/16" White  
(carrier size 1" x 2")




MSL-4MT-PB-7  
7/8" x 2" Pink  
(carrier size 1" x 2 1/8")

## MICROSCOPE SLIDE LASER LABELS


Cat. No.	Label Size	Labels/Sheet	Sheets/Carton
SLS-LASER	1 5/16" x 1 5/16"	88	57
SLS-LASER-RC	1 5/16" x 1 5/16"	80	63
SES-LASER	1 5/32" x 1 5/16"	168	30



SLS-LASER-1  
1 5/16" x 1 5/16" White  
(square corner)



SLS-LASER-RC-1  
1 5/16" x 1 5/16" White  
(round corner)




SES-LASER-6  
1 5/32" x 1 5/16" Blue  
(end labels)

# Specialty Laser Labels for Cryogenic Procedures


## DURA-STICK LASER LABELS FOR -54°C

Glove Friendly • Smudge Resistant


Cat. No.	Label Size	Labels/Sheet	Sheets/Carton
LASER-KDS-43	7/16" Dia.	212	24
LASER-KDS-59	1/2" x 9/10"	176	29
LASER-KDS-5125	1/2" x 1 1/4"	120	42
LASER-KDS-515	1/2" x 1 1/2"	100	500
LASER-KDS-759	3/4" x 9/10"	104	48
LASER-KDS-7512	3/4" x 1 1/4"	84	60
LASER-KDS-7515	3/4" x 1 1/2"	70	72




LASER-KDS-43-30  
7/16" Dia. Fl. Chartreuse




LASER-KDS-59-1  
1/2" x 9/10" White




LASER-KDS-5125-11  
1/2" x 1 1/4" Rose




LASER-KDS-515-2  
1/2" x 1 1/2" Yellow



LASER-KDS-759-1  
3/4" x 9/10" White



LASER-KDS-7512-1  
3/4" x 1 1/4" White



LASER-KDS-7515-3  
3/4" x 1 1/2" Green

SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM




## LIQUI-STICK LASER LABELS FOR -196°C

Withstand Liquid or Vapor Phase Liquid Nitrogen.  
Resistant to UV and Gama Irradiation • Smudge Resistant Surface

Cat. No.	Label Size	Description	Labels/Sheet	Sheets/Carton
* LASER-CL12T	.5" x .94"	7 Color Choices	212	16
* LASER-CL23T	.5" x 1.28"	7 Color Choices	176	16
* LASER-CL48T	.4333" Dia.	7 Color Choices	120	16
LASER-CL60	.94" x 1.125" & .44" Dia.	White Only	100	16
SLS-LASER-CL32T	.78" x .94"	White Only	104	16
LASER-TRCL12	.5" x .94"	Transparent	84	16
LASER-TRCL23	.5" x 1.28"	Transparent	70	16

Printing Templates  
are available for  
all Laser Labels.


### \* LABELS AVAILABLE IN 7 COLORS



\* LASER-CL12T-6  
.5" x .94"  
Blue



\* LASER-CL23T-1  
.5" x 1.28"  
White




\* LASER-CL48T-2  
.4333" Dia.  
Yellow



\*To indicate label color when ordering, place the color number after the catalog number above.  
For instance, LASER-CL12T-1 in yellow would be ordered as LASER-CL12T-1-2.

### OPAQUE WHITE LABELS




LASER-CL60  
.94" x 1.125" For Tube  
.44" Dia. for Tube Tops



SLS-LASER-CL32T  
.78" x .94"  
For Microscope Slides

### TRANSPARENT LABELS



LASER-TRCL12  
.5" x .94"

LASER-TRCL23  
.5" x 1.28"

Clear, See Through Labels

# Specialty Laser Label Materials

## STANDARD LASER MATERIAL

**Laser Bright White** Smudge resistant material with general purpose, permanent adhesive providing excellent initial adhesion to most surfaces including plastic, glass and paper. May be applied at 25°F (-4°C), holds in temperatures ranging from -40°F to 300°F (-40°F to 149°C).

## LASER MATERIALS FOR SPECIAL APPLICATIONS

**Laser Removable** Standard weight stock for precise print quality, removable adhesive ideal for temporary labeling applications on most surfaces. May be applied at 32°F (0°C), holds in temperatures ranging from -40°F to 200°F (-40°F to 93°C).

**Laser Smudgeproof Kimdura®** Glove free, solvent resistant, aggressive permanent adhesive. Its surface accepts a clear imprint that won't smear or flake. May be applied at 35°F (2°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C).

**Laser Piggyback** Permanent adhesive will adhere to moist glass bottles. After 10 minute set-up time, will withstand immersion in ice water, provide UV and solvent resistance. May be applied at 33°F (1°C), holds in temperatures ranging from -40°F to 280°F (-40°C to 138°C).

# Specialty Labels for Dot Matrix Printers

Available in a wide range of sizes. Specialty labels run through your dot matrix printer without jamming, accept a clear imprint, stick tightly to most surfaces and adapt to a variety of labeling needs. Labels are supplied on computer litho stock. Special application label materials are available for all pinfed labels.



## GENERAL PURPOSE LABORATORY LABELS

General purpose pinfed labels are available in other sizes and in formats of three to five labels across. Contact Shamrock for complete information.

### BLANK LABELS -- ONE WIDE



SNO-52 1<sup>15</sup>/<sub>16</sub>" x 2<sup>3</sup>/<sub>4</sub>"

Cat. No.	Size	Perf Every	Cat. No.	Size	Perf Every
SNO-12	7/16" x 2 <sup>3</sup> / <sub>4</sub> "	6"	SNO-44	1 <sup>7</sup> / <sub>16</sub> " x 4"	1 <sup>1</sup> / <sub>2</sub> "
SNO-13	7/16" x 3 <sup>1</sup> / <sub>2</sub> "	6"	SNO-44-R (Removable)	1 <sup>7</sup> / <sub>16</sub> " x 4"	1 <sup>1</sup> / <sub>2</sub> "
SNO-14	7/16" x 4"	6"	SNO-45	1 <sup>7</sup> / <sub>16</sub> " x 5"	1 <sup>1</sup> / <sub>2</sub> "
SNO-15	7/16" x 5"	6"	SNO-701	1 <sup>1</sup> / <sub>2</sub> " x 3"	1 <sup>3</sup> / <sub>4</sub> "
SNO-17	7/16" x 7 <sup>3</sup> / <sub>4</sub> "	6"	SNO-420	1 <sup>11</sup> / <sub>16</sub> " x 3"	1 <sup>3</sup> / <sub>4</sub> "
SNO-105	3/4" x 3"	1"	SNO-52	1 <sup>15</sup> / <sub>16</sub> " x 2 <sup>3</sup> / <sub>4</sub> "	2"
SNO-559	3/4" x 4 <sup>1</sup> / <sub>2</sub> "	1 <sup>3</sup> / <sub>16</sub> "	SNO-153	1 <sup>15</sup> / <sub>16</sub> " x 3"	2"
SNO-110	1 <sup>5</sup> / <sub>16</sub> " x 2"	1"	SNO-53	1 <sup>15</sup> / <sub>16</sub> " x 3 <sup>1</sup> / <sub>2</sub> "	2"
SNO-31	1 <sup>5</sup> / <sub>16</sub> " x 2 <sup>1</sup> / <sub>2</sub> "	1"	SNO-54	1 <sup>15</sup> / <sub>16</sub> " x 4"	2"
SNO-32	1 <sup>5</sup> / <sub>16</sub> " x 2 <sup>3</sup> / <sub>4</sub> "	1"	SNO-55	1 <sup>15</sup> / <sub>16</sub> " x 5"	2"
SNO-700	1 <sup>5</sup> / <sub>16</sub> " x 3"	1"	SNO-373	2 <sup>7</sup> / <sub>16</sub> " x 3"	2 <sup>1</sup> / <sub>2</sub> "
SNO-33	1 <sup>5</sup> / <sub>16</sub> " x 3 <sup>1</sup> / <sub>2</sub> "	1"	SNO-64	2 <sup>7</sup> / <sub>16</sub> " x 4"	2 <sup>1</sup> / <sub>2</sub> "
SNO-34	1 <sup>5</sup> / <sub>16</sub> " x 4"	1"	SNO-175	2 <sup>7</sup> / <sub>16</sub> " x 4 <sup>5</sup> / <sub>8</sub> "	2 <sup>1</sup> / <sub>2</sub> "
SNO-35	1 <sup>5</sup> / <sub>16</sub> " x 5"	1"	SNO-60	2 <sup>15</sup> / <sub>16</sub> " x 4"	3"
SNO-37	1 <sup>5</sup> / <sub>16</sub> " x 7 <sup>3</sup> / <sub>4</sub> "	1"	SNO-65	2 <sup>15</sup> / <sub>16</sub> " x 5"	3"
SNO-42	1 <sup>7</sup> / <sub>16</sub> " x 2 <sup>3</sup> / <sub>4</sub> "	1 <sup>1</sup> / <sub>2</sub> "	SNO-465	3 <sup>1</sup> / <sub>4</sub> " x 1 <sup>15</sup> / <sub>16</sub> "	3 <sup>1</sup> / <sub>2</sub> "
SNO-43	1 <sup>7</sup> / <sub>16</sub> " x 3 <sup>1</sup> / <sub>2</sub> "	1 <sup>1</sup> / <sub>2</sub> "	SNO-193	3 <sup>15</sup> / <sub>16</sub> " x 6"	4"

### BLANK LABELS -- TWO WIDE



SNO-335 1<sup>7</sup>/<sub>16</sub>" x 3"

Cat. No.	Size	Perf Every
SNO-443	1 <sup>5</sup> / <sub>16</sub> " x 1 <sup>5</sup> / <sub>16</sub> "	1"
SNO-33-2	1 <sup>5</sup> / <sub>16</sub> " x 3 <sup>1</sup> / <sub>2</sub> "	6"
SNO-34-2	1 <sup>5</sup> / <sub>16</sub> " x 4"	6"
SNO-335	1 <sup>7</sup> / <sub>16</sub> " x 3"	1 <sup>1</sup> / <sub>2</sub> "
SNO-44-2	1 <sup>7</sup> / <sub>16</sub> " x 4"	6"
SNO-305	1 <sup>15</sup> / <sub>16</sub> " x 3"	2"


# Specialty Labels for Dot Matrix Printers



## LABELS FOR LABORATORY SOFTWARE




BDP-3 1/4" x 1" (Perf every 11")



BDP-2 7/16" x 1"

Most stock labels are white. Color may be added to your specifications.



ACTG-SCPF-1 5/8" x 1" (Perf every 11")



BDP-5 15/16" x 15/16"



BDP-1 15/16" x 15/16"



BDP-T-20 15/16" x 2 1/2"




BDP-52 15/16" x 3"



BDP-M-50 15/16" x 3 1/2"



BDP-CHC-30  
1 1/16" x 1 11/16"



BDP-337 1 7/8" x 8 1/16"

# Specialty Labels for Dot Matrix Printers

## STANDARD DOT MATRIX MATERIAL

**EDP Permanent** Permanent adhesive offers high initial tack and ultimate adhesion on most surfaces, will not flag or fall off curved surfaces down to 1/2" diameter. Will hold in temperatures ranging from -65°F to 220°F (-54°C to 104°C).

## DOT MATRIX MATERIALS FOR SPECIAL APPLICATIONS

**EDP Removable** General purpose removable adhesive provides a good initial tack and is easily removed on most surfaces including plastic and glass. May be applied at 40°F (4°C), holds in temperatures ranging from 0°F to 120°F (-18°C to 49°C).

**Pliaprint** General purpose permanent adhesive features high initial tack and ultimate adhesion. Excellent for labeling small containers, irregular and curved surfaces. May be applied at 40°F (4°C), holds in temperatures ranging from -65°F to 160°F (-54°C to 71°C).

**Dupli-Copy** Everything printed on label automatically copies onto the label's liner paper. Strong permanent adhesive is solvent resistant. May be applied at 35°F (1°C), holds in temperatures ranging from -65°F to 200°F (-54°C to 93°C).

**Twin Stock** Two matte finish labels in one. Printing on top label copies onto an underlying, second label. Permanent adhesive holds label firmly in place. May be applied at 25°F (-4°C), holds in temperatures ranging from -75°F to 200°F (-60°C to 93°C).


# Specialty Labels for Every Printer, All Software

Shamrock can supply labels blank or individualized to meet your special requirements, tinted, die cut to exactly match your software, on materials specific to your thermal, laser, or dot matrix printer.


If you don't see the exact label needed for your software, contact Shamrock at 800-323-0249. We offer a wide variety of stock and custom materials plus design assistance to meet your individual labeling requirements.




### THERMAL



### LASER



### DOT MATRIX



SHAMROCK LABELS 34 DAVIS DRIVE, BELLWOOD, ILLINOIS 60104

TOLL FREE: 1-800-323-0249 FAX: 1-800-248-1907 EMAIL: SALES@SHAMROCKLABELS.COM SHOP ONLINE: SHAMROCKLABELS.COM

## Terms and Conditions of Sale

All prices subject to change without notice.

A \$25.00 minimum order is required.

Prices are F.O.B. Bellwood, Illinois, unless otherwise specified. Foreign import duties and exchange rates, if any, are additional.

Terms: Net 30 days.

### RETURN POLICY

No returns will be accepted without prior authorization. Shamrock's numbered return authorization form must accompany all merchandise returned.

**Defective Merchandise:** Should any Shamrock product prove to be defective, a claim for return must be placed within 60 days of shipment. Products shown to be defective and returned within this period will be replaced.

**Stock Items:** No return of stock items will be accepted after six months. A 5% restocking charge is imposed on return of all stock items.

**Custom Items:** Custom orders may not be returned.

### QUANTITIES

Over runs or under runs, not to exceed 10% of quantities ordered, will constitute acceptable delivery.

### PLATE CHARGES

Due to natural deterioration of printing plates, and in order to provide only top-quality products, Shamrock must remake its plates every three years. If an order is not placed within that time, a new plate charge must be made upon receipt of the next order.



814





shamrock  
labels

**SHAMROCK LABELS**

34 Davis Drive, Bellwood, Illinois 60104

National Toll Free: 800-323-0249

California Toll Free: 877-461-6908

Fax Toll Free: 800-248-1907

Email: [sales@shamrocklabels.com](mailto:sales@shamrocklabels.com)

Shop Online: [shamrocklabels.com](http://shamrocklabels.com)